Deutsch

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English

1. Prior to start-up

You have acquired a valuable RIESTER diagnostic set manufactured in compliance with Directive 93/42/EC for medical products and subject to continuous stringent quality control, whose excellent quality will ensure reliable diagnoses.

These Instructions describe the use of RIESTER uni[®], I, II, III instruments and economy[®] and accessories. Please read these Instructions carefully prior to start-up and keep in a safe place.

Should you have any queries, please contact the Company or your RIESTER Agent who will be pleased to assist you. For addresses see last page of these Instructions. The address of your authorised RIESTER Agent will be supplied to you on request.

Please note that any instruments described in these Instructions are only suited for application by trained operators.

Please also note that correct and safe operation of instruments will only be guaranteed when RIESTER instruments and accessories are used throughout.

2. Battery handle and accessories

2.1. Purpose

The RIESTER battery handles described are used for supplying the instrument heads with power (the lamps are included in appropriate instrument heads), also serving as a bracket.

2.2. Battery handle range Important information
Any instrument heads described in these operating instructions will fit the following battery handles and may therefore be individually combined.

a) **Type C battery handle incl. a 2.5 V**
order to operate these battery handles, you will require two commercial type C alkaline baby batteries (IEC standard designation LR14) or one 2.5 V rechargeable battery (RI-accu[®]). The handle and the RIESTER rechargeable battery (ri-accu[®]) may only be charged by a RIESTER ri-charger[®] charging base.

b) **Type C battery handle incl. a 3.5 V**
This battery handle may only be operated with a RIESTER rechargeable battery of 3.5 V (ri-accu[®]). This may only be charged in a RIESTER ri-charger[®] charging base.

c) **Type C rechargeable battery handle of 2.5 V or 3.5 V with a rheostat for charging in a mains socket of 230 V or 120 V**
This handle is available both as a 2.5 V or 3.5 V version and may be ordered for 230 V or 120 V mains operation. The handle is supplied with a suitable rechargeable battery (ri-accu[®]). Please note that the handle must be operated with a RIESTER rechargeable battery (ri-accu[®])

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uni[®], economy[®]

2.1. Insertion and removal of battery handle types (2.2. a) and b)
Turn handle over on bottom handle section.
Depending on what handle has been purchased for what mains voltage (see 2.2), either insert 2 commercial alkaline type "C" baby batteries of 1.5 V (IEC standard designation LR 14) or RIESTER rechargeable battery (ri-accu[®]) of 2.5 or 3.5 V into the sleeve of the handle, with the plus pole pointing towards the top section of the handle. There is an additional arrow on the rechargeable battery next to the plus sign, indicating the direction for insertion into the handle. Reighten handle cover firmly over the handle.
Remove batteries by first releasing the battery handle cover, followed by shaking the handle lightly.
Prior to operation, insert rechargeable batteries (ri-accu[®]) into handle and charge with a RIESTER ri-charger[®] charging base. Each charging base is supplied with separate operating instructions which should be complied with.

Handle types (2.2. c)
Prior to first start-up of the mains socket handle, charge for 24 hours in a mains socket.
Attention: The plug-in handle should not be charged for longer than 24 hours.
Leave rechargeable battery (ri-accu[®]) in the handle as supplied. Rotate bottom handle section in counter-clockwise direction to release, thus exposing the socket contacts. Round contacts are used for 220 V mains operation, flat contacts for 110 V. Now insert bottom handle section into socket for charging.
Warning: Never leave handle in socket whilst the battery is replaced!
When changing the rechargeable battery (ri-accu[®]), rotate battery handle cover at handle bottom section in counter-clockwise direction to release. Remove rechargeable battery (ri-accu[®]) from battery handle by shaking the handle lightly downwards. Insert rechargeable battery (ri-accu[®]) into battery handle. Ensure that rechargeable batteries of 2.5 V are inserted into the handle with the polished metal head first. For 3.5 V rechargeable batteries, the direction of insertion is unimportant. Finally tighten battery handle cover in clockwise direction on handle. Mains voltage: Optionally 230 or 120 V

Note:
• When the unit is not being used for an extended period of time or when travelling, remove batteries from battery compartment.
• Insert new batteries when the light intensity of the instrument is reduced, thus possibly affecting examination.
• For maximum light it is recommended to insert two new high-quality batteries (as described in 2.2).
• Should you suspect that fluid or condensation has penetrated into the handle, do not recharge. This may cause lethal electrocution, in particular in handles connected to a mains socket.

Disposal
Please note that batteries must be disposed of separately. For information ask your local council and/or environmental officer.

2.4. Attaching instrument heads
Ensure that the knurled screw on the battery handle is released so far that the tip of the screw can no longer be seen in the interior of the handle.
Attach required instrument head to the support on the battery handle, ensuring that the notch is directed towards the knurled screw. Fix by the locking screw.

2.5. Switching the instrument on/off
Switch instrument on by pressing the red On/Off switch on the black knurled plastic ring (rheostat) and turning it to the left, away from the "0", whilst still keeping it pressed.
Switch instrument off by keeping the On/Off switch pressed and turning it to the right to position "0".

2.6. Rheostat for control of light intensity
The rheostat allows to set light intensity. Depending on how far the switch including the black knurled ring is being screwed in in counter-clockwise or clockwise direction, light intensity will be increased or decreased. For orientation, use marking below the rheostat.

2.7. Spares
Item No. 10426 Top handle section comprising a rheostat and plug-in locking system

Note: Comply with operating instructions!

Type B application section

3. The Otolaryngoscope and accessories

3.1. Purpose

The RIESTER otolaryngoscope described in these Instructions has been produced with examination of the auditory canal in mind, combined with a RIESTER ear speculum.

3.2. Application and removal of ear speculum
Position the required ear speculum on the head of the otoscope, with the recess in the metal section of the ear speculum fitting into the guide pin of the head. Lock by turning the speculum in clockwise direction. For removal of the speculum, turn fully in opposite direction, then remove from otoscope.

3.3. Swivel lens for enlargement
The otoscope comprises an enlarging lens, swivelling by 360° for approximately 4X enlargement. The swivel lens may be easily removed by withdrawal.

3.4. Introduction of external instruments into the ear
When intending to introduce external instruments into the ear (such as a pair of tweezers), remove glass cover with the black ring and the connection for pneumatic testing by turning in counter-clockwise direction and withdrawing.
Reposition glass. Ensure that the recess in the black plastic ring fits into the guide pin of the otoscope. Lock ring by rotating in clockwise direction.

3.5. Pneumatic testing
In order to perform a pneumatic test (= examination of the ear drum), you will require a ball not included in the standard set supplied, which may, however, be ordered as an option (see 3.8 Spares and accessories).

The otoscope head comprises a black outer ring and a connection. Attach the tube end of the ball to this connection, then carefully blow the required air volume into the auditory canal.

3.6. Replacement of lamp
Detach ear speculum from otoscope (see 3.2.). Release lamp by turning in counter-clockwise direction.
Insert new lamp by screwing in in clockwise direction and re-attach required ear speculum (see 3.2.).

3.7. Spares and accessories
Ear speculum
Item No. 10460 Reusable ear speculum, 2 mm
Item No. 10461 Reusable ear speculum, 3 mm
Item No. 10462 Ear speculum for multiple use, 4 mm
Item No. 10463 Ear speculum for multiple use, 5 mm
Item No. 10464 Ear speculum for multiple use, 9 mm

Replacement lamps (vacuum)
Item No. 10421 Pack of 6 lamps each of 2.7 V for otoscope

Replacement lamps (halogen)
Item No. 10590 Pack of 6 lamps each of 2.5 V for otoscope

Replacement lamps (xenon)
Item No. 10592 Pack of 6 lamps each of 3.5 V for otoscope

4. The ophthalmoscope and accessories

4.1. Purpose

The RIESTER May ophthalmoscope described in these Instructions has been designed for the examination of the eye and its background.

4.2. Lens wheel and correcting lenses
The correcting lenses may be adjusted on the lens wheel.
The following correcting lenses are available: diopters 0 to +20 and 0 to -20. Readings will be displayed on a lit panel. Plus values are displayed on a black background, minus values on a red background.

4.3. Aperture
There is one aperture (a round circle for standard fundus examinations).

4.4. Lamp replacement
Rotate knurled screw on ophthalmoscope in counter-clockwise direction and withdraw the socket together with the lamp. Rotate lamp in counter-clockwise direction to remove.
Fully screw in new lamp in clockwise direction. Insert socket into ophthalmoscope, ensuring that the screw fits into the recess provided below the knurled screw. Then fully screw in knurled screw.

4.5. Spares and accessories
Replacement lamps (vacuum)
Item No. 10424 Pack of 6 lamps of 2.7 V for May ophthalmoscopes
Replacement lamps (halogen)
Item No. 10590 Pack of 6 lamps of 2.5 V for May ophthalmoscopes
Replacement lamps (xenon)
Item No. 10593 Pack of 6 lamps of 3.5 V for May ophthalmoscopes

5. Bent arm illuminator and accessories

5.1. Purpose

The bent arm illuminator described in these Instructions is manufactured for lighting the oral cavity and the pharyngeal area.

5.2. Replacement of lamp
Screw down lamp at front of bent arm illuminator in counter-clockwise direction and screw new lamp into bent arm illuminator in clockwise direction.

5.3. Spares and accessories
Replacement lamps (vacuum)
Item No. 10421 Pack of 6 lamps of 2.7 V for bent arm illuminator
Replacement lamps (halogen)
Item No. 10590 Pack of 6 lamps of 2.5 V for bent arm illuminator
Replacement lamps (xenon)
Item No. 10592 Pack of 6 lamps of 3.5 V for bent arm illuminator

5. Bent arm illuminator and accessories

5.1. Purpose
The bent arm illuminator described in these Instructions is manufactured for lighting the oral cavity and the pharyngeal area.

5.2. Replacement of lamp
Screw down lamp at front of bent arm illuminator in counter-clockwise direction and screw new lamp into bent arm illuminator in clockwise direction.

5.3. Spares and accessories
Replacement lamps (vacuum)
Item No. 10421 Pack of 6 lamps of 2.7 V for bent arm illuminator
Replacement lamps (halogen)
Item No. 10590 Pack of 6 lamps of 2.5 V for bent arm illuminator
Replacement lamps (xenon