Standard-Instrumente /

Standard instruments / Instruments standard / Instrumentos estándar / Strumenti standard / Стандартни инструменти / Standardinstrumenter / Standardinstrumentid / Vakioinstrumentti / Πρότυπα όργανα / Standardni instrumenti / Standarta instrumenti / Standartiniai instrumentai / Standaardinstrumenten / Standardowe instrumenty / Instrumentos padrão / Instrumente standard / Standardinstrument / Štandardné nástroje / Standardni instrumenti / Standardní nástroje / Standard eszközök

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Seite 5

Operating Manual

Page 29

Mode d'emploi

Page 51

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Página 75

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Página 97

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

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Table of contents

1	Symbol descriptions	30
2	Introduction	31
3	Scope	31
3.1	Intended Use	34
3.2	Contraindications	41
4	Warnings	42
5	Handling	43
6	Preparation	43
6.1	Reutilization restrictions	43
6.2	Information on instrument preparation	43
6.3	Preparation at the place of use	43
6.4	Ultrasound bath (optional)	44
6.5	Manual cleaning	44
6.6	Mechanical cleaning	44
6.7	Drying	46
7	Maintenance, inspection	47
7.1	Functional check	47
8	Sterilization	48
8.1	Packaging	48
9	Lifetime	48
10	Storage	49
11	Warranty / Repair	49
12	Service and manufacturer address	49

1 Symbol descriptions

Symbols	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
  Hinweis auf eIFU	(Electronic) instruction for use

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

3 Scope

The use of Standard Instruments is mainly in the following areas:

Comedones Extractors; Plastic Spatula	The instrument is used for treatments. It may only be used by trained and qualified professionals.
Eye probes; Foreign Body instruments	The instrument is used for treatments in ophthalmology. It may only be used by trained and qualified professionals.
Rectal director	The instrument is used in gastroenterological and urological examinations. The examination may only be carried out by suitably trained and qualified specialist personnel.
Tongue depressor; Measuring probes	The instrument is used for ENT examinations and treatments. It may only be used by trained and qualified professionals.
Uterine Dilators	The instrument is used during examinations in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Dilatation bougies	The instrument is used during examinations in urology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Screwdriver; Hand Fixation Device; Bending forceps; Slat cutters; Hammers; Spatula	The instrument is used in various procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Marking instruments	The instrument is used in various examinations. The examination may only be carried out by suitably trained and qualified specialist personnel.
Cuticle- and Nail instruments; Cuticle Nippers; Pedicure nippers; Nail Scissors; Nail Nippers	The instrument is used for nail and foot care.
Areolamarkers	The instrument is used in plastic surgery. It may only be used by trained and qualified professionals.
Uterine depressants	The product is used in examinations and treatments in gynecology. It may only be used by trained and qualified personnel.

Laryngeal Mirrors; Laryngeal mirror handles; Middle Ear Mirror; Nasal specula (simple); ear specula; Probes; Cotton Applicators; Holder, Speculum; Lacrimal duct Dilators	The product is used in ENT examinations and treatments. It may only be used by trained and qualified professionals.
Extension Bows	The product is used in orthopedics. It may only be used by trained and qualified personnel.
Stainless steel bowls	The use of vessels takes place in almost every medical field. It may only be used by trained and qualified professionals.
Laminectomy Spreaders; Lamina Spreaders	The instrument is used for surgical procedures on the spine. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Gall Stone Probes	The instrument is used in surgical procedures of the gallbladder and bile duct. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Bile duct dilators	The instrument is used in surgical procedures of the bile duct. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Micro Instruments; Sharp Hooks; Myoma Screw; Biopsy curettes; Uterine curettes	The instrument is used for surgical procedures in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Eye Hooks; Strabismus Hooks; Exenteration Scoops; Chalazion Curettes; Lid Retractor; Eye Specula; Muscle and fixation hooks	The instrument is used for surgical procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Curettes; curettes, titanium; sharp spoons; ring curettes; Drill and guide protection forceps; Bending irons, bones	The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Lung Spatulas	The instrument is used for thoracic surgery procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Bladder spatula; Gallstone spoon; Uterine spoon	The instrument is used in gastroenterological and urological surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Trocars; Breast Enlargement Sets; Lupus curettes; lupus spoon; lupus scraper; Tissue spatula; Hooks; Exploration- and Coagulation Hooks; Root Retractor; Spreaders; Tissue hook; Ligature Guides; Ligature Conductors; Tendon passing instrument; Wire Strut Guide; Skin Hook; Micro Hooks, reparatory, spoons, forks; Intestinal Needles; Muscular and Skin Sutures; Reverdin Needles; Ligature needle; Guide needles; Laser Instruments, Hook	The instrument is used in various surgical procedures. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Abdominal and Intestinal Spatulas	The instrument is used in various surgical procedures in the abdominal area. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Elevators; Meniscus retractors	The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Tissue dissectors; Separator	The instrument is used for various procedures in general or plastic surgery. It is not intended for use in direct contact with the heart, central circulatory system or central nervous system. The procedure must be performed by trained and qualified personnel.
Vascular dilators	The instrument is used in various procedures in vascular surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Endarterectomy strippers; Varicose Vein Probe Set;	The instrument is used in various procedures in vascular surgery. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Vein instruments, Vein retractor, vein extractor	The instrument is used during surgical ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Septum and Periosteal Elevators; Hooks for posterior nares; Ears Curettes; Sinus Curettes; Ear Loops; Nasal Polypus Snares; Tonsil Snares; Ear Hooks; Paracentesis Needles; Foot plate instruments; Tympanum Needle	

3.1 Intended Use

Plaster Spatula	An instrument, usually made of stainless steel, for applying plaster. It is a non-invasive, reusable instrument.
Comedones Extractors	A hand instrument for removing sebaceous deposits and comedones. It may be a single- or double-ended instrument with blunt working end and a handle positioned centrally or at the proximal end. It is usually made of stainless steel. It is a non-invasive, reusable instrument.
Foreign Body instruments	An ophthalmic instrument to remove foreign bodies from inside the eye or from the surface of the eyeball. It is typically made of high quality stainless steel and has a fine handle at the proximal end and a thin concave blade with a rounded nose at the distal end. The ophthalmologist uses it to remove the foreign body / object. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Eye probes	A slender rod-shaped instrument made of flexible metal, with a blunt tip of tuber, used to examine eye tissue. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Rectal director	An instrument used for exploring the tissues and structures of the rectum during a gastroenterological / urological examination. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Tounge depressor	A hand-held tongue removal instrument to facilitate examination of the surrounding organs / tissues. It is typically a double-ended, flat product of various sizes with a slightly angled working end to better manipulate the tongue. It is made of stainless steel. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Measuring probes	A product / instrument used for comparative measurement, e.g. to measure inner and outer diameters, lengths, depths or thicknesses. The instrument ist not calibrated. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Uterine- Dilators	A round rod-like massive instrument that is inserted through the cervix and dilates the cervical canal. Commonly known as a uterine dilator, it is typically made of metal or plastic and has a tapered rounded distal tip. It can have one or two different diameters, one at each end of the instrument, and can be part of a set of size-increasing diameters which are inserted successively. The instrument allows insertion of a curette or other instrument into the uterus and / or expands the uterus by pushing against the walls of the uterus from the inside. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.

Dilatation bougies	A round, rod-like instrument for examining and dilating the urethra. It has a semi-rigid, conical shaft that merges into a narrow, bulbous finish. The instrument can be made hollow to pass other, narrow instrument can. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Spatula	A hand instrument, usually in spoon shape or rounded at the end of work for preparatory activities in various procedures. The instrument is made of high quality stainless steel. It is a non-invasive, reusable instrument.
Hammers	A manual instrument for applying force (e.g., by beating) to another product (e.g., a chisel, pry bar, driver) during a procedure. It is available in different designs and materials. The shaft is usually made of metal (chromed iron or stainless steel), but may also be made of synthetic materials. The hammer head is usually solid, annular, mostly double-ended, and is made of solid (iron / steel) or a shock-absorbing material (e.g., rubber or plastic). It is a non-invasive, reusable instrument.
Bending forceps; Slat top cutters	An instrument for external bending of wires or plates during various procedures. It is a non-invasive, reusable instrument.
Hand Fixation Device	A device for fixing the hand. It is a non-invasive, reusable product.
Screwdriver	A tool with a shaft that fits into a screw head to tighten or remove this screw. The working end may e.g. for a slotted screw, Phillips screw, Phillips, Hexagon. It can have an indication of the torque. It is a non-invasive, reusable product
Marking instruments	An instrument, usually in the form of a suitable pen containing non-toxic ink, used to demarcate areas on a patient's skin (e.g., incisions and for dermatological applications). It is a non-invasive, reusable product.
Cuticle- and Nail instruments; Cuticle Nippers; Pedicure nippers; Nail Scissors; Nail Nippers	An instrument for cutting fingernails and toenails. It can be bent to match the natural curvature of human nails. The toenail models are usually made more robust to cut thick toenails. The instrument consists of two movable branches, with holes for fingers and thumbs in the handle, and cuts by closing the leaves over the nail. It is a non-invasive, reusable instrument
Areolamarkers	An instrument used to mark on the surface of the skin, e.g. the areola. It is a non-invasive, reusable product.
Uterine depressants	A hand-held instrument for holding off the uterus in order to facilitate an examination of the surrounding organs / tissues. It is typically a double-ended product in various designs (e.g., curved, spoon-shaped, or fenestrated) with a slightly angled working end to better manipulate the organs. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.

Nasal specula (simple); ear specula	An instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used in insertion and opening to stretch or stretch tissue in the nostril. It is used for examination or during procedures. It is usually made of stainless steel with a self-opening mechanism or an adjustable mechanism. It is a reusable instrument intended for transient use. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Holder, Speculum	A device for fixing an ear, nose and throat (ENT) speculum in a fixed position during examinations and treatments. It is a non-invasive, reusable instrument.
Probes; Cotton Applicators	A narrow, rod-like, flexible metal instrument, usually with a blunt, beaded, or pointed tip, for examining fistulas, lumens during an ENT procedure. The distal end, in front of the bead / tip, may be bent at right angles to the shaft. It is a reusable instrument intended for transient use.
Laryngeal Mirrors; Laryngeal mirror handles; Middle Ear Mirror	A hand instrument with a polished surface, for reflecting sufficient undiffused light to obtain a visible image of an object during an ENT procedure. This mirror is mounted on a long, narrow handle, which is held by the doctor in the vicinity of the object to be observed. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Lacrimal duct Dilators	An ophthalmological instrument for use during ophthalmological examinations to dilate a cavity, canal or opening. Eye dilators also include instruments used for lacrimal or lacrimal canal surgery. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Extension Bows	The product is part of a traction unit. It consists of metal bars with a curved frame, for attaching hooks or pins for exerting a pull on the patient. The traction unit exerts a tensile force to allow the body parts. It is a non-invasive, reusable product.
Stainless steel bowls	A container, usually in the form of a bowl, that is used to store medical materials. It is a non-invasive, reusable product.
Laminectomy Spreaders; Lamina Spreaders	A self-locking surgical instrument for separating the edges or withdrawing the wound edges of a surgical incision to permit the removal of a lamina of the vertebral arch. This may be necessary to remove a tumor, detachment, or tension caused by fracture or disc protrusion. The instrument resembles a large scissors with bulging branches with fenestrated, curved hooks that often end sharply to grab tissue. It is usually made of stainless steel and is available in different sizes and designs. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.

Gall Stone Probes	A surgical instrument consisting of a slender rod of flexible metal with a blunt buttoned tip used to examine (bite) the bile duct, e.g. is used on the union of the ductus sympatheticus with the ductus cysticus. It is used during surgical procedures to examine or dilate bile duct narrowing. It is a reusable instrument intended for transient use.
Bile duct dilators	A surgical instrument consisting of a slender hollow or solid body of metal or plastic in a cylindrical shape in a variety of sizes for dilating the bile duct, e.g. at the union of the ductus sympatheticus with the ductus cysticus. It is a reusable instrument intended for transient use.
Micro Instruments; Sharp Hooks	A shaft-like surgical instrument whose distal end has the form of a hook or more parallel hooks for fixing or pulling uterine fibroids during a gynaecological procedure. It is typically a one-piece stainless-steel instrument with a handle in various designs at the proximal end. It is a reusable instrument intended for transient use.
Biopsy curettes; Uterine curettes	A surgical instrument with an elliptically windowed end on a long shaft for scraping out the uterine lining, scraping tissue during a gynecological procedure, and removing waste from the uterus. The instrument is usually made of stainless steel and can be either rigid or flexible. The working end is prepared sharp or blunt in different versions. The general product group includes uterine, placental and fallopian tube curettes. It is a reusable instrument intended for transient use.
Myoma Screw	A manual surgical hand instrument for detecting and holding a fibroid (i.e., a tumor of muscle tissue, usually in the uterus) during its surgical removal. It is typically a one-piece instrument, similar to a corkscrew with a pointed, spindled working end (for screwing into the myoma) and a long, slender shaft and a T-shaped handle. It is usually made of stainless steel. It is a reusable instrument intended for transient use.
Eye Hooks; Strabismus Hooks	An ophthalmic surgical instrument with a shaft-like handle in various versions that tapers towards the working end. The distal part of the instrument is blunt or pointed bent or angled to better manipulate the ocular tissue or to remove debris. It is a reusable instrument intended for transient use.
Lid Retractor; Eye Specula; Muscle and fixation hooks	A surgical eye instrument used to hold off the edges and associated tissue of an incision in the eye during a surgical procedure. It usually has a self-holding arm with curved hooks or blades at the distal end that grip and hold the tissue away. The proximal end is usually connected and has a lockable mechanism to hold the working end in an expanded position. It is made of stainless steel and is available in various sizes and designs. It is a reusable instrument intended for transient use.
Exenteration Scoops; Chalazion Curettes	An ophthalmic surgical instrument with a windowed, spoon-like or ring-shaped working end that is either sharp or blunt. It is used to scrape out eye tissue. It is usually made of stainless steel. It is a reusable instrument intended for transient use.

Curettes; curettes, titanium; sharp spoons; ring curettes	A hand-held surgical instrument used to cut and excise bone tissue, typically during orthopedic or plastic surgery procedures. It can have a spoon, ring or fenestrated working part. It is used to remove bone tissue. It is a reusable instrument intended for transient use.
Drill and guide protection forceps	A surgical hand instrument in the form of a cannula for simultaneously placing and anchoring a rotating drill (i.e., a surgical twist drill) through hard tissue. Surrounding soft tissue is protected while the drill is in operation. It usually has ridges at the distal end to give a secure placement on the tissue surface. At the proximal end is a handle for holding and moving. The product is usually adjusted to the drill size and is made of stainless steel. It is a reusable instrument intended for transient use.
Bending irons, bones	A surgical product for bending bone plates. It is a reusable instrument intended for transient use.
Lung Spatulas	A surgical instrument, usually made of stainless steel, which is used to manipulate lung tissue surfaces. The instrument consists of a proximal handle and a flat blade with no sharp edges at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. It is a reusable instrument intended for transient use.
Bladder spatula; Gallstone spoon; Uterine spoon	A spoon-shaped surgical instrument used to remove tissue during gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Guide needles	A product for connecting a catheter to another product (e.g., a plug) or for introducing a flushing or drainage catheter into the body. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Intestinal Needles; Mucular and Skin Sutures; Reverdin Needles; Ligature needle	A cylindrically shaped, solid surgical instrument of various diameters, whose distal tip can be made in various geometric shapes. It is used for introducing and pulling suture through superficial tissue, e.g. Skin or muscle, used. The suture serves to piece two wound edges together. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Trocars; Breast Enlargement Sets	A surgical instrument with a sharp pyramidal or conical tip used to pierce body cavity. It is usually assembled and used with a compatible sleeve, thereby allowing insertion of the instrument. After piercing, the trocar is withdrawn and provides a working channel into the body cavity. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.

Ligature Guides; Ligature Conductors; Tendon passing instrument; Wire clips or claws for holding the material while passing through the tissue; Strut Guide	A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, passing instrument; Wire clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Laser Instruments, Hook	A surgical instrument used in conjunction with a suitable endoscope during endotherapy. It is used for mechanical work, e.g. grasp tissue or foreign objects or width of an area / lumen. It works without electricity, including e.g. High-frequency, electromagnetic, ultrasonic or laser energy. It is a non-active product that is reusable and intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Skin Hook; Micro Hooks, raparatory, spoons, forkshandle	A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks that apply tension to the skin. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Hooks; Exploration- and Coagulation Hooks; Root Retractor; Spreaders; Tissue hook	A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks to hold tissue. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Tissue spatula	A surgical instrument, usually made of stainless steel, that is used to manipulate tissue or to remove material from a surface. The instrument consists of a proximal handle and a flat blade with no sharp edges at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Lupus curettes; lupus spoon; lupus scraper	A spoon-shaped surgical instrument used to remove tissue during gastroenterological and urological surgical procedures. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Abdominal and Intestinal Spatulas	A surgical hand instrument for manipulating the intestine and abdominal organs during an abdominal surgical procedure. It typically has a large, flat, blunt blade (without sharp edges) at the distal end and a strong handle at the proximal end. It is usually made of stainless steel and is available in various designs, e.g. curved and flat, or straight and flat with a rounded end, or right angled. Some models have a ridged leaf surface to better hold the slippery abdominal tissue. The handle is typically designed so that it can be kept non-slip. It is a reusable instrument intended for transient use.

Elevators; Meniscus retractors	A hand-held surgical instrument used to lift, position, or lever open bone structures, other anatomical structures or surgical material during an orthopedic surgery. It is made of stainless steel and is available in a variety of sizes, styles and working ends. It is a reusable instrument intended for transient use.
Tissue dissectors; Sepa-rator	A surgical hand instrument, usually in spoon shape or rounded at the end of work. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. It is a reusable instrument intended for transient use. Es ist nicht zur Verwendung in direktem Kontakt mit dem Herz, dem zentralen Kreislaufsystem oder dem zentralen Nervensystem vorgesehen.
Vascular dilators	A surgical instrument that is inserted into a vessel and used to dilate or adjust a vessel during a surgical procedure. The instrument is not intended for use on the blood vessels arteriae pulmonales, aorta ascendens, arcus aortae, descending aorta to bifurcation aortae, coronary arteries, common carotid artery, external carotid artery, internal carotid artery, cerebral arteries, brachiocephalic trunk, venae cordis, vena pulmonary, superior vena cava, and inferior vena cava. It is a rod-like instrument with a long, slender solid body and a distal tip in various shapes and sizes. It is available in a variety of sizes and degrees of softness and is typically made of metal. It is a reusable instrument intended for transient use.
Endarterectomy strips; Varicose Vein Probe Set; Vein instruments, Vein retractor, vein extractor	A surgical instrument to cut a piece out of a vessel. A stripper can be used to excise part of a vein or artery except for the arterial arteries pulmonary, ascending aorta, aortic arc, descending aorta, aortic bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, arteriae cerebrales, brachiocephalic trunk, venae cordis, vena pulmonales, superior vena cava, and inferior vena cava. Vessel Stripper is available in two different versions: 1) a flexible stainless steel cable with a peeling dome or disc at one end and a guide tip at the other end; and 2) a rigid rod terminating in a closed loop or loop, e.g., an external stripper. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Ear Loops; Nasal Polypus Snares; Tonsil Snares	A surgical ENT instrument, usually consisting of a canal or cannula and two fixed and sliding finger rings. A sling positions a circle of round, flexible wire around tissue in the ENT area and holds it in place. By moving the sliding finger ring you can control how much wire slides through the channel. The instrument is commonly used to exfoliate depleted foreign bodies. It is a reusable instrument intended for transient use.
Hooks for posterior nares	A shaft-like surgical instrument whose proximal end forms a handle that can vary in shape. The shaft tapers toward the distal end toward a hook or hooks used in plastic nose surgery. In the multi-hook variant, the hooks are arranged in parallel. It is a reusable instrument intended for transient use.

Ear Hooks; Paracentesis A fine surgical instrument that consists of a handle that tapers into a shaft and ends in a fine, 90 ° bent tip. It is used during ENT surgical procedures to repair structures in the middle ear, e.g. to edit the stapes. It is a reusable instrument intended for transient use.

Ears Currettes; Sinus Cu-A surgical hand instrument having an annular, fenestrated, or spoon-shaped working portion which may be either sharp or blunt to provide tissue by scraping or scraping during an ENT procedure. It is a curette for general use in ENT surgery, it is made of stainless steel and is available in various designs and shapes. It is a reusable instrument intended for transient use.

Septum and Periosteal Elevators A surgical hand instrument for lifting, positioning and rupturing anatomical structures or surgical material in the ENT area. It can be a single or double instrument with a blunt end (leaf or hook in different versions) and a handle in the middle. Size and design of intended use, e.g. a dull impeller elevator, a septum elevator, a stronger nose artery elevator, a special viable ear elevator. The instrument is usually made of stainless steel. It is a reusable instrument intended for transient use.

3.2 Contraindications

No contraindications are known.

4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Please observe the additional information enclosed with the products.
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use.

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do **not** use hard brushes (e.g. metal brushes or metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

V The drying time during validation was 1 hour.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through. Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

V To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection. Medicine products which are possessing a bad ultrasound transmission, e.g. soft materials are not usable for the ultrasound bath.

6.5 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

6.6 Mechanical cleaning

Due to international standards (EN ISO 15883) and national guidelines, only validated mechanical cleaning and disinfection procedures should be used. For automated cleaning, we recommend a standard program for surgical instruments, e.g. instruments from Miele.

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

Flexible (complex) instruments with non-visible surfaces should be pre-cleaned manually before machine cleaning.



We recommend manual pre-cleaning for all instruments where the surfaces lie on top of each other during cleaning (e.g. translated bone splitter and gouge forceps) in order to achieve an optimal, residue-free cleaning result.

Observe the following by loading:

- Place the disassembled or opened instruments securely in the tray.
- Instruments with openings and recesses must be placed with the open side facing downwards so that they can be cleaned and no water from the cleaning process can collect in them.
- If available, use a coordinated rinsing device.

- Do not overload trays, avoid creating any overlaps

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable.

After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between **80 and 95°C**, with an **exposure time** as outlined in **EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V	Drying was omitted in the validation (worst case condition).
----------	--

7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 7.1.

Instruments that are stained, blunt, bent, no longer functional or damaged in any other way must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement. Carry out a function check in accordance with the intended application of the instrument.

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132°C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

V Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	Drying time

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

V	During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.
----------	---

9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Website: www.HEBUmedical.de

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