

Air Compressors ONE - TWO SURGERY UNITS

Promed now offer a range of air compressors from market leaders such as Durr, Cattani and MGF. Please ask your account manager for further information.



If you are not sure which air compressor to choose, call us, our Service Engineers will be happy to advise you!

MGF COMPRESSORS SHOWN HERE ARE FOR ONE - TWO SURGERY UNITS



OIL FREE COMPRESSOR 30/10 GENESI SCE

- Bycylindrical oil-free compressor 1,15 kW – 230V 50Hz
- FAD@5bar: 120 l/min
- Regulated pressure (min-max): 5 -7 bar
- Air tank: 30 l.
- Noise level: 69 dB/lm
- Weight: 34 Kg
- Size: 43x57x62 cm
- Filtration: 10 µm + drying system with -20°C dew point value

OIL FREE COMPRESSOR SKY 30/10 GENESI SCE

- Bycylindrical oil-free compressor 1,15 kW – 230V 50Hz
- FAD@5bar: 120 l/min
- Regulated pressure (min-max): 5 -7 bar
- Air tank: 30 l.
- Noise level: 62 dB/lm
- Weight: 42 Kg
- Size: 54X65X76 cm
- Filtration: 10 µm + drying system with -20°C dew point value

OIL FREE COMPRESSOR CS 30/15 PRIME SCE

- Bycylindrical oil-free compressor 1,5 kW – 230V 50Hz
- FAD@5bar: 152 l/min
- Regulated pressure (min-max): 5 -7 bar
- Air tank: 30 l.
- Noise level: 54 dB/lm
- Weight: 89 Kg
- Size: 45X70X70 cm
- Filtration: 10 µm + drying system with -20°C dew point value



K100 ONE-CYLINDER COMPRESSOR FOR 1 SURGERY

Oil-free compressor with dry air system mounted on a 25ltr tank

Performance
3.8 CFM F.A.D. 67N l/min at 5 real bar

Dimensions
Width = 505mm
Depth=550mm
Height=735mm

Noise level
70dB (A)



K200 TWO-CYLINDER COMPRESSOR FOR 2-3 SURGERIES

Oil-free compressor with dry air system mounted on a 50ltr tank

Performance
8.2 CFM F.A.D. 160N l/min at 5 real bar

Dimensions
Width = 560mm
Depth=590mm
Height=760mm

Noise level
71.8dB (A)

CATTANI AND DURR DENTAL COMPRESSORS SHOWN ARE FOR ONE - TWO SURGERY UNITS

This information is just to give an overview of the range of equipment Promed can provide. For further information on our range of compressors and suction units please contact your account manager on freephone 1800 619 619



TORNADO 70 - FOR 1 SURGERY

The Tornado 70 is the smallest variant as an entry-level model with one cylinder and an output of 70 l/min.



PRIMO - FOR 1 SURGERY

The Primo is the high performance unit with one cylinder, with an output of 55 l/min specially developed for the one workstation practice.



DUO FOR 2 SURGERIES

The Duo with its 2-cylinder V-aggregate is the classic among the compact compressors with a 20-litre tank.

TORNADO 130 - FOR 2 SURGERIES

The Tornado 130 with dry air is the economic entry-level model in the Silver Airline family.