

## Dental Suction Units ONE - TWO SURGERY UNITS

Promed now offer a range of suction units from market leaders such as Durr, Cattani and MGF. Please ask your account manager for further information.



Promed's Service Engineers will be able to best advise you on which suction unit would best suit your requirements, and offer tips on best practice for their maintenance. If you are undertaking a complete surgery set up, please call one of the Promed Account Managers as they will be able to advise you of our Introductory offers and any available discounted pricing.

All units will be installed by Promed Service Engineers who can provide equipment training to staff, if required.

**Dürr Dental has stood for progress and innovation in dental technology for over six decades.**

The Dürr dental range of suction motors has both semi wet and dry suction motors from 1 surgery units to 3 surgery units.

This information is just to give an overview of the range of equipment Promed can provide. For further information on our range of compressors and suction units available please contact your account manager on freephone 1800 619 619

## ONE - TWO SURGERY UNITS



### VSA 300S - SUITABLE FOR ONE SURGERY

A compact and powerful combined vacuum pump and amalgam separator suitable for one surgery, compliant with regulations and with a proven track record. Easy to install in surgeries configured to traditional wet line systems or if an upgrade to amalgam separation is required.

- Over 99% amalgam separation rate
- Cleaning system/spare collection canister included as standard
- Easy to maintain when used with Dürr Orotol



### VS 600 - SUITABLE FOR TWO SURGERIES

VS 600 The VS 600 is a two surgery suction motor with combined automatic separation. Two-in-one is the magic formula. Dürr Dental has once again been successful in combining two machines: a suction motor and an auto separator, which separates the mixture of secretions and solids. For a maximum of two surgeries at the same time, even with high liquid flow rates. The VS 600 can be easily upgraded to VSA 600.

- Dimensions: (H) 500 x (W) 420 x (D) 410mm

The success of CATTANI around the world can be attributed to superior quality; from design, to the selection of materials, to the final production. No aspect of their manufacture is neglected; pistons and cylinders are individually checked to ensure conformity to parameters within microns. It is for reasons such as these that they give a 5-year warranty on every compressor they sell.



### CATTANI TURBOSMART - SUITABLE FOR TWO SURGERIES

- Consistent suction with variable speed
- High airflow and controlled heat output
- Easy upgrade from a 2 to a 4 surgery system
- Small footprint and energy saving
- Smart motor starting

### TWO TURBOSMART MODELS ARE AVAILABLE:

**TURBOSMART A**  
A maximum speed of 4800 rpm (for two surgeries running simultaneously)



All units will be installed by Promed Service Engineers who can provide equipment training to staff, if required.

## Dental Suction Units THREE -FOUR SURGERY UNITS

Promed now offer a range of suction units from market leaders such as Durr, Cattani and MGF. Please ask your account manager for further information.



### CATTANI TURBOSMART B - SUITABLE FOR FOUR SURGERIES

- Consistent suction with variable speed
- High airflow and controlled heat output
- Easy upgrade from a 2 to a 4 surgery system
- Small footprint and energy saving
- Smart motor starting

A maximum speed of 6200 rpm (for three to four surgeries running simultaneously) This machine is identical to the Turbosmart A, the only difference being that the models have different maximum speeds



### VS 900 - SUITABLE FOR UP TO THREE SURGERIES

Semi-wet suction motor like the VS 600 however for up to three surgeries.

- Dimensions: (H) 550mm x (W) 5420mm x (L) 410mm
- Noise: 69 db(A)
- Weight: 35kg (77lbs) without housing

### CA 4 - SUITABLE FOR UP TO FOUR SURGERIES

As an integral component of our state-of-the-art suction systems, the new Durr CA 4 satisfies the greatest demands of a modern dental practice.

- Centrifugal separator with hydro-dynamic self cleaning central operation with up to 4 patients simultaneously.
- Rate of flow of up to 16 litres per minute
- Optical and acoustic displays of the operating status
- Convenient removal of the amalgam collector

MGF COMPRESSORS SRL specializes in compressed air applications. Their production is composed of industrial, dental, medical and professional equipment that offers a wide range of solutions for applications where top silence & high reliability is required."

If you are not sure which dental suction unit to choose, call us, our Service Engineers will be happy to advise you!



Promed has a wide range of dental suction units from MGF, with varying options, with or without acoustic cabinets and canisters. Some variations are shown below.



All units will be installed by Promed Service Engineers who can provide equipment training to staff, if required.



COMPLETE SUCTION PUMPS  
WITHOUT ACOUSTIC CABINET



COMPLETE SUCTION  
PUMPS WITH ACOUSTIC  
CABINET WITHOUT  
CANISTER



COMPLETE SUCTION PUMPS  
WITH ACOUSTIC CABINET  
AND CANISTER