

sphygmotensiophone



Extreme precision and refined functionality – our two-tube sphygmomanometer. The chromium-plated pressure gage housing fits firmly on the cuff, is constantly in the field of view and is neat and tidy. The high-quality cotton cuff finished in artificial leather is provided with a metal clip through which to pull the cuff. This ensures safe and accurate measurement at all times.

- › Precision air release valve with fine adjustment and virtually wear-free
- › Microfilter protects air release valve and measuring system
- › Specially tempered copper-beryllium diaphragm almost non-aging
- › Diaphragm pressure-loading capacity up to 600g mmHg
- › Maximum error tolerance of +/- 3 mmHg
- › Easy to read face with scale of Ø 49 mm from 0 to 300 mmHg
- › latex-free version available (please enter 'LF' when ordering for latex free version)
- › Includes vinyl zipper bag

No. 1380-100	sphygmotensiophon manometer, 2-tube, w/o cuff
No. 1381	sphygmotensiophon manometer, 2-tube, velcro cuff (children) pull through strap
No. 1380	sphygmotensiophon manometer, 2-tube, velcro cuff (adult) pull through strap

110035-50 Rev. C 2023-01-23

Please visit riester.de for more information