

sanaphon®

Sphygmomanometer for Self-measurement



sanaphon® is a classic palm-style sphygmomanometer with integrated stethoscope for the precise self-measurement of blood pressure.

- › Sturdy and light aluminium casing, corrosion protected
- › Plug-in tube connector for cuff tube
- › Precision-release valve can be finely adjusted without wear and tear
- › Micro-filter protects pressure-release valve and measurement system
- › Spoon-shaped handle made of stainless steel
- › Special hardened copper-beryllium membranes, virtually non-wearing
- › Membrane can withstand overpressures of up to 600 mm Hg
- › Easy-to-read linear scale diam. 64 mm (2.5 in) from 0 - 300 mm Hg
- › Maximum error tolerance: +/- 2 mm Hg
- › In padded nylon zipper bag
- › Adult, large adult and child velcro cuff sizes with pull-through strap and integrated chest piece
- › Latex free stethoscope with soft ear-tips

No. 1442-100	sanaphon® manometer, w/o cuff, w/o stethoscope
No. 1442	sanaphon®, 1-tube with stethoscope, velcro cuff (adult) pull through strap
No. 1442-142	sanaphon®, 1-tube with stethoscope, velcro cuff (large adult) pull through strap
No. 1442-144	sanaphon®, 1-tube with stethoscope, velcro cuff (children) pull through strap

* Sphygmomanometers for self-measurement are alternatively available with **precisa®** N ABS, **e-mega®** or **ri-san®** gauges.

Please visit riester.de for more information