

## SPECIFICATION

### Silcoat Balloon Catheter



The latex balloon catheter with a silicone coating for medium indwelling time. The silcoat catheter has a pure latex core with an additional silicone-elastomer coating both inside and outside to minimise encrustation, at comparatively low costs. The catheter can remain in place for up to 14 days which is longer than associated with regular latex catheters.

- With optimised colour-coded plastic valve.
- With smoothly rounded eyelets
- Easily inflatable balloon (3 – 5 ml / 5 - 15 ml)
- Available in sizes Ch. 8 - 26

## SPECIFICATION

### Silcoat Balloon Catheter

Type	Size	PZN	REF	Unit
Type Paediatric - Balloon size 3 - 5 ml	CH 08	04639118	I7 180-08	10 pcs.
Type Paediatric - Balloon size 3 - 5 ml	CH 10	04639124	I7 180-10	10 pcs.
Type Adult - Balloon size 5 - 15 ml	CH 12	00368987	I7 605-12	10 pcs.
Type Adult - Balloon size 5 - 15 ml	CH 14	00368993	I7 605-14	10 pcs.
Type Adult - Balloon size 5 -15 ml	CH 16	00369001	I7 605-16	10 pcs.
Type Adult - Balloon size 5- 15 ml	CH 18	00369018	I7 605-18	10 pcs.
Type Adult - Balloon size 5 - 15 ml	CH 20	00369024	I7 605-20	10 pcs.
Type Adult - Balloon size 5 - 15 ml	CH 22	00369030	I7 605-22	10 pcs.
Type Adult - Balloon size 5 - 15 ml	CH 24	00369047	I7 605-24	10 pcs.
Type Adult - Balloon size 5 - 15 ml	CH 26	00369053	I7 605-26	10 pcs.