

**Actimove**<sup>®</sup>

Feel good. Move better.

# *Walker Range*

Focused Protection and Stabilisation  
for the Lower Leg, Ankle and Foot



# High quality support for every step

Actimove® Standard Walker,  
Actimove® Standard Air Walker and  
Actimove® Closed Shell Air Walker

The Actimove Walker range ensures protection of the foot, ankle and lower leg. Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort. Walking with Actimove Walker products feels natural due to the rocker bottom which simulates a physiological gait.

## Standard Walker

Secure and lightweight stabilisation for the lower leg, ankle and foot



## Standard Air Walker

Adjustable stabilisation with integrated pump and valve system to easily alter desired compression level



## Closed Shell Air Walker

Designed to offer the highest protection for the lower leg, ankle and foot



### Indications:

- Fractures of the foot and ankle
- Lower extremity soft-tissue injuries
- Severe ankle and foot sprains
- Post-operative and rehabilitation use

### For Closed Shell Air Walker only:

- Achilles tendon ruptures (only in combination with separately available Actimove Heel Wedge for Walkers)

## Actimove Walker range advanced features for safety and patient care



### **Circumferential strapping**

For critical stability



### **Shock absorbing sole and rocker bottom**

Minimises impact and simulates a natural gait



### **Comfortable breathable liner**

## Integrated pump system of Actimove Standard Air Walker and Actimove Closed Shell Air Walker



### **Integrated pump and valve system**

Designed to easily adjust desired compression level



## Additional adjustable features of the Actimove Closed Shell Air Walker



### **Heel wedge**

Adjustable heel wedge for Achilles tendon injuries




### **Adjustable toe guard**

Provides protection in the toe area







## Actimove Walker

	Essity Product Code		Size	Measurements UK Shoe Size	NHS SC Code	
	Low Walker	High Walker			Low Walker	High Walker
	76270-30	76270-35	X Small	< 3.5	ERE270	ERE275
	76270-31	76270-36	Small	4 - 6.5	ERE271	ERE276
	76270-32	76270-37	Medium	7 - 9.5	ERE272	ERE277
	76270-33	76270-38	Large	10 - 12	ERE273	ERE278
	76270-34	76270-39	X Large	> 12	ERE274	ERE279

## Actimove Walker Air

	Essity Product Code		Size	Measurements UK Shoe Size	NHS SC Code	
	Low Walker	High Walker			Low Walker	High Walker
	76272-30	76272-35	X Small	< 3.5	ERE328	ERE334
	76272-31	76272-36	Small	4 - 6.5	ERE333	ERE335
	76272-32	76272-37	Medium	7 - 9.5	ERE326	ERE330
	76272-33	76272-38	Large	10 - 12	ERE332	ERE329
	76272-34	76272-39	X Large	> 12	ERE331	ERE327

## Actimove Walker Air Closed Shell

	Essity Product Code		Size	Measurements UK Shoe Size	NHS SC Code	
	Low Walker	High Walker			Low Walker	High Walker
	76276-10	76276-15	X Small	< 3.5	ERE314	ERE319
	76276-11	76276-16	Small	4 - 6.5	ERE315	ERE320
	76276-12	76276-17	Medium	7 - 9.5	ERE316	ERE321
	76276-13	76276-18	Large	10 - 12	ERE317	ERE322
	76276-14	76276-19	X Large	> 12	ERE318	ERE323

Code	Description	Size
76277-94	Heel wedge	Universal