

Actimove[®]

Feel good. Move better.

Walker Range

Focused Protection and Stabilisation
for the Lower Leg, Ankle and Foot



High quality support for every step

Actimove® Standard Walker,
Actimove® Standard Air Walker and
Actimove® Closed Shell Air Walker

The Actimove Walker range ensures protection of the foot, ankle and lower leg. Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort. Walking with Actimove Walker products feels natural due to the rocker bottom which simulates a physiological gait.

Standard Walker

Secure and lightweight stabilisation for the lower leg, ankle and foot



Standard Air Walker

Adjustable stabilisation with integrated pump and valve system to easily alter desired compression level



Closed Shell Air Walker

Designed to offer the highest protection for the lower leg, ankle and foot

Indications:

- Fractures of the foot and ankle
- Lower extremity soft-tissue injuries
- Severe ankle and foot sprains
- Post-operative and rehabilitation use

For Closed Shell Air Walker only:

- Achilles tendon ruptures (only in combination with separately available Actimove Heel Wedge for Walkers)

Actimove Walker range advanced features for safety and patient care



Circumferential strapping

For critical stability



Shock absorbing sole and rocker bottom

Minimises impact and simulates a natural gait



Comfortable breathable liner

Integrated pump system of Actimove Standard Air Walker and Actimove Closed Shell Air Walker



Integrated pump and valve system

Designed to easily adjust desired compression level



Additional adjustable features of the Actimove Closed Shell Air Walker



Heel wedge

Adjustable heel wedge for Achilles tendon injuries



Adjustable toe guard

Provides protection in the toe area





Actimove Walker

Essity Product Code					NHS SC Code	
	Low Walker	High Walker	Size	Measurements UK Shoe Size	Low Walker	High Walker
	76270-30	76270-35	X Small	< 3.5	ERE270	ERE275
	76270-31	76270-36	Small	4 - 6.5	ERE271	ERE276
	76270-32	76270-37	Medium	7 - 9.5	ERE272	ERE277
	76270-33	76270-38	Large	10 - 12	ERE273	ERE278
	76270-34	76270-39	X Large	> 12	ERE274	ERE279

Actimove Walker Air

Essity Product Code					NHS SC Code	
	Low Walker	High Walker	Size	Measurements UK Shoe Size	Low Walker	High Walker
	76272-30	76272-35	X Small	< 3.5	ERE328	ERE334
	76272-31	76272-36	Small	4 - 6.5	ERE333	ERE335
	76272-32	76272-37	Medium	7 - 9.5	ERE326	ERE330
	76272-33	76272-38	Large	10 - 12	ERE332	ERE329
	76272-34	76272-39	X Large	> 12	ERE331	ERE327

Actimove Walker Air Closed Shell

Essity Product Code					NHS SC Code	
	Low Walker	High Walker	Size	Measurements UK Shoe Size	Low Walker	High Walker
	76276-10	76276-15	X Small	< 3.5	ERE314	ERE319
	76276-11	76276-16	Small	4 - 6.5	ERE315	ERE320
	76276-12	76276-17	Medium	7 - 9.5	ERE316	ERE321
	76276-13	76276-18	Large	10 - 12	ERE317	ERE322
	76276-14	76276-19	X Large	> 12	ERE318	ERE323
	Code	Description	Size			
	76277-94	Heel wedge	Universal			