



Personalized PDF Catalog
Catalogue Generated 23 January 2025

Disposal bin, SHARPSGUARD®, orange

Supplier: DANIELS HEALTHCARE LIMITED

PP.





SHARPSGUARD® orange is a comprehensive range of containers for the disposal of sharps, excluding those contaminated with medicinal products and their residues. These general purpose sharps containers provide maximum capacity in minimum working space. A robust side handle prevents the container from rolling for safe transport and ease of disposal.

For the SHARPSGUARD® orange 22 XA bin (129-0721) the 'XA' denotes extra access. This container is designed for the disposal of larger instruments. It has a hands-free hinged door, which can stay open at 45 and 90 degrees, providing a large aperture for the disposal of large instruments with minimum handling.

The SHARPSGUARD® yellow 35L theatre bin (129-1046) is designed to be used upright to ease the disposal of long instruments such as staple guns as well as chest drains. It is ideal for theatres, acute medical areas and X-ray departments but is not recommended for high public access areas.

Certifications: They are manufactured to the most stringent quality and safety standards: BS EN ISO 23907, UN3291, AFNOR NFX 30-500 and ISO 9002.

ORDER

VWR Catalog Number	Price	Quantity	Description	Colour	Capacity	W×D×H	Ø×H	P
129-0450	£468.00		SHARPSGUARD® orange 5	Orange/yellow	5 l	- mm	176×285 mm	4
129-0726	£180.00		SHARPSGUARD® orange 11.5	Orange/yellow	11,5 l	- mm	280×250 mm	2
129-0725	£128.00		SHARPSGUARD® orange 22	Orange/yellow	22 l	- mm	330×345 mm	1
129-0721	£117.00		SHARPSGUARD® orange 22 xa	Orange/yellow	22 l	- mm	350×345 mm	7

Product(s) marked with this symbol are discontinued - sold till end of stock. Alternatives may be available by searching with the VWR Catalogue Number listed above. Contact VWR Customer Service on 0800 22 33 44.