


Technical Specification

AD-VANCE VISION SINGLE USE FACE VISOR

<p style="text-align: center;">Part Number AD-V 001</p> <p>Part number and Date of manufacture are moulded into headband for traceability</p> <p style="text-align: center;">CE marked</p> <p style="text-align: center;">PPE category 1</p> <p style="text-align: center;">Applications</p> <p>Protection from projections of particulate matter in good visibility conditions</p>	
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------

Components	Main Headband (adjustable)	Food Grade Polypropylene with REACH certification
	PET shield	0.5mm PET sheet Anti-fog Visor suitable for use with prescription lenses and protective masks Optically correct-can be used for prolonged periods as it does not create distortion or fatigue
Cleaning/ sterilisation instructions	<p>Where required, the PP headband can be cleaned using boiling water, and both parts can be cleaned effectively with Isopropyl alcohol of similar appropriate disinfectant.</p> <p>If the item will be re-used it should only be by the same person.</p>	
Size adjustment range	<p>45.9cm to 68.8 cm head circumference</p> <p>Mask is offset from face by 50mm to allow space for glasses wearers and prevent fogging</p>	
Compliance and certification	UK/ EU testing	PPE Class I of REGULATION (EU) 2016/425 CE conformity
Other material safety information	<p>PVC free</p> <p>Phthalate free</p> <p>BPA free</p>	