

iHealth® PT3

Non contact infrared thermometer



iHealth PT3 - Features

- LED display with backlight
- Easy to use
- 1 second measurement
- 1 measurement mode: Forehead
- °C or °F mode
- Automatic switch off



Technical specifications:

- Model: PT3
- Commercial name: iHealth PT3
- Dimensions: 137mm x 34mm x 39mm
- Weight: 90g (without batteries)
- Power: 2x 1.5V AAA batteries
- Measurement distance: ≤ 3 cm
- Measurement duration: 1s
- Accuracy:
 - $\pm 0.2^{\circ}\text{C}/0.4^{\circ}\text{F}$ during $35.5\sim 42.0^{\circ}\text{C}$ ($95.0\sim 107.6^{\circ}\text{F}$)
 - $\pm 0.3^{\circ}\text{C}/0.5^{\circ}\text{F}$ for other ranges of temperatures
- Measurement range:
 $32^{\circ}\text{C}\sim 42.9^{\circ}\text{C}$ ($89.6^{\circ}\text{F}\sim 109.2^{\circ}\text{F}$)
- Operating environment:
Temp. : $15^{\circ}\text{C}\sim 40^{\circ}\text{C}$, humidity : $\leq 85\%$ RH
- Conditions for transport and storage: $-20^{\circ}\text{C}\sim 55^{\circ}\text{C}$ ($-4^{\circ}\text{F}\sim 131^{\circ}\text{F}$)



In the box :

- 1 Non contact infrared Thermometer iHealth PT3
- 1 User Manual
- 1 Quick Start Guide
- 2 AAA Batteries



www.iHealthLabs.eu

