



## MATERIAL SAFETY DATA SHEET

### STEROGRIP

PRODUCT CODE: VARIOUS

PRODUCT NAME: STEROGRIP TUBULAR SUPPORT BANDAGE

USE: For support to limbs sprains, dislocations, etc may be re-used. Can be cut to size. For external use only.

#### Physical Description / Properties

Appearance and Odour:	Off white ribbed stockinette, constructed from cotton with rubber cores covered by polyester.
Melting Point/Boiling Point:	Not determined
Vapour Pressure:	N/A
Specific Gravity:	N/A
Flashpoint:	N/A
Solubility in Water (g/L):	N/A
Percent volatile by volume:	N/A

---

#### Components

Cotton

Rubber cores coated with polyester – Contains Latex, but not directly to the skin.

---

#### Health Hazard Information

Health Effects

Acute

Swallowed:

Not a likely route of exposure. If swallowing did occur, it may cause gastrointestinal irritation, nausea and dizziness.

Eye:

Not a likely route of exposure.

Skin

Prolonged skin contact or repeated exposure not likely to cause significant skin irritation unless there was a pre-existing skin disorder. Alternatively irritation may occur if an unsuitable width is used.

Inhaled: Symptoms may include redness and itchiness.  
Chronic: Not a likely route of exposure.  
Not likely.

---

### **First Aid**

Swallowed: In the unlikely event that a bandage piece was swallowed. DO NOT induce vomiting. Give water. Get medical attention.  
Eye: Not likely  
Skin: Wash thoroughly with soap and water  
Inhaled: Not likely  
First-Aid facilities: Eye wash plus normal washroom facilities nearby.  
Advice to Doctor: Treat symptomatically

---

### **Precautions for Use**

Wearing: Not recommended to wear overnight due to circulation problems.  
Incorrect size used: May result in circulation problems.  
Flammability: The bandage is non-flammable. Packaging is a combustible solid.

---

### **Safe Handling Information**

Storage and Transport: Store in a cool (25 C) dry place.  
Spills and Disposal: Dispose of in accordance with local waste regulations. May be incinerated in an approved facility.  
Used: soiled bandages may be required to be disposed of in accordance with the local authority's infectious waste requirements.  
Fire or Explosion hazard: No fire or explosion hazard.  
Packaging material: combustible solid.  
Extinguishing Media: According to the surrounding area. Any standard fire extinguisher may be used.

---

## Other Information

Other Information: Stable at room temperature.  
Conditions to Avoid: Heat.  
Incompatible Materials: Avoid contact with strong oxidising agents.  
Hazardous Decomposition Products: Carbon monoxide.  
Classification: Class I – Medical Device  
Environmental Issues: No environmental issues expected

---

Date of First Issue: April 2008  
Reviewed: October 2016

Prepared by: Quality Manager

This Data Sheet summarises our best knowledge of the health and safety hazard information of the product and how to safely handle and use the product in the workplace. Each user should read this data sheet and consider the information in the context of how the product will be handled and used in the workplace including in conjunction with other products. If clarification or further information is needed to ensure that an appropriate risk assessment can be made, the user should contact this company.

Steroplast Healthcare Ltd  
Alpha Point, Bradnor Road  
Manchester, M22 4TE  
Tel: 0161 902 3030  
Fax: 0161 902 3040  
Email: [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk)

