

iHealth Track

Connected Arm Blood Pressure Monitor

Compatible with iPhone, iPad and Android smartphones



Blood pressure monitoring on iPhone, iPad and Android smartphones

For iOS and Android version 4.4 +.

See the list of compatible devices on

ihealthlabs.eu/support

iHealth Track features

- Clinically validated blood pressure (ESH)
- Systolic, diastolic and heart rate (number of beats per minute)
- Detection of arrhythmia
- Universal cuff (22 - 42cm)
- Evaluation of measurements as recommended by the WHO*
- Colour changing display
- Easy to use and comfortable
- 99 memories
- One-Click : Save all your measurements in your MyVitals app account

iHealth MyVitals App

- Track data, graphs and history on the iHealth MyVitals App
- Add measurement reminders
- Free and secure data storage
- Share results with your coach and friends or healthcare professionals

iHealth



*World Health Organisation

Compatible with many other great iHealth products:
Activity trackers, scales and SpO2 monitors



Technical specifications

- Connectivity: Bluetooth 4.0
- Dimensions: 98 mm x 98 mm x 46mm
- Cuff circumference: 22-42cm
- Weight: 348g (Cuff+batteries included)
- Measuring method: Oscillometric, with automatic inflation and deflation
- Power : 4x1.5V SIZE AAA
- AAA battery life : 250 measurements
- Range: Sys: 60-260mmHg / Dia: 40-199mmHg
- Pulse: 40-180bts/min,
- Operating Environment:
Temp: 5 °C ~ 40 °C, humidity: <85% RH
- Conditions for transport and storage:
Temp: -20 °C ~ 55 °C, humidity <85% RH



In the box:

- 1 iHealth BP Track blood pressure monitor
- 1 Cuff
- 1 user manual
- 4 AAA Batteries



www.iHealthLabs.eu

iHealth™
Connected to my health