



The exclusive Irish distributor of


COOLMED<sup>®</sup>

When health  
matters,  
we're there



# Contents

- 4 **CMS29** Solid Door Small Tabletop Medical Refrigerator
- 5 **CMG29** Glass Door Small Tabletop Medical Refrigerator
- 6 **CMS59** Solid Door Small Tabletop /Under Counter Medical Refrigerator
- 7 **CMG59** Glass Door Small Tabletop /Under Counter Medical Refrigerator
- 8 **CMS125** Solid Door Medium Under Counter Medical Refrigerator
- 9 **CMG125** Glass Door Medium Under Counter Medical Refrigerator
- 10 **CMS300** Solid Door Large Medical Refrigerator
- 11 **CMG300** Glass Door Large Medical Refrigerator
- 12 **CMS400** Solid Door Large Medical Refrigerator
- 13 **CMG400** Glass Door Large Medical Refrigerator





# COOLMED<sup>®</sup>

Say hello to award-winning refrigeration, that doesn't cost the earth.

Born from a mission to bring reliable, and high-quality refrigeration, CoolMed is the UK's fastest growing medical refrigeration brand. We're product serve 80% of the UK pharmacy market, including TESCO and Wells, with high-quality fridges at unbeatable prices.

Now, we're delighted to the bring the same quality products to Ireland via our partnership via our distribution partner, Promed. Their commitment to exceptional customer service means our customers can expect the same high-quality medical fridges, backed by a 5 year warranty and unbeatable technical support.

From local practices to major healthcare providers, trust CoolMed and Promed to supply award-winning reliability.

**Because when health matters, we're there.**

# CMS29 | Solid Door Small Tabletop Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally friendly refrigerant
- Rear brackets for wall mounting if required



Climate Class	4
Gross Capacity	30L
Net Capacity	29L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	2+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	445*485*535
Depth with door open (mm)	875
Number of fans	1
Energy Consumption (kW/year)	219
Noise Level (dbA)	50

# CMG29 | Glass Door Small Tabletop Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally friendly refrigerant
- Rear brackets for wall mounting if required



Climate Class	4
Gross Capacity	30L
Net Capacity	29L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	2+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	445*485*535
Depth with door open (mm)	875
Number of fans	1
Energy Consumption (kW/year)	219
Noise Level (dbA)	50

# CMS59 | Solid Door Small Tabletop/Under Counter Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant
- Rear brackets for wall mounting if required



Climate Class	4
Gross Capacity	62L
Net Capacity	59L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	3+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	445*485*535
Depth with door open (mm)	875
Number of fans	1
Energy Consumption (kW/year)	219
Noise Level (dbA)	50

# CMG59 | Glass Door Small Tabletop/Under Counter Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant
- Rear brackets for wall mounting if required



Climate Class	4
Gross Capacity	62L
Net Capacity	59L
Coolant	R600A
Temperature Range (°C)	+2-8 °C
No. of Shelves	3+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	445*485*535
Depth with door open (mm)	875
Number of fans	1
Energy Consumption (kW/year)	219
Noise Level (dbA)	50

# CMS125 | Solid Door Medium Under Counter Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant

Climate Class	4
Gross Capacity	148L
Net Capacity	145L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	3+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	595*595*820
Depth with door open (mm)	1145
Number of fans	1
Energy Consumption (kW/year)	237.25
Noise Level (dbA)	50



# CMG125 | Glass Door Medium Under Counter Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant

Climate Class	4
Gross Capacity	148L
Net Capacity	145L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	3+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	595*595*820
Depth with door open (mm)	1145
Number of fans	1
Energy Consumption (kW/year)	237.25
Noise Level (dbA)	50



# CMS300 | Solid Door Large Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant

Climate Class	4
Gross Capacity	308L
Net Capacity	300L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	6+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	595*695*1845
Depth with door open (mm)	1235
Number of fans	1
Energy Consumption (kW/year)	547.50
Noise Level (dbA)	50



# CMG300 | Glass Door Large Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant

Climate Class	4
Gross Capacity	308L
Net Capacity	300L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	6+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	595*695*1845
Depth with door open (mm)	1235
Number of fans	1
Energy Consumption (kW/year)	547.50
Noise Level (dbA)	50



# CMS400 | Solid Door Large Medical Refrigerator

## Key Features:

- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant

Climate Class	4
Gross Capacity	408L
Net Capacity	400L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	6+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	595*695*1845
Depth with door open (mm)	1190
Number of fans	1
Energy Consumption (kW/year)	657
Noise Level (dbA)	50



# CMG400 | Glass Door Large Medical Refrigerator

## Key Features:


- Core temperature display represented by air temperature probe
- External digital controller and temperature display with battery back up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock with two keys
- R600a environmentally-friendly refrigerant

Climate Class	4
Gross Capacity	408L
Net Capacity	400L
Coolant	R600A
Set Temperature Range (°C)	+2-8 °C
No. of Shelves	6+1 (bottom shelf)
Unit Dimensions WxDxH (mm)	595*695*1845
Depth with door open (mm)	1190
Number of fans	1
Energy Consumption (kW/year)	657
Noise Level (dbA)	50



When health matters,  
we're are there.

Sales enquiries:

 Promed,  
Tulligmore  
Killorglin, Co Kerry  
Ireland  
V93 YD9D

 [info@promed.ie](mailto:info@promed.ie)

 1800 619 619

General enquiries:

 CoolMed Ltd,  
8 Coronet Way  
Salford  
Greater Manchester  
M50 1RE

 [sales@coolmed.co.uk](mailto:sales@coolmed.co.uk)

 0161 7725 666

