

MEDI-GRIP® Latex Standard

SURGICAL GLOVES

Powder-free latex surgical gloves, offering reliable comfort and protection

- **Dependable comfort and dexterity:** The soft latex formulation and optimized thickness of these latex surgical gloves ensure surgical staff get the hand protection, comfort and dexterity they need when performing everyday procedures
- **Ergonomic performance:** Ansell's proprietary anatomic-shaped formers are designed to offer the wearer excellent fit and feel with minimal hand fatigue
- **Effortless donning:** MEDI-GRIP® Latex Standard powder-free surgical gloves have a polymer coating applied to the inside of the glove to make them easy to don
- **Safe chemo handling:** Tested in accordance with ASTM D6978, these reliable surgical gloves can be safely used with a range of chemotherapy drugs**

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

**Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Sterile procedures in the O.R.
- Enhanced grip in wet environments
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES



MEDI-GRIP®
Everyday Protection you can
Depend on

Trusted around the world to deliver the every day protection surgeons, O.R. nurses and surgical technicians can depend on to perform their daily tasks safely and comfortably.

MEDI-GRIP® Latex Standard

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Beaded
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330200755	330200760	330200765	330200770	330200775	330200780	330200785	330200790

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com