



Ethicon Wound Closure Catalog

Sutures - Adhesives

Johnson&Johnson
MedTech

Table of Contents

General Information 4

Needle Index11

True-to-Scale Needle Diagram 29

DERMABOND™ Portfolio 34

STRATAFIX™ Knotless Tissue Control Devices 36

Sutures: Absorbable 56

Sutures: Nonabsorbable 134

Primary Suture Storage Carts, Carrying Racks, and
Other Wound Closure Devices 208

Suture Chart214

Index 217

General Information

Explanation of a Typical Selection and Page

The needle type and size/shape are shown only for needle/suture combinations. A double needle image indicates a "double-armed" suture, that is, a suture with a needle attached to each end.

ATRALOC™ Surgical Needles

ATRALOC Needles are supplied in a variety of sizes, shapes and chord lengths. Surgeons should select the type of needle that, in their experience, is appropriate for the specific procedure and tissue.

Needle shapes are generally classified according to the degree of curvature of the body — 5/8, 1/2, 3/8, or 1/4 circle and straight — with taper, cutting, blunt points, or TAPERCUT™ Surgical Needle tips.

In general, the same size of needle may be made from finer gauge wire for use in soft or delicate tissues and from heavier gauge wire for use in tough or fibrosed tissues (the surgeon's choice).

Needle Type/Needle Shape Symbol

Needle type and needle shape are shown in this subhead.

Material

Suture materials are designated by color and construction characteristics (eg, Undyed Braided, Violet Braided, etc.).

Length

Suture length is shown in both inches and the approximate metric equivalent. Multi-strand (MS) means there is more than 1 strand of suture per package. For example: MS-2 means 2 strands per pack, MS-4 means 4 strands per pack..

The **CR** refers to CONTROL RELEASE™ Needle Sutures.

Product Name

Sutures and other Ethicon products are designated by trademarked and generic names.

Size

Suture size, denoting the diameter of the material, is in accordance with the United States Pharmacopeia (USP).

Needle Silhouette

Needle silhouettes are shown in ACTUAL SIZE. The designated number above the needle, within the Needle section, is the actual length measured in millimeters.

Code

Product codes ending in suffix N are packaged six eaches, or one-half dozen packets, to a box. Product codes ending in suffix G are packaged one dozen packets to a box. Product codes ending in suffix T are packaged two dozen packets to a box. Product codes ending in suffix H are packaged three dozen to a box. In the few instances where no alphabetic suffix is in the product code, the number of packets per box is identified by a footnote. The number of strands per multi-strand packet appears on the left of suture length. Primary suture packets do not have suffix G, H, or T.

Needle Length

Needle length is the measurement in millimeters of the length of the needle along the curvature. Ethicon has applied a rounding convention to all needle lengths. Needles less than 10 mm are rounded to the nearest tenth. Needles greater than 10 mm are rounded to the nearest whole number.

Explanation of a Typical Selection and Page

MULTIPASS™ Needles


Elements of Needle Design

Needle design involves analyzing a surgical procedure and the density of the tissue involved in great detail. Ethicon engineers work continuously to improve upon their needle line.


The anatomy of the ideal surgical needle includes three factors: Alloy, Geometry, and Coating.

The combination of coating and alloy make up the MULTIPASS Needle Technology, ensuring excellent strength, tissue penetration, and control.¹ MULTIPASS Needles are the highest performing needles that Ethicon offers.¹

CONTROL RELEASE Needle Sutures






Needle is held securely in the holder. Suture is grasped securely just below needle, pulling strand taut.



Needle is released with a straight tug of the needleholder.

Selection and Use of Needleholders¹



- The needleholder must be an appropriate size instrument for the size needle selected
- It must be made from good quality steel with a secure jaw design
- Needles should be grasped in an area about one-third to one-half the distance from the swaged area to the point. Avoid placement on or near the swaged area
- The needle should be placed securely in the tip of the needleholder jaws
- Do not damage taper points or cutting edges when using the needleholder to pull the needle out through the tissue. Grasp as far back as possible
- When placing the needle in tissue, any force applied should be in the direction following the curve of the needle
- Do not take excessively large bites of tissue with small needles
- Do not force or twist the needle in an effort to bring the point out through the tissue. Withdraw the needle and replace tissue
- Do not force a dull needle through tissue; obtain a new one
- Avoid using the needle to bridge or approximate tissues for suturing
- If the needle is held too tightly in a sharp, hard-jawed, or defective needleholder, the needle may be damaged or notched in such a manner that it will have more of a tendency to bend or break on successive passes through tissue
- When the needle holder with needle is passed to the surgeon, the needle should be pointing in the right direction in which it will start to be used without need for readjustments
- In some patients, the tissues may be tougher or fibrosed more than normal and require the use of a heavier gauge needle
- In deep or confined areas, ideal positioning of the needle may not be possible. The surgeon should proceed with caution or use a heavier gauge needle


Note: The American College of Surgeons recommends universal adoption of blunt-tip suture needles for the closure of fascia and muscle in order to reduce needlestick injuries.²

An ATRALOC™ Surgical Needle is an instrument that complements the surgeon's skills. The variety of Ethicon needles offers the surgeon a choice of needles that are dependable for consistent quality for practically all surgical applications.


Steps to Create Performance



Conventional Needle



Super Smooth



Needle Swaging

Controlled Manufacturing

Ethicon makes its own needle/suture combinations. Each step, from raw material to finished needle, is carefully controlled.

Strong, Sharp Micrograin Stainless Steel and Tungsten-Rhenium Alloy

Ethicon makes its needles from premium grade stainless steel and tungsten-rhenium alloy formulated specifically for surgical use.

The Ethicon process controls the orientation of the micrograins so that ductility and strength are closely controlled — for strong needles that start sharp and stay sharp, thus assuring easy passage through tissue.

Super Smooth

Virtually all ATRALOC™ Surgical Needles have a micro-thin super-smooth finish that allows the needle to penetrate and pass through the toughest tissue with an almost fluid-like flow, with minimal resistance and drag, giving the surgeon a feeling of control.


Needle Swaging

Needles are permanently attached to the suture materials, eliminating the need for threading and providing a new, sharp needle each time. ATRALOC Needle swaging gives the surgeon a needle with an unbroken surface from end to end; without double strands, it draws through tissue easier with less tissue trauma.

Ribbed for Security and Control

Longitudinal ribbing on the inside and outside curvatures of most needles has a cross-locking action in the needleholder, which minimizes turning, rocking, or twisting, and provides unsurpassed needle control.

Shapes



Types



Conventional Cutting Needle

Two opposing cutting edges, with a third on inside curve. Change in cross-section from a triangle cutting tip to a flattened body.



Conventional Spatula Needle

Visibility of the point at bottom provides control of depth penetration.



ETHIGUARD™ Blunt Point Needle

Taper body. For blunt dissection and suturing friable tissue.



MICRO-POINT™ Reverse Cutting Needle

Cutting edge on outer curve. Extremely smooth. Extremely sharp for ophthalmic surgery.



MICRO-POINT™ Surgical Needle (Spatula)

Thin, flat profile. Specially designed for ophthalmic anterior segment surgery.



Precision Cosmetic - Conventional Cutting PRIME™ Needle

For delicate plastic or cosmetic surgery. Conventional cutting tip and PRIME geometry for increased sharpness in the tip leading to a flattened body.



Precision Point - Reverse Cutting PRIME Needle

For plastic or cosmetic surgery, reverse cutting tip and PRIME geometry for increased sharpness in the tip.



Reverse Cutting Needle

Cutting edge on outer curve. For tough, difficult-to-penetrate tissues..



SABRELOC™ Spatula Needle

Side-cutting spatula shaped edges. For layers of scleral or corneal tissue. Spatulated needle point is centered for maximum needle stability in thin sclera. Four equidistant and properly defined edges provide greater control.



TAPERCUT™ Surgical Needle

Cutting tip, taper body. For tough tissue, like two needles in one.



Taperpoint Needle

For soft, easily penetrated tissues.



ULTIMA™ Ophthalmic Needle (Spatula)

Reduced edge-angles provide better penetration. Readily facilitates knot rotation in ophthalmic surgery.



VISI-BLACK™ Surgical Needle

Slim taper point needles with a black finish for improved visibility.

Trademarks

- ATRALOC™ Surgical Needle
- CONTROL RELEASE™ Needle or Needle Suture **CR**
- DERMABOND ADVANCED™ Topical Skin Adhesive
- DERMABOND™ Mini Topical Skin Adhesive
- DERMABOND™ PRINEO™ Skin Closure System
- ETHALLOY™ Needle Alloy
- ETHIBOND™ Extra Polyester Suture
- ETHIBOND EXCEL™ Polyester Suture
- ETHIGUARD™ Blunt Point Needle
- ETHILON™ Nylon Suture
- ETHI-PACK™ Pre-cut Sutures
- EVERPOINT™ Cardiovascular Needle **EVP**
- HEMO-SEAL™ Needle Suture **HS**
- LAPRA-TY™ Suture Clip
- LIGAPAK™ Dispensing Reel
- LIGAPAK™ Ligature
- MERSILENE™ Polyester Fiber Suture
- MERSILENE™ Polyester Fiber Strip
- MICRO-POINT™ Surgical Needle
- MONOCRYL™ (poliglecaprone 25) Suture
- MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture
- MULTIPLASS™ Needles 
- NUROLON™ Nylon Suture
- P, PS PRIME™ Needles
- PC PRIME™ Needle
- PDS™ (polydioxanone) Suture
- PDS™ II (polydioxanone) Suture
- PDS™ Plus Antibacterial (polydioxanone) Suture
- PERMA-HAND™ Silk Suture
- PROLENE™ Polypropylene Suture
- PRONOVA™ Poly (Hexafluoropropylene – VDF) Suture
- RELAY Suture Delivery System
- RELAY Suture Package
- SABRELOC™ Spatula Needle
- STRATAFIX™ Knotless Tissue Control Devices
- STRATAFIX™ Spiral Knotless Tissue Control Device
- STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device
- STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device
- STRATAFIX™ Spiral Polypropylene Knotless Tissue Control Device
- STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device
- SUTUPAK™ Pre-cut Sterile Sutures
- TAPERCUT™ Surgical Needle
- ULTIMA™ Ophthalmic Needles
- VICRYL™ (polyglactin 910) Knitted Mesh
- VICRYL™ (polyglactin 910) Woven Mesh
- VICRYL™ (polyglactin 910) Mesh Bag
- VICRYL™ (polyglactin 910) Perirenal Mesh Bag
- VICRYL™ (polyglactin 910) Perisplenic Mesh Bag
- VICRYL™ (polyglactin 910) Perihepatic Mesh Bag
- VICRYL™ (polyglactin 910) Perihepatic Right Lobe Mesh Bag
- VICRYL™ (polyglactin 910) Perihepatic Left Lobe Mesh Bag
- Coated VICRYL™ (polyglactin 910) Suture
- Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture **VB**
- VICRYL RAPIDE™ (polyglactin 910) Suture
- VISI-BLACK™ Surgical Needle

Needle Index

1/4 Circle

TAPERCUT™ Surgical Needle



Needles	Needle Length	Available Products
V75-3	3.5 mm	ETHILON
V100-3	3.5 mm	ETHILON
CIF-4	13 mm	PROLENE

SABRELOC™ Spatula Needle



Needles	Needle Length	Available Products
S-4	8 mm	GUT CHROMIC, MERSILENE
S-29	7.6 mm	COATED VICRYL
S-24	8 mm	MERSILENE
S-14	8 mm	COATED VICRYL, GUT CHROMIC, GUT PLAIN, ETHILON, MERSILENE, PERMA HAND SILK

1/4 Circle

ULTIMA™ Ophthalmic Needle



Needles	Needle Length	Available Products
CS90-6	7 mm	VICRYL
CTC-6L	13 mm	PROLENE
CTC-6	13 mm	PROLENE

3/8 Circle

Taperpoint



Needles	Needle Length	Available Products
TE	32 mm	TEMPORARY PACING WIRE
BB-1	22 mm	TEMPORARY PACING WIRE
BB	17 mm	STRATAFIX SPIRAL MONOCRYL PLUS, GUT CHROMIC, ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK, TEMPORARY PACING WIRE
C-1	13 mm	PDS PLUS, GUT CHROMIC, PDS II, ETHIBOND, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK
BV	11 mm	PDS PLUS, PDS II, PROLENE, PRONOVA
BV-1	9.3 mm	PDS PLUS, PDS II, ETHIBOND, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK
BV175-8	9.3 mm	PROLENE, PRONOVA
BV175-7	8 mm	PROLENE, PRONOVA

3/8 Circle

Taperpoint



Needles	Needle Length	Available Products
BV175-6	8 mm	PROLENE, PRONOVA
BV130-5	6.5 mm	COATED VICRYL, ETHILON, PROLENE, PRONOVA
BV130-4	5 mm	COATED VICRYL, ETHILON
BV130-3	3.8 mm	ETHILON
BV100-4	5.1 mm	ETHILON, PROLENE
BV100-3	3.8 mm	VICRYL
BV75-4	5.1 mm	ETHILON
BV75-3	3.8 mm	ETHILON, PROLENE
BV50-3	3.8 mm	ETHILON

TAPERCUT™ Surgical Needle



Needles	Needle Length	Available Products
V-4	17 mm	ETHIBOND, PROLENE, PERMA HAND SILK
CC-1	13 mm	ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK
CC	9.3 mm	PROLENE
V130-5	5.8 mm	COATED VICRYL
V130-4	4.2 mm	COATED VICRYL
V130-3	3.8 mm	ETHILON

3/8 Circle

Reverse Cutting



Needles	Needle Length	Available Products
LR	75 mm	ETHIBOND, ETHILON, MERSILENE, PROLENE, STEEL, PERMA HAND SILK
FSLX	36 mm	ETHILON, PROLENE, STEEL, PERMA HAND SILK
FSL	30 mm	COATED VICRYL PLUS, COATED VICRYL, ETHILON, NUROLON, PROLENE, PERMA HAND SILK
FS	26 mm	STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, COATED VICRYL, ETHIBOND, ETHILON, MERSILENE, NUROLON, PROLENE, PERMA HAND SILK
FS-1	24 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, ETHILON, MERSILENE, PROLENE, STEEL, PERMA HAND SILK
FS-2	19 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, ETHILON, MERSILENE, PROLENE, STEEL, PERMA HAND SILK
M-2	13 mm	GUT CHROMIC
C-3	13 mm	COATED VICRYL, GUT CHROMIC, ETHILON, PERMA HAND SILK
C-2	13 mm	ETHILON, PERMA HAND SILK

3/8 Circle

Conventional Cutting



Needles	Needle Length	Available Products
PC-5	19 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, ETHILON, PROLENE, PERMA HAND SILK
PC-3	16 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, MONOCRYL, PDS II, ETHILON, PROLENE, PERMA HAND SILK
PC-1	13 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT PLAIN, MONOCRYL, PDS II, ETHILON, NUROLON, PROLENE

Precision Point Reverse Cutting



Needles	Needle Length	Available Products
PSLX	36 mm	ETHILON, PERMA HAND SILK
PSL	30 mm	STRATAFIX SPIRAL MONOCRYL PLUS, COATED VICRYL PLUS, COATED VICRYL, ETHILON, PERMA HAND SILK
PS	26 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, GUT CHROMIC, ETHILON, PERMA HAND SILK
PS-1	24 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT PLAIN, MONOCRYL, PDS II, ETHILON, PROLENE, PERMA HAND SILK
PS-2	19 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX, SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK

3/8 Circle

Precision Point Reverse Cutting



Needles	Needle Length	Available Products
PS-3	16 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK
P-3	13 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, MERSILENE, PROLENE, PERMA HAND SILK
P-1	11 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS, PLUS, COATED VICRYL, VICRYL RAPIDE, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK
P-6	8 mm	ETHILON, PROLENE

ULTIMA™ Ophthalmic Needle



Needles	Needle Length	Available Products
CS140-8	6.5 mm	ETHILON
CS140-6	6.5 mm	VICRYL, ETHILON, MERSILENE, PROLENE


MICRO-POINT™ Surgical Needle (Spatula)



Needles	Needle Length	Available Products
TG140-8	6.5 mm	COATED VICRYL, VICRYL, GUT CHROMIC, GUT PLAIN, ETHILON, PROLENE, PERMA HAND SILK
TG140-6	6.5 mm	COATED VICRYL, ETHILON
TG100-8	6.5 mm	PROLENE
VAS100-4	5 mm	ETHILON

1/2 Circle

Blunt Point



Needles	Needle Length	Available Products
BP	80 mm	SURGICAL GUT SUTURE - CHROMIC MERSILENE SUTURE
BP-1	65 mm	MERSILENE SUTURE
BP-2	48 mm	SURGICAL GUT SUTURE - CHROMIC
TPB-1	65 mm	COATED VICRYL PLUS SUTURE PDS PLUS SUTURE
CTXB	48 mm	STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, PROLENE
CTB	40 mm	COATED VICRYL PLUS, PDS PLUS
CTB-1	36 mm	STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, PROLENE
CTB-2	26 mm	COATED VICRYL PLUS, PDS PLUS
SHB	26 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS COATED VICRYL, MONOCRYL

1/2 Circle

Taperpoint

**Needles****Needle Length****Available Products**

XLH

70 mm

COATED VICRYL PLUS, PDS PLUS,
COATED VICRYL, GUT CHROMIC,
GUT PLAIN, PDS II, ETHIBOND,
PROLENE



LH

40 mm

GUT CHROMIC



MH

36 mm

COATED VICRYL PLUS, COATED
VICRYL, GUT CHROMIC,
MONOCRYL, ETHIBOND, PROLENE,
PRONOVA, PERMA HAND SILK



SH

26 mm

STRATAFIX SPIRAL MONOCRYL
PLUS, STRATAFIX SPIRAL PDS PLUS,
STRATAFIX SYMMETRIC PDS PLUS,
COATED VICRYL PLUS, MONOCRYL
PLUS, PDS PLUS, COATED VICRYL,
VICRYL RAPIDE, GUT CHROMIC,
GUT PLAIN, MONOCRYL, PDS II,
ETHIBOND, MERSILENE, NUROLON,
PROLENE, PROLENE HS, PRONOVA,
PERMA HAND SILK, TEMPORARY
PACING WIRE



SH-1

22 mm

STRATAFIX SPIRAL POLYPRO-
PYLENE, COATED VICRYL PLUS,
PDS PLUS, COATED VICRYL, GUT
CHROMIC, MONOCRYL, PDS II,
ETHIBOND, NUROLON, PROLENE,
PRONOVA, PERMA HAND SILK



SH-2

20 mm

ETHIBOND, PROLENE HS

1/2 Circle

Taperpoint



Needles	Needle Length	Available Products
RB-1	17 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATA- FIX SPIRAL, PDS PLUS, COATED VICRYL, PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBO- ND, NUROLON, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK, TEMPO- RARY PACING WIRE
RB-2	13 mm	COATED VICRYL PLUS, PDS PLUS, PDS II, PROLENE, PRO- LENE HS, PROVONA, PERMA HAND SILK
UCL	13 mm	PDS PLUS, PDS II, ETHIBOND
TF	13 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, ETHIBO- ND, NUROLON, PROLENE, PERMA HAND SILK
RB-3	11 mm	PROLENE
BVH100-3	5 mm	ETHILON

1/2 Circle

TAPERCUT™ Surgical Needle



Needles	Needle Length	Available Products
V-56	60 mm	STEEL
V-40	48 mm	ETHIBOND, STEEL
V-37	40 mm	ETHIBOND, STEEL
V-34	36 mm	STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, ETHIBOND, PROLENE
V-7	26 mm	STRATAFIX SPIRAL PDS PLUS, GUT PLAIN, ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK
V-5	17 mm	PDS PLUS, GUT CHROMIC, PDS II, ETHIBOND, PROLENE, PERMA HAND SILK

1/2 Circle

Reverse Cutting



Needles	Needle Length	Available Products
LS-1	45 mm	STEEL
CPX	48 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, PDS II
CP	40 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, PDS II, NUROLON
CP-1	36 mm	STRATAFIX SPIRAL MONOCRYL PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, STEEL, PERMA HAND SILK
CP-2	26 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, PDS II, NUROLON
OS-8	40 mm	STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, COATED VICRYL, ETHIBOND, NUROLON
OS-6	36 mm	STRATAFIX SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, PDS II, ETHIBOND
OS-4	22 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, PDS II, ETHIBOND
OS-2	17 mm	ETHIBOND
X-1	22 mm	COATED VICRYL PLUS, COATED VICRYL, GUT PLAIN, ETHILON, NUROLON, PERMA HAND SILK

1/2 Circle

Reverse Cutting



Needles	Needle Length	Available Products
J-1	17 mm	COATED VICRYL, GUT CHROMIC, PERMA HAND SILK
M-1	13 mm	GUT CHROMIC, GUT PLAIN
PS-4C	16 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC
PS-4	16 mm	STRATAFIX SPIRAL MONOCRYL PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK
PS-5	13 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, PROLENE
PS-6	10 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL
P-2	8 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II
G-2	8 mm	GUT CHROMIC
G-7	8 mm	PERMA HAND SILK

Conventional Cutting



Needles	Needle Length	Available Products
PC-12	16 mm	COATED VICRYL, MONOCRYL, ETHILON

SABRELOC™ Spatula Needle



Needles	Needle Length	Available Products
S-28	7.6 mm	COATED VICRYL
S-2	8 mm	GUT CHROMIC, MERSILENE
S-22	8 mm	ETHILON

1/2 Circle

ULTIMA™ Ophthalmic Needle



Needles	Needle Length	Available Products
CS175-8	7 mm	GUT CHROMIC, ETHILON
CS160-8	5.5 mm	VICRYL, ETHILON, PROLENE
CS175-6	7 mm	ETHILON
CS160-6	5.5 mm	VICRYL, ETHILON, MERSILENE, PROLENE
CSM-6	4.3 mm	ETHILON

MICRO-POINT™ Surgical Needle (Spatula)



Needles	Needle Length	Available Products
TG175-8	7 mm	ETHILON
TG160-8	5.5 mm	COATED VICRYL
TG160-6	5.5 mm	COATED VICRYL, ETHILON
TG160-4	5.3 mm	ETHILON, MERSILENE

5/8 Circle

Taperpoint



Needles	Needle Length	Available Products
UR-4	40 mm	COATED VICRYL, GUT CHROMIC
UR-5	36 mm	COATED VICRYL PLUS, MONOCRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL
UR-6	26 mm	STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL
TF-5	8 mm	PERMA HAND SILK

Straight

Taperpoint



Needles	Needle Lengths	Available Products
ST-1	45 mm	PERMA HAND SILK
ST-3	22 mm	COATED VICRYL, ETHIBOND, PDS II
ST-4	19 mm	COATED VICRYL, ETHIBOND, PDS II

ST100-3 3 mm ETHILON

TAPERCUT™ Surgical Needle



Needles	Needle Lengths	Available Products
V-47	36 mm	PROLENE

Straight Cutting



Needles	Needle Lengths	Available Products
SKS-3	60 mm	TEMPORARY PACING WIRE
SKS	60 mm	TEMPORARY PACING WIRE
KS	60 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, PROLENE, STEEL, PERMA HAND SILK
TS	48 mm	PROLENE
SC-1	13 mm	GUT PLAIN

SABRELOC™ Spatula Needle



Needles	Needle Lengths	Available Products
STC-6	16 mm	PROLENE

Trocar

Trocar

Needles	Needle Lengths	Available Products
STP-10	254 mm	PDS II, ETHIBOND

Compound Curve

ULTIMA™ Ophthalmic Needle



Needles	Needle Lengths	Available Products
CSC-6	5 mm	ETHILON, PROLENE

Bi-Curve

ULTIMA™ Ophthalmic Needle



Needles	Needle Lengths	Available Products
CSB-6	5 mm	VICRYL, ETHILON

Sternum

Conventional Cutting - Sternum



Needles	Needle Lengths	Available Products
CCS-1	60 mm	COATED VICRYL, PDS II, STEEL

CCS	48 mm	ETHIBOND, STEEL
-----	-------	-----------------

Center Point Spatula

SABRELOC™ Spatula Needle



Needles	Needle Lengths	Available Products
RD-1	8 mm	COATED VICRYL, ETHILON, MERSILINE

Endoscopy

Taperpoint



Needles	Needle Lengths	Available Products
EN-S	19 mm	PDS
V-38	30 mm	PDS II

True-to-Scale Needle Diagram

True-to-Scale Needle Diagram

Cutting Edge


Precision Point 

Reverse Cutting

3/8 Circle



1/2 Circle




Reverse Cutting 

3/8 Circle



1/2 Circle



Precision Cosmetic 

Conventional Cosmetic PRIME™

3/8 Circle



1/2 Circle



MICRO-POINT™ 

Reverse Cutting

3/8 Circle



1/2 Circle



ULTIMA™ 

Spatula


1/8 Circle



1/4 Circle



3/8 Circle



1/2 Circle



ULTIMA™ 

Compound Curve



SPATULA 

1/4 Circle



3/8 Circle




1/2 Circle



MICRO-POINT™ SPATULA 

3/8 Circle



1/2 Circle



SABERLOC™ 

Conventional Spatula

1/4 Circle



1/2 Circle



Spatula

1/8 Circle



1/4 Circle



1/2 Circle




Center Point Spatula

1/4 Circle



Straight Transchamber

STC-6

STRAIGHT CUTTING 

SC-1

SKS

TS

SKS-3

KS

STERNUM 




Taperpoint




3/8 Circle



5/8 Circle




1/2 Circle



STRAIGHT 




TAPERCUT™




Trocar

STP-10

Blunt Point



Endoscopy



DERMABOND™

Portfolio

DERMABOND™ Portfolio

DERMABOND ADVANCED™ Topical Skin Adhesive



Product Code	Description	Adhesive Volume	EA/BX
DNX12	2-octyl cyanoacrylate (2-OCA) topical skin adhesive with an easy-to-use applicator for precise control for wide and fine line applications	0.7 mL	12
DNX6	2-OCA topical skin adhesive with an easy-to-use applicator for precise control for wide and fine line applications	0.7 mL	6

DERMABOND™ PRINEO™ Skin Closure System



Product Code	Description	Adhesive Volume	EA/BX
CLR222US	2-OCA topical skin adhesive and 22 cm self-adhering mesh creates an optimal wound healing environment ¹⁻³ and leads to greater overall patient satisfaction when compared to staples ^{4*†}	3.8 mL	1 × 22 cm mesh patch
CLR422US	2-OCA topical skin adhesive and 42 cm self-adhering mesh creates an optimal wound healing environment ¹⁻³ and leads to greater overall patient satisfaction when compared to staples ^{4*†}	3.8 mL	2 × 22 cm mesh patches
CLR602US	2-OCA topical skin adhesive and 60 cm self-adhering mesh creates an optimal wound healing environment ¹⁻³ and leads to greater overall patient satisfaction when compared to staples ^{4*†}	3.8 mL	1 × 60 cm mesh dispenser

DERMABOND™ Mini Topical Skin Adhesive



Product Code	Description	Adhesive Volume	EA/BX
DHVM12	2-OCA high viscosity formulation for small incisions and lacerations	0.36 mL	12

*Double-blinded quantitative market research study comparing surgeon experience with DERMABOND PRINEO System and skin staples in total knee arthroplasty. N=83 orthopaedic surgeons. 90% c.i. Fielded June/July 2017.

†In a TKA

References: 1. Su W. Study report for in vitro evaluation of microbial barrier properties of DERMABOND ProTape. 06TR071. December 4, 2006. Ethicon, Inc. 2. Kumar A. Study to Compare the tissue holding strength of PRINEO™ skin closure system with conventional wound closure techniques. AST-2012-0290. October 11, 2012. Ethicon, Inc. 3. Kannon GA, Garrett AB. Moist wound healing with occlusive dressings. *Dermatol Surg*. 1995;21:583-590. 4. PRINEO Claims Research Quant Detail. August 16, 2017. Ethicon, Inc.

STRATAFIX™ Knotless Tissue Control Devices

STRATAFIX™ Knotless Tissue Control Devices

STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device is the only barbed suture appropriate for closing fascia.^{1,2*}

STRATAFIX™ Spiral Knotless Tissue Control Devices provide smooth tissue passage and a secure hold for excellent tissue approximation.^{1,3-9}



*Refers to STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device only.

References: 1. Time Zero Tissue Holding - Competitive Claims Comparisons for STRATAFIX Knotless Tissue Control Devices vs Various Products. 100326296. May 20, 2015. Ethicon, Inc. 2. V-Loc 180 Absorbable Wound Closure Device. Instructions for Use. Medtronic. 3. Moran ME, Marsh C, Perrotti M. Bidirectional-barbed sutured knotless running anastomosis v classic Van Velthoven suturing in a model system. *J Endourol.* 2007;21(10):1175-1178. 4. Vakil JJ, O'Reilly MP, Sutter EG, Mears SC, Belkoff SM, Khanuja HS. Knee arthroscopy repair with a continuous barbed suture: a biomechanical study. *J Arthroplasty.* 2011;26(5):710-713. 5. Eickmann T, Quane E. Total knee arthroplasty closure with barbed sutures. *J Knee Surg.* 2010;23(3):163-167. 6. Einarsson JI, Chavan NR, Suzuki Y, Jonsdottir G, Vellinga TT, Greenberg JA. Use of bidirectional barbed suture in laparoscopic myomectomy: evaluation of perioperative outcomes, safety, and efficacy. *J Minim Invasive Gynecol.* 2011;18(1):92-95. 7. Levine BR, Ting N, Della Valle CJ. Use of a barbed suture in the closure of hip and knee arthroplasty wounds. *Orthopedics.* 2011;34(9):e473-e475. doi:10.3928/01477447-20110714-35 8. Ethicon study AST-2012-0331. Tissue gapping under tension of porcine cadaveric skin incisions closed with Stratafix Spiral in comparison to Monocryl in both interrupted and continuous stitching patterns. Approved on August 24, 2012. Ethicon, Inc. 9. Nawrocki J. Technical memo: Explanation on the science behind the anchors on knotless tissue control devices. 2017. Ethicon, Inc.

STRATAFIX™ Spiral Polypropylene Knotless Tissue Control Device






Taperpoint 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 SH-1 22mm 1/2 Circle	2-0	Undyed	12"/30cm	12	SXPL1B400
 CT-3 22mm 1/2 Circle	0	Undyed	5"/14cm x 5"/14cm	12	SXPL2B401
 CT-1 36mm 1/2 Circle	1	Undyed	9"/23cm x 9"/23cm	12	SXPL2B400

Bidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device








Taperpoint 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	4-0	Undyed	6"/16cm + 6"/16cm	12	SXMP2B403
	3-0	Violet	3"/7cm + 3"/7cm	12	SXMP2B415
	3-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B416
	3-0	Undyed	6"/16cm + 6"/16cm	12	SXMP2B402
 SH-1 22mm 1/2 Circle	3-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B422

Bidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 <p>SH 26mm 1/2 Circle</p>	3-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B421	
	2-0	Violet	3"/7cm + 3"/7cm	12	SXMP2B419	
	2-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B420	
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Undyed	5"/14cm + 5"/14cm	12	SXMP2B400	
	2-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B423	
	2-0	Undyed	14"/36cm + 14"/36cm	12	SXMP2B401	
 <p>UR-6 26mm 5/8 Circle</p>	3-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B417	
	2-0	Violet	6"/16cm + 6"/16cm	12	SXMP2B418	
 <p>ST-4 19mm Straight</p>	3-0	Undyed	12"/30cm x 12"/30cm	12	SXMP2B427	


Bidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device



PRIME Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>PS-2 19mm 3/8 Circle</p>	4-0	Undyed	3"/7cm + 3"/7cm	12	SXMP2B405
	4-0	Undyed	5"/14cm + 5"/14cm	12	SXMP2B407
	4-0	Undyed	12"/30cm + 12"/30cm	12	SXMP2B409
	4-0	Undyed	18"/45cm + 18"/45cm	12	SXMP2B151
	3-0	Undyed	3"/7cm + 3"/7cm	12	SXMP2B404
	3-0	Undyed	5"/14cm + 5"/14cm	12	SXMP2B406
	3-0	Undyed	12"/30cm + 12"/30cm	12	SXMP2B408
	3-0	Undyed	18"/45cm + 18"/45cm	12	SXMP2B152
 <p>PS-1 24mm 3/8 Circle</p>	4-0	Undyed	12"/30cm + 12"/30cm	12	SXMP2B150
	3-0	Undyed	12"/30cm + 12"/30cm	12	SXMP2B410
 <p>PS 26mm 3/8 Circle</p>	4-0	Undyed	12"/30cm + 12"/30cm	12	SXMP2B413
	3-0	Undyed	5"/14cm + 5"/14cm	12	SXMP2B411
	3-0	Undyed	12"/30cm + 12"/30cm	12	SXMP2B412

Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>CP-2 26mm 1/2 Circle</p>	2-0	Undyed	5"/14cm + 5"/14cm	12	SXMP2B414

Unidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device




Taperpoint							
Needle	Suture USP	Color	Length	EA/BX	Product Code		
 RB-1 17mm 1/2 Circle	4-0	Undyed	6"/15cm	12	SXMP1B434		
	4-0	Undyed	12"/30cm	12	SXMP1B435		
	3-0	Undyed	6"/15cm	12	SXMP1B424		
	3-0	Undyed	9"/23cm	12	SXMP1B438		
 CT-2 26mm 1/2 Circle	3-0	Undyed	12"/30cm	12	SXMP1B425		
	3-0	Violet	6"/15cm	12	SXMP1B449		
	3-0	Violet	9"/23cm	12	SXMP1B450		
	2-0	Undyed	6"/15cm	12	SXMP1B415		
	2-0	Undyed	12"/30cm	12	SXMP1B416		
	2-0	Undyed	18"/45cm	12	SXMP1B417		
	2-0	Violet	6"/15cm	12	SXMP1B443		
	2-0	Violet	12"/30cm	12	SXMP1B444		
 SH 26mm 1/2 Circle	2-0	Violet	18"/45cm	2	SXMP1B445		
	4-0	Undyed	27"/70cm	12	SXMP1B437		
	3-0	Undyed	9"/23cm	12	SXMP1B427		
	3-0	Undyed	27"/70cm	12	SXMP1B428		
	3-0	Violet	6"/15cm	12	SXMP1B451		
	3-0	Violet	9"/23cm	12	SXMP1B452		
	2-0	Undyed	6"/15cm	12	SXMP1B408		
	2-0	Undyed	9"/23cm	12	SXMP1B409		
	2-0	Undyed	12"/30cm	12	SXMP1B410		
	2-0	Violet	9"/23cm	2	SXMP1B447		
 CT-1 36mm 1/2 Circle	3-0	Undyed	18"/45cm	12	SXMP1B429		
	2-0	Undyed	6"/15cm	12	SXMP1B411		
	2-0	Undyed	12"/30cm	12	SXMP1B412		
	2-0	Undyed	18"/45cm	12	SXMP1B413		
	2-0	Undyed	27"/70cm	12	SXMP1B414		
	2-0	Violet	6"/15cm	12	SXMP1B440		
	2-0	Violet	9"/23cm	12	SXMP1B441		
	2-0	Violet	12"/30cm	12	SXMP1B442		
	 ST-4 19mm Straight	3-0	Undyed	12"/30cm	12	SXMP1B433	

Unidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device





Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Product Code
 UR-6 26mm 5/8 Circle	2-0	Violet	9"/23cm	12	SXMP1B448
	3-0	Violet	6"/15cm	12	SXMP1B453
	3-0	Violet	9"/23cm	12	SXMP1B454
	3-0	Violet	12"/30cm	12	SXMP1B455

Cutting Edge Reverse




Needle	Suture USP	Color	Length	EA/BX	Product Code
 PSL 30mm 3/8 Circle	2-0	Undyed	27"/70cm	12	SXMP1B421
 CP-2 26mm 1/2 Circle	2-0	Undyed	12"/30cm	12	SXMP1B420
 CP-1 36mm 1/2 Circle	2-0	Undyed	27"/70cm	12	SXMP1B419

PRIME Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Product Code
 PS-2 19mm 3/8 Circle	4-0	Undyed	12"/30cm	12	SXMP1B117
	4-0	Undyed	18"/45cm	12	SXMP1B118
	4-0	Undyed	27"/70cm	12	SXMP1B119
	3-0	Undyed	6"/15cm	12	SXMP1B105
	3-0	Undyed	12"/30cm	12	SXMP1B106
	3-0	Undyed	18"/45cm	12	SXMP1B107
	3-0	Undyed	24"/60cm	12	SXMP1B108
	3-0	Undyed	27"/70cm	12	SXMP1B109
 PS-1 24mm 3/8 Circle	4-0	Undyed	12"/30cm	12	SXMP1B114
	4-0	Undyed	18"/45cm	12	SXMP1B115
	4-0	Undyed	27"/70cm	12	SXMP1B116
	3-0	Undyed	12"/30cm	12	SXMP1B101
	3-0	Undyed	18"/45cm	12	SXMP1B102
	3-0	Undyed	24"/60cm	12	SXMP1B103
	3-0	Undyed	27"/70cm	12	SXMP1B104





Unidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device



PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 PS 26mm 3/8 Circle	3-0	Undyed	12"/30cm	12	SXMP1B110	
	3-0	Undyed	18"/45cm	12	SXMP1B111	
	3-0	Undyed	27"/70cm	12	SXMP1B113	
 PS-4 16mm 1/2 Circle	4-0	Undyed	18"/45cm	12	SXMP1B120	

Bidirectional STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 RB-1 17mm 1/2 Circle	3-0	Violet	6"/16cm + 6"/16cm	12	SXPP2B431	
 SH-1 22mm 1/2 Circle	3-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B430	
 SH 26mm 1/2 Circle	3-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B428	
	2-0	Violet	4"/10cm + 4"/10cm	12	SXPP2B413	
	2-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B414	
	2-0	Violet	14"/36cm + 14"/36cm	12	SXPP2B427	

Bidirectional STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device



Taperpoint


Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>CT-2 26mm 1/2 Circle</p>	0	Violet	5"/14cm + 5"/14cm	12	SXPP2B404
	2-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B407
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Violet	9"/24cm + 9"/24cm	12	SXPP2B408
	2-0	Violet	14"/36cm + 14"/36cm	12	SXPP2B412
	0	Violet	5"/14cm + 5"/14cm	12	SXPP2B409
	0	Violet	9"/24cm + 9"/24cm	12	SXPP2B410
	0	Violet	14"/36cm + 14"/36cm	12	SXPP2B411
	0	Violet	18"/45cm + 18"/45cm	12	SXPP2B434
	1	Violet	5"/14cm + 5"/14cm	12	SXPP2B401
	1	Violet	9"/24cm + 9"/24cm	12	SXPP2B402
	1	Violet	12"/30cm + 12"/30cm	12	SXPP2B403
	 <p>MO-4 36mm 1/2 Circle</p>	1	Violet	14"/36cm + 14"/36cm	12

Bidirectional STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device





Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 CT 40mm 1/2 Circle	1	Violet	14"/36cm + 14"/36cm	12	SXPP2B435	
 CTX 48mm 1/2 Circle	1	Violet	18"/45cm + 18"/45cm	12	SXPP2B405	
 UR-6 26mm 5/8 Circle	3-0	Violet	6"/16cm + 6"/16cm	12	SXPP2B432	
	2-0	Violet	6"/16cm + 6"/16cm	12	SXPP2B433	





Bidirectional STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device




PRIME Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Product Code
	4-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B423
PS-2 19mm 3/8 Circle	3-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B100
	3-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B417
PS 26mm 3/8 Circle	3-0	Violet	9"/24cm + 9"/24cm	12	SXPP2B419


Cutting Edge Reverse ▾

Needle	Suture USP	Color	Length	EA/BX	Product Code
	2-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B418
FS 26mm 3/8 Circle	2-0	Violet	9"/24cm + 9"/24cm	12	SXPP2B420
	3-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B102
X-1 22mm 1/2 Circle	0	Violet	9"/24cm + 9"/24cm	12	SXPP2B415
	0	Violet	9"/24cm + 9"/24cm	12	SXPP2B415
CP-2 26mm 1/2 Circle	1	Violet	9"/24cm + 9"/24cm	12	SXPP2B200
	1	Violet	12"/30cm + 12"/30cm	12	SXPP2B201
OS-6 36mm 1/2 Circle					



Bidirectional STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device




Cutting Edge Reverse ▾

Needle	Suture USP	Color	Length	EA/BX	Product Code
 OS-8 40mm 1/2 Circle	1	Violet	14"/36cm + 14"/36cm	12	SXPP2B202
	1	Violet	18"/45cm + 18"/45cm	12	SXPP2B203

TAPERCUT™ ▾


Needle	Suture USP	Color	Length	EA/BX	Product Code
 V-7 26mm 1/2 Circle	2-0	Violet	5"/14cm + 5"/14cm	12	SXPP2B426
 V-34 36mm 1/2 Circle	0	Violet	5"/14cm + 5"/14cm	12	SXPP2B425
	0	Violet	9"/24cm + 9"/24cm	12	SXPP2B424






ETHIGUARD™ Blunt Point Needle ●

Needle	Suture USP	Color	Length	EA/BX	Product Code
 CTB-1 36mm 1/2 Circle	0	Violet	5"/14cm + 5"/14cm	12	SXPP2B300
	1	Violet	5"/14cm + 5"/14cm	12	SXPP2B301

Undirectional STRATAFIX™ PDS™ Plus Knotless Tissue Control Device












Taperpoint 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	4-0	Violet	6"/15cm	12	SXPP1B427
	4-0	Violet	12"/30cm	12	SXPP1B428
	3-0	Violet	6"/15cm	12	SXPP1B422
 CT-2 26mm 1/2 Circle	3-0	Violet	27"/70cm	12	SXPP1B425
	2-0	Violet	6"/15cm	12	SXPP1B413
	2-0	Violet	9"/23cm	12	SXPP1B432
	2-0	Violet	27"/70cm	12	SXPP1B414
	0	Violet	12"/30cm	12	SXPP1B405
 SH 26mm 1/2 Circle	4-0	Violet	27"/70cm	12	SXPP1B426
	3-0	Violet	6"/15cm	12	SXPP1B420
	3-0	Violet	27"/70cm	12	SXPP1B421
	2-0	Violet	6"/15cm	12	SXPP1B415
	2-0	Violet	9"/23cm	12	SXPP1B433
	2-0	Violet	12"/30cm	12	SXPP1B416
	2-0	Violet	27"/70cm	12	SXPP1B417
 SH-1 22mm 1/2 Circle	4-0	Violet	9"/23cm	12	SXPP1B454
	3-0	Violet	9"/23cm	12	SXPP1B453
 CT-1 36mm 1/2 Circle	2-0	Violet	6"/15cm	12	SXPP1B409
	2-0	Violet	9"/23cm	12	SXPP1B456
	2-0	Violet	12"/30cm	12	SXPP1B410
	2-0	Violet	18"/45cm	12	SXPP1B411
	2-0	Violet	27"/70cm	12	SXPP1B412
	0	Violet	6"/15cm	12	SXPP1B406
	0	Violet	9"/23cm	12	SXPP1B455
	0	Violet	12"/30cm	12	SXPP1B450
	0	Violet	18"/45cm	12	SXPP1B407
	0	Violet	27"/70cm	12	SXPP1B408
	1	Violet	12"/30cm	12	SXPP1B429
	1	Violet	18"/45cm	12	SXPP1B430
	1	Violet	27"/70cm	12	SXPP1B431

Undirectional STRATAFIX™ PDS™ Plus Knotless Tissue Control Device








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 <p>MO-4 36mm 1/2 Circle</p>	1	Violet	27"/70cm	12	SXPP1B205	
 <p>CTX 48mm 1/2 Circle</p>	1	Violet	27"/70cm	12	SXPP1B401	
 <p>UR-6 26mm 5/8 Circle</p>	2-0	Violet	12"/30cm	12	SXPP1B419	
 <p>EN-S 19mm Ski</p>	0	Violet	12"/30cm	12	SXPP1B452	
 <p>ST-4 19mm Straight</p>	3-0	Violet	6"/15cm	12	SXPP1B424	

TAPERCUT™ 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 <p>V-7 26mm 1/2 Circle</p>	2-0	Violet	12"/30cm	12	SXPP1B418	
 <p>V-34 36mm 1/2 Circle</p>	2-0	Violet	9"/23cm	12	SXPP1B458	
	0	Violet	9"/23cm	12	SXPP1B457	
	0	Violet	12"/30cm	12	SXPP1B404	


Undirectional STRATAFIX™ PDS™ Plus Knotless Tissue Control Device







Cutting Edge Reverse ▼					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 FS 26mm 3/8 Circle	2-0	Violet	18"/45cm	12	SXPP1B101
	2-0	Violet	27"/70cm	12	SXPP1B102
 CP-2 26mm 1/2 Circle	0	Violet	18"/45cm	12	SXPP1B202
	1	Violet	18"/45cm	12	SXPP1B403
 OS-6 36mm 1/2 Circle	1	Violet	18"/45cm	12	SXPP1B201
	1	Violet	27"/70cm	12	SXPP1B203
 OS-8 40mm 1/2 Circle	1	Violet	27"/70cm	12	SXPP1B204
 X-1 22mm 1/2 Circle	3-0	Violet	27"/70cm	12	SXPP1B120


Undirectional STRATAFIX™ PDS™ Plus Knotless Tissue Control Device





PRIME Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 PS-4 16mm 1/2 Circle	4-0	Violet	18"/45cm	12	SXPP1B119
	4-0	Violet	6"/15cm	12	SXPP1B111
 PS-2 19mm 3/8 Circle	4-0	Violet	12"/30cm	12	SXPP1B112
	4-0	Violet	18"/45cm	12	SXPP1B113
	4-0	Violet	27"/70cm	12	SXPP1B114
	3-0	Violet	6"/15cm	12	SXPP1B103
	3-0	Violet	12"/30cm	12	SXPP1B104
	3-0	Violet	18"/45cm	12	SXPP1B105
	3-0	Violet	27"/70cm	12	SXPP1B106
	4-0	Violet	18"/45cm	12	SXPP1B116
 PS-1 24mm 3/8 Circle	4-0	Violet	23"/60cm	12	SXPP1B117
	3-0	Violet	27"/70cm	12	SXPP1B110
	4-0	Violet	18"/45cm	12	SXPP1B118
 PS 26mm 3/8 Circle	3-0	Violet	12"/30cm	12	SXPP1B107
	3-0	Violet	18"/45cm	12	SXPP1B108
	3-0	Violet	27"/70cm	12	SXPP1B109
	3-0	Violet	27"/70cm	12	SXPP1B109





Unidirectional STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device




ETHIGUARD™ Blunt Point Needle ●






Needle	Suture USP	Color	Length	EA/BX	Product Code
 CTB-1 36mm 1/2 Circle	0	Violet	23"/60cm	12	SXPP1A303
	1	Violet	18"/45cm	12	SXPP1A301
	1	Violet	23"/60cm	12	SXPP1A304
 CTXB 48mm 1/2 Circle	0	Violet	18"/45cm	12	SXPP1A302
	0	Violet	23"/60cm	12	SXPP1A305
	1	Violet	18"/45cm	12	SXPP1A300
	1	Violet	23"/60cm	12	SXPP1A306

Taperpoint ○


Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	3-0	Violet	6"/15cm	12	SXPP1A411
	3-0	Violet	12"/30cm	12	SXPP1A427
 SH-1 22mm 1/2 Circle	2-0	Violet	12"/30cm	12	SXPP1A430
 CT-2 26mm 1/2 Circle	2-0	Violet	6"/15cm	12	SXPP1A415
	2-0	Violet	9"/23cm	12	SXPP1A422
	2-0	Violet	18"/45cm	12	SXPP1A408
	0	Violet	9"/23cm	12	SXPP1A446
	0	Violet	18"/45cm	12	SXPP1A407
	1	Violet	6"/15cm	12	SXPP1A419
 SH 26mm 1/2 Circle	3-0	Violet	6"/15cm	12	SXPP1A413
	3-0	Violet	12"/30cm	12	SXPP1A428
	3-0	Violet	18"/45cm	12	SXPP1A410
	2-0	Violet	6"/15cm	12	SXPP1A414
	2-0	Violet	9"/23cm	12	SXPP1A421
	2-0	Violet	12"/30cm	12	SXPP1A431
	2-0	Violet	18"/45cm	12	SXPP1A409

Unidirectional STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device





Taperpoint 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	3-0	Violet	12"/30cm	12	SXPP1A429
	3-0	Violet	18"/45cm	12	SXPP1A438
	2-0	Violet	6"/15cm	12	SXPP1A416
	2-0	Violet	9"/23cm	12	SXPP1A423
	2-0	Violet	12"/30cm	12	SXPP1A432
	2-0	Violet	18"/45cm	12	SXPP1A403
	2-0	Violet	24"/60cm	12	SXPP1A439
	0	Violet	6"/15cm	12	SXPP1A418
	0	Violet	9"/23cm	12	SXPP1A425
	0	Violet	12"/30cm	12	SXPP1A433
	0	Violet	18"/45cm	12	SXPP1A401
	0	Violet	24"/60cm	12	SXPP1A440
	1	Violet	6"/15cm	12	SXPP1A420
	1	Violet	12"/30cm	12	SXPP1A435
	1	Violet	18"/45cm	12	SXPP1A404
1	Violet	24"/60cm	12	SXPP1A443	
 <p>CT 40mm 1/2 Circle</p>	0	Violet	12"/30cm	12	SXPP1A412
	0	Violet	18"/45cm	12	SXPP1A406
	0	Violet	24"/60cm	12	SXPP1A441
	1	Violet	18"/45cm	12	SXPP1A405
	1	Violet	24"/60cm	12	SXPP1A444
 <p>CTX 48mm 1/2 Circle</p>	0	Violet	12"/30cm	12	SXPP1A434
	0	Violet	18"/45cm	12	SXPP1A402
	0	Violet	24"/60cm	12	SXPP1A442
	1	Violet	12"/30cm	12	SXPP1A436
	1	Violet	18"/45cm	12	SXPP1A400
	1	Violet	24"/60cm	12	SXPP1A445
 <p>UR-6 26mm 5/8 Circle</p>	0	Violet	6"/15cm	12	SXPP1A447
	2-0	Violet	6"/15cm	12	SXPP1A417


Unidirectional STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device




TAPERCUT™ 




Needle	Suture USP	Color	Length	EA/BX	Product Code
 V-7 26mm 1/2 Circle	2-0	Violet	9"/23cm	12	SXPP1A424
 V-34 36mm 1/2 Circle	0	Violet	9"/23cm	12	SXPP1A426

Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 CP-2 26mm 1/2 Circle	3-0	Violet	18"/45cm	12	SXPP1A437
 OS-6 36mm 1/2 Circle	0	Violet	18"/45cm	12	SXPP1A200
	0	Violet	24"/60cm	12	SXPP1A202
	1	Violet	18"/45cm	12	SXPP1A201
	1	Violet	24"/60cm	12	SXPP1A203
 OS-8 40mm 1/2 Circle	0	Violet	24"/60cm	12	SXPP1A204
	1	Violet	24"/60cm	12	SXPP1A205

Unidirectional STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device



PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 PS-2 19mm 3/8 Circle	3-0	Violet	18"/45cm	12	SXPP1A101	
 PS-1 24mm 3/8 Circle	3-0	Violet	18"/45cm	12	SXPP1A100	

Sutures: Absorbable




VICRYL RAPIDE™ (polyglactin 910) Suture





VICRYL RAPIDE Suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide.¹



VICRYL RAPIDE™ (polyglactin 910) Suture




Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	4-0	Undyed	1 × 27"/70cm	12		VR214
 SH 26mm 1/2 Circle	3-0	Undyed	1 × 27"/70cm	12		VR416
	2-0	Undyed	1 × 27"/70cm	12		VR417
 CT-1 36mm 1/2 Circle	3-0	Undyed	1 × 36"/90cm	12		VR944
	2-0	Undyed	1 × 36"/90cm	12		VR945
 CTX 48mm 1/2 Circle	0	Undyed	1 × 36"/90cm	12		VR1695


TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-34 36mm 1/2 Circle	3-0	Undyed	1 × 36"/90cm	12		VR916
	2-0	Undyed	1 × 36"/90cm	12		VR917
	0	Undyed	1 × 36"/90cm	12		VR2218












VICRYL RAPIDE™ (polyglactin 910) Suture




Taperpoint , Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1 36mm/ CP-1 36mm 1/2 Circle/ 1/2 Circle	3-0	Undyed	1 × 54"/135cm	12		VR931
	2-0	Undyed	1 × 54"/135cm	12		VR932







PRIME Reverse Cutting 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12		VR490
 P-3 13mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12		VR493
	4-0	Undyed	1 × 18"/45cm	12		VR494
 PS-2 19mm 3/8 Circle	4-0	Undyed	1 × 27"/70cm	12		VR426
	4-0	Undyed	1 × 18"/45cm	12		VR496
	3-0	Undyed	1 × 18"/45cm	12		VR497
 PS-1 24mm 3/8 Circle	3-0	Undyed	1 × 27"/70cm	12		VR935

VICRYL RAPIDE™ (polyglactin 910) Suture



PRIME Conventional Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass™ NEEDLES	VR834
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass™ NEEDLES	VR835
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass™ NEEDLES	VR844
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass™ NEEDLES	VR845



Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture

VICRYL Suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide.¹

VICRYL Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.¹

VICRYL Plus Suture is indicated for use in general soft tissue approximation and/or ligation, except for ophthalmic, cardiovascular and neurological tissues.¹


- Shown in vitro to inhibit bacterial colonization of the suture for 7 days or more²







Reference: 1. Ethicon Wound Closure Manual. 2021. Ethicon, Inc. 2. Rothenburger S, Spangler D, Bhende S, Burkley D. In vitro antimicrobial evaluation of coated Vicryl Plus Antibacterial Suture (coated polyglactin 910 with triclosan) using zone of inhibition assays. *Surg Infect (Larchmt)*. 2002;3(suppl):S79-S87.

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTB-2 26mm 1/2 Circle	2-0	Undyed	1 × 27"/70cm	36	Plus	VCPB269H
	2-0	Violet	8 × 18"/45cm	12	Plus CR	VCPB726D
	0	Violet	8 × 18"/45cm	12	Plus CR	VCPB727D
 SHB 26mm 1/2 Circle	3-0	Undyed	8 × 18"/45cm	12	Plus CR	VCPB864D*
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCPB416H
	3-0	Violet	8 × 18"/45cm	12	Plus CR	VCPB774D*
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCPB417H
 CTB-1 36mm 1/2 Circle	3-0	Undyed	1 × 36"/90cm	36	Plus	VCPB944H
	2-0	Violet	8 × 18"/45cm	12	Plus CR	VCPB739D*
	2-0	Undyed	8 × 18"/45cm	12	Plus CR	VCPB839D*
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCPB259H
	2-0	Undyed	1 × 36"/90cm	36		JB945H
	2-0	Undyed	1 × 36"/90cm	36		VCPB945H
	0	Undyed	1 × 36"/90cm	36	Plus	VCPB946H
	0	Undyed	8 × 18"/45cm	12	Plus CR	VCPB840D*

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTB-1 36mm 1/2 Circle</p>	0	Violet	8 × 18"/45cm	12	Plus CR	VCPB740D*
	0	Undyed	1 × 27"/70cm	36	Plus	VCPB260H
	0	Undyed	8 × 18"/45cm	12	Plus CR	VCPB41D*
	0	Violet	1 × 27"/70cm	36	Plus	VCPB340H
	0	Violet	8 × 18"/45cm	12	Plus CR	VCPB31D*
	1	Undyed	8 × 18"/45cm	12	CR	JB841H
	1	Undyed	1 × 36"/90cm	36		JB947H
	1	Undyed	1 × 36"/90cm	36	Plus	VCPB947H
	1	Violet	1 × 27"/70cm	36	Plus	VCPB341H
	1	Undyed	8 × 18"/45cm	12	Plus CR	VCPB841D*
 <p>CTB 40mm 1/2 Circle</p>	0	Undyed	1 × 36"/90cm	36	Plus	VCPB958H
	 <p>CTXB 48mm 1/2 Circle</p>	0	Undyed	8 × 18"/45cm	12	Plus CR
0		Undyed	1 × 36"/90cm	36		JB978H
0		Undyed	1 × 36"/90cm	36	Plus CR	VCPB978H
1		Undyed	8 × 18"/45cm	12	Plus CR	VCPB725D*
1		Undyed	1 × 36"/90cm	36		JB977H
1		Undyed	1 × 36"/90cm	36	Plus	VCPB977H
1		Violet	8 × 18"/45cm	12	CR	VCPB765D*

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture








Coated VICRYL™ (polyglactin 910) Suture



ETHIGUARD™ Blunt Point Needle


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TPB-1 65mm 1/2 Circle	1	Violet	4 X 27"/70cm	12	Plus	VCPB650G








Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-4 5mm 3/8 Circle	8-O	Violet	1 x 5"/13cm	12	 MultiPass™ NEEDLES	J405G
 BV130-5 6.5mm 3/8 Circle	8-O	Violet	1 x 5"/13cm	12	 MultiPass™ NEEDLES	J401G
 TF 13mm 1/2 Circle	5-0	Undyed	1 x 27"/70cm	36		J433H
 RB-2 13mm 1/2 Circle	4-0	Violet	8 x 18"/45cm	12	Plus CR	VCP743D*
	5-0	Undyed	1 x 27"/70cm	36	Plus	VCP433H
 RB-1 17mm 1/2 Circle	4-0	Undyed	1 x 27"/70cm	36	Plus	VCP434H
	6-0	Undyed	1 x 27"/70cm	36		J212H
	5-0	Undyed	1 x 27"/70cm	36	Plus	VCP213H
	5-0	Violet	1 x 27"/70cm	36	Plus	VCP303H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	4-0	Undyed	8 × 18"/45cm	12		VCP714D*
	4-0	Undyed	1 × 27"/70cm	36		J214H
	4-0	Undyed	1 × 27"/70cm	36	Plus	VCP214H
	4-0	Violet	1 × 27"/70cm	36		J304H
	4-0	Violet	1 × 27"/70cm	36	Plus	VCP304H
	3-0	Violet	8 × 18"/45cm	12		VCP713D*
	3-0	Undyed	1 × 27"/70cm	36		VCP215H
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP305H
	2-0	Violet	1 × 27"/70cm	36		J306H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP306H
 CT-3 22mm 1/2 Circle	3-0	Violet	1 × 27"/70cm	36		J327H
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP327H
	2-0	Violet	1 × 27"/70cm	36		J328H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP328H
	O	Violet	1 × 27"/70cm	36		J329H
	O	Violet	1 × 27"/70cm	36	Plus	VCP329H



*Safety organizer tray

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH-1 22mm 1/2 Circle	4-0	Violet	8 × 18"/45cm	12	Plus	VCP771D*
	4-0	Undyed	1 × 27"/70cm	36	Plus CR	VCP218H
	4-0	Violet	1 × 27"/70cm	36		J310H
	4-0	Violet	1 × 27"/70cm	36	Plus	VCP310H
	3-0	Violet	8 × 18"/45cm	12		VCP772D*
	3-0	Undyed	1 × 27"/70cm	36		J219H
	3-0	Undyed	1 × 27"/70cm	36	Plus CR	VCP219H
 CT-2 26mm 1/2 Circle	3-0	Violet	1 × 27"/70cm	36	Plus	VCP311H
	3-0	Undyed	1 × 27"/70cm	36		J232H
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP232H
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP332H
	2-0	Violet	8 × 18"/45cm	12		J726D*
	2-0	Undyed	8 × 18"/45cm	12		VCP726D*
	2-0	Undyed	1 × 27"/70cm	36		J269H
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP269H
	2-0	Violet	1 × 27"/70cm	36		J333H
	2-0	Violet	1 × 27"/70cm	36	Plus CR	VCP333H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture








Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-2 26mm 1/2 Circle	0	Violet	8 × 18"/45cm	12		J726D*
	0	Violet	8 × 18"/45cm	36	Plus 	VCP727D
	0	Undyed	1 × 27"/70cm	36		VCP270H
	0	Violet	1 × 27"/70cm	36		VCP334H
	1	Violet	8 × 18"/45cm	36		J335H
	1	Violet	1 × 27"/70cm	36		VCP335H
 SH 26mm 1/2 Circle	4-0	Violet	8 × 18"/45cm	12	Plus	VCP773D*
	4-0	Undyed	1 × 27"/70cm	36		J415H
	4-0	Undyed	1 × 27"/70cm	36	Plus	VCP415H
	4-0	Violet	1 × 27"/70cm	36		J315H
	4-0	Violet	1 × 27"/70cm	36	Plus	VCP315H
	4-0	Violet	8 × 18"/45cm	12		J783G
	3-0	Undyed	8 × 18"/45cm	12		J864D
	3-0	Undyed	8 × 18"/45cm	12	Plus 	VCP864D*
	3-0	Violet	8 × 18"/45cm	12		J774D*
	3-0	Violet	8 × 18"/45cm	12	Plus 	VCP774D*
	3-0	Undyed	1 × 27"/70cm	36		VCP416H
	3-0	Violet	1 × 27"/70cm	36		J316H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture
















Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH 26mm 1/2 Circle	3-0	Violet	1 × 27"/70cm	36	Plus	VCP316H
	3-0	Violet	8 × 27"/70cm	12	Plus CR	VCP784D*
	2-0	Violet	8 × 18"/45cm	12	CR	J775D*
	2-0	Violet	8 × 18"/45cm	12	Plus CR	VCP775D*
	2-0	Undyed	1 × 27"/70cm	36		J417H
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP417H
	2-0	Violet	1 × 27"/70cm	36		J317H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP317H
	2-0	Violet	8 × 27"/70cm	12	Plus CR	VCP785D*
	0	Undyed	1 × 27"/70cm	36		J418H
	0	Undyed	1 × 27"/70cm	36	Plus	VCP418H
 SH 26mm 1/2 Circle  MO-6 26mm 1/2 Circle	3-0	Violet	1 × 36"/90cm	36	Plus	VCP527H
	2-0	Violet	8 × 18"/45cm	12	CR	VCP789D*

*Safety organizer tray

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture





Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 MO-5 31mm 1/2 Circle	0	Violet	8 × 18"/45cm	12		J801D*
	0	Violet	8 × 18"/45cm	12	Plus 	VCP801D*
	0	Violet	8 × 27"/70cm	12	Plus 	VCPP81D*
	1	Violet	8 × 27"/70cm	12	Plus 	VCPP80D*
 CT-1 36mm 1/2 Circle	4-0	Violet	1 × 27"/70cm	36		J337H
	3-0	Undyed	8 × 18"/45cm	12	Plus 	VCP838D*
	3-0	Violet	8 × 18"/45cm	12	Plus 	VCP738D*
	3-0	Undyed	1 × 27"/70cm	36		J258H
	3-0	Violet	1 × 27"/70cm	36		J338H
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP338H
	3-0	Undyed	1 × 36"/90cm	36		J944H
	3-0	Undyed	1 × 36"/90cm	36	Plus	VCP944H
	3-0	Violet	1 × 36"/90cm	36		J344H
	3-0	Violet	1 × 36"/90cm	36	Plus	VCP344H
	2-0	Undyed	8 × 18"/45cm	12		J839D*
	2-0	Undyed	8 × 18"/45cm	12	Plus 	VCP839D*
	2-0	Violet	8 × 18"/45cm	12		J739D*

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture





Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Violet	8 × 18"/45cm	12	Plus CR	VCP739D*
	2-0	Undyed	1 × 27"/70cm	36		J259H
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP259H
	2-0	Violet	1 × 27"/70cm	36		J339H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP339H
	2-0	Undyed	8 × 27"/70cm	12		VCPP42D*
	2-0	Undyed	1 × 36"/90cm	36	Plus CR	J945H
	2-0	Undyed	1 × 36"/90cm	36	Plus	VCP945H
	2-0	Violet	1 × 36"/90cm	36		J345H
	2-0	Violet	1 × 36"/90cm	36	Plus	VCP345H
	O	Undyed	8 × 18"/45cm	12		VCP840D*
	O	Violet	8 × 18"/45cm	12	CR	J740D*
	O	Violet	8 × 18"/45cm	12	Plus CR	VCP740D*
	O	Undyed	1 × 27"/70cm	36		J260H
	O	Undyed	1 × 27"/70cm	36	Plus CR	VCP260H
	O	Violet	1 × 27"/70cm	36		J340H
	O	Violet	1 × 27"/70cm	36	Plus	VCP340H
	0	Undyed	8 × 27"/70cm	12		VCPP41D*

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture







Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	0	Violet	8 × 27"/70cm	36	Plus CR	VCPP31D*
	0	Undyed	1 × 36"/90cm	36		J946H
	0	Undyed	1 × 36"/90cm	36	Plus	VCP946H
	0	Violet	1 × 36"/90cm	36		J346H
	0	Violet	1 × 36"/90cm	36		VCP346H
	1	Undyed	8 × 18"/45cm	12	Plus CR	VCP841D*
	1	Violet	8 × 18"/45cm	12	CR	J741D*
	1	Violet	8 × 18"/45cm	12	Plus CR	VCP741D*
	1	Undyed	1 × 27"/70cm	36		J261H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP261H
	1	Violet	1 × 27"/70cm	36	Plus	VCP341H
	1	Undyed	8 × 27"/70cm	12	Plus CR	VCPP40D*
	1	Undyed	1 × 36"/90cm	36		J947H
	1	Undyed	1 × 36"/90cm	36	Plus	VCP947H
1	Violet	1 × 36"/90cm	36		J347H	
1	Violet	1 × 36"/90cm	36	Plus	VCP347H	
 <p>MH 36mm 1/2 Circle</p>	3-0	Violet	1 × 27"/70cm	36		J322H
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP322H
	2-0	Violet	1 × 27"/70cm	36		J323H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP323H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>MO-4 36mm 1/2 Circle</p>	2-0	Violet	8 × 18"/45cm	12	Plus CR	VCP700D*
	0	Violet	8 × 27"/70cm	12	Plus CR	VCPP71D*
	0	Violet	8 × 18"/45cm	12	Plus CR	VCP701D*
	0	Violet	1 × 27"/70cm	36		J436H
	0	Violet	1 × 27"/70cm	36	Plus	VCP436H
	1	Violet	8 × 18"/45cm	12	Plus CR	VCP702D*
	1	Violet	1 × 27"/70cm	36		J437H
	1	Violet	1 × 27"/70cm	36	Plus	VCP437H
 <p>CT 40mm 1/2 Circle</p>	3-0	Undyed	1 × 36"/90cm	36	Plus	VCP954H
	3-0	Violet	1 × 36"/90cm	36	Plus	VCP356H
	2-0	Violet	8 × 18"/45cm	12	Plus CR	VCP751D*
	2-0	Undyed	1 × 27"/70cm	36		J275H
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP275H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP351H
	2-0	Undyed	1 × 36"/90cm	36	Plus	VCP957H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture




Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT 40mm 1/2 Circle</p>	2-0	Violet	1 × 36"/90cm	36	Plus	VCP357H
	0	Violet	8 × 18"/45cm	12	Plus CR	VCP752D*
	0	Undyed	1 × 27"/70cm	36	Plus	VCP280H
	0	Violet	1 × 27"/70cm	36		J352H
	0	Violet	1 × 27"/70cm	36	Plus	VCP352H
	0	Undyed	1 × 36"/90cm	36		J958H
	0	Undyed	1 × 36"/90cm	36	Plus	VCP958H
	0	Violet	1 × 36"/90cm	36	Plus	VCP358H
	1	Violet	8 × 18"/45cm	12	Plus CR	VCP753D*
	1	Undyed	1 × 27"/70cm	36		J281H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP281H
	1	Violet	1 × 27"/70cm	36		J353H
	1	Violet	1 × 27"/70cm	36	Plus	VCP353H
	1	Undyed	1 × 36"/90cm	36		J959H
	1	Undyed	1 × 36"/90cm	36	Plus	VCP959H
	1	Violet	1 × 36"/90cm	36		J359H
	1	Violet	1 × 36"/90cm	36	Plus	VCP359H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture Coated VICRYL™ (polyglactin 910) Suture










Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	3-0	Undyed	1 × 36"/90cm	36		J980H
	3-0	Undyed	1 × 36"/90cm	36	Plus	VCP980H
	2-0	Undyed	8 × 18"/45cm	12	Plus CR	VCP723D*
	2-0	Violet	1 × 27"/70cm	36		J363H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP363H
	2-0	Undyed	1 × 36"/90cm	36		J979H
	2-0	Undyed	1 × 36"/90cm	36	Plus	VCP979H
	2-0	Violet	1 × 36"/90cm	36		J369H
	2-0	Violet	1 × 36"/90cm	36	Plus	VCP369H
	0	Undyed	8 × 18"/45cm	12	Plus CR	VCP724D*
	0	Violet	8 × 18"/45cm	36	Plus	VCP764D*
	0	Violet	1 × 27"/70cm	36		J364H
	0	Violet	1 × 27"/70cm	36	Plus	VCP364H
	0	Undyed	1 × 36"/90cm	36		J978H
	0	Undyed	1 × 36"/90cm	36	Plus	VCP978H
	0	Violet	1 × 36"/90cm	36		J370H
	0	Violet	1 × 36"/90cm	36	Plus	VCP370H
	1	Violet	8 × 18"/45cm	12	Plus CR	VCP765D*
	1	Violet	1 × 27"/70cm	36		J365H
	1	Violet	1 × 27"/70cm	36	Plus	VCP365H
1	Undyed	1 × 36"/90cm	36		J977H	
1	Undyed	1 × 36"/90cm	36	Plus	VCP977H	
1	Violet	1 × 36"/90cm	36		J371H	
1	Violet	1 × 36"/90cm	36	Plus	VCP371H	

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>TP-1 65mm 1/2 Circle</p>	1	Violet	4 × 27"/70cm	12		J650G
	2	Undyed	2 × 27"/70cm	12	Plus	VCP849G
	2	Violet	4 × 27"/70cm	12		J649G
	2	Violet	4 × 27"/70cm	12	Plus	VCP649G
	2	Undyed	1 × 54"/135cm	24	Plus	VCP880T
 <p>XLH 70mm 1/2 Circle</p>	2-0	Violet	1 × 27"/70cm	12		J581G
	2-0	Violet	1 × 27"/70cm	12	Plus	VCP581G
	0	Violet	1 × 27"/70cm	12		J582G
	1	Violet	1 × 27"/70cm	12		J583G
	1	Violet	1 × 27"/70cm	12	Plus	VCP583G
 <p>UR-6 26mm 5/8 Circle</p>	2-0	Violet	1 × 27"/70cm	36	Plus	VCP602H
	0	Violet	1 × 27"/70cm	36	Plus	VCP603H
 <p>UR-5 36mm 5/8 Circle</p>	2-0	Violet	1 × 27"/70cm	36	Plus	VCP375H
	0	Violet	1 × 27"/70cm	36	Plus	VCP376H
 <p>UR-4 40mm 5/8 Circle</p>	0	Violet	1 × 27"/70cm	36		J381H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture



Coated VICRYL™ (polyglactin 910) Suture



TAPERCUT™ 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V130-4 4.2mm 1/4 Circle	8-0	Violet	1 × 9"/23cm	12		J408G
 V130-5 5.8mm 1/4 Circle	8-0	Violet	1 × 9"/23cm	12		J409G
 V-34 36mm 1/2 Circle	3-0	Undyed	1 × 36"/90cm	36		J916H
	3-0	Undyed	1 × 36"/90cm	36	Plus	VCP915H
	3-0	Violet	1 × 36"/90cm	36	Plus	VCP516H
	2-0	Undyed	1 × 36"/90cm	36		J917H
	2-0	Undyed	1 × 36"/90cm	36	Plus	VCP917H
	2-0	Violet	1 × 36"/90cm	36		J517H
	2-0	Violet	1 × 36"/90cm	36	Plus	VCP517H
	0	Violet	1 × 36"/90cm	36		J518H
	0	Violet	1 × 36"/90cm	36	Plus	VCP518H
	1	Violet	1 × 36"/90cm	36		J519H
1	Violet	1 × 36"/90cm	36	Plus	VCP519H	

Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-3 13mm 3/8 Circle	5-0	Violet	1 × 18"/45cm	36		J385H
	4-0	Violet	1 × 18"/45cm	36		J386H
 FS-2 19mm 3/8 Circle	5-0	Violet	1 × 18"/45cm	36	Plus	VCP391H
	5-0	Undyed	1 × 27"/70cm	36		J421H
	5-0	Undyed	1 × 27"/70cm	36	Plus	VCP421H
	4-0	Violet	1 × 18"/45cm	36		J392H
	4-0	Violet	1 × 18"/45cm	36	Plus	VCP392ZH
	4-0	Undyed	1 × 27"/70cm	36	Plus	VCP422H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture






Reverse Cutting ▾							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
	3-0	Violet	1 × 18"/45cm	36	Plus	VCP393H	
FS-2 19mm 3/8 Circle	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP423H	
	4-0	Undyed	1 × 27"/70cm	36	Plus	VCP441H	
	4-0	Violet	1 × 27"/70cm	36		J451H	
	4-0	Violet	1 × 27"/70cm	36	Plus	VCP451H	
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP452H	
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP443H	
	3-0	Violet	8 × 18"/45cm	12	Plus CR	VCP787D*	
	3-0	Violet	8 × 18"/45cm	12	Plus CR	VCP787Z	
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP588H	
	FSL 30mm 3/8 Circle						
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP502H	
	PSL 30mm 3/8 Circle						
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP596H	
	PSL 30mm 3/8 Circle						
	4-0	Undyed	1 × 27"/70cm	36		J455H	
	J-1 17mm 1/2 Circle						
	2-0	Undyed	3 × 18"/45cm	24	Plus CR	VCP748T	
	2-0	Undyed	3 × 18"/45cm	24	Plus CR	VCP706T	
	2-0	Undyed	1 × 27"/70cm	36		J693H	
OS-4 22mm 1/2 Circle							

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture





Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 OS-4 22mm 1/2 Circle	0	Undyed	3 × 18"/45cm	24	Plus CR	VCP5749T
	0	Violet	3 × 18"/45cm	24	CR	J707T
	0	Undyed	1 × 27"/70cm	36	Plus	VCP694H
	1	Violet	3 × 18"/45cm	24	Plus CR	VCP708T
	1	Undyed	1 × 27"/70cm	36	Plus	VCP695H
 X-1 22mm 1/2 Circle	3-0	Violet	8 × 18"/45cm	12	Plus CR	VCP790D*
	3-0	Undyed	1 × 27"/70cm	36		J458H
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP458H
	3-0	Violet	1 × 27"/70cm	36		J460H
	3-0	Violet	1 × 27"/70cm	36	Plus	VCP460H
	2-0	Undyed	1 × 27"/70cm	36		J459H
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP459H
	2-0	Violet	1 × 27"/70cm	36		J461H
 CP-2 26mm 1/2 Circle	3-0	Undyed	8 × 18"/45cm	12	Plus	VCP761D*
	3-0	Undyed	1 × 27"/70cm	36		J868H
	3-0	Undyed	8 × 18"/45cm	12	Plus	VCP761D*
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP868H
	2-0	Undyed	8 × 18"/45cm	12	Plus CR	VCP762D*

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture






Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-2 26mm 1/2 Circle	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP869H
	0	Undyed	1 × 27"/70cm	36		J870H
	0	Undyed	1 × 27"/70cm	36	Plus	VCP870H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP871H
	1	Undyed	1 × 27"/70cm	36		J871H
 CP-1 36mm 1/2 Circle	3-0	Violet	1 × 36"/90cm	36		J471H
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP265H
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP266H
	2-0	Violet	1 × 27"/70cm	36	Plus	VCP466H
	2-0	Violet	1 × 36"/90cm	36	Plus	VCP472H
	0	Undyed	1 × 27"/70cm	36		J267H
	0	Undyed	1 × 27"/70cm	36	Plus	VCP267H
	0	Violet	1 × 27"/70cm	36		J467H
	0	Violet	1 × 27"/70cm	36	Plus	VCP467H
	0	Violet	1 × 36"/90cm	36		J473H
	1	Undyed	1 × 27"/70cm	36		J268H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP268H
	1	Violet	1 × 27"/70cm	36		J468H
	1	Violet	1 × 27"/70cm	36	Plus	VCP468H
	1	Violet	1 × 36"/90cm	36		J474H
1	Violet	1 × 36"/90cm	36	Plus	VCP474H	

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture





Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>OS-6 36mm 1/2 Circle</p>	2-0	Undyed	3 × 18"/45cm	24	Plus CR	VCP750T
	2-0	Violet	3 × 18"/45cm	24	Plus CR	VCP710T
	2-0	Undyed	1 × 27"/70cm	36	Plus	VCP533H
	0	Undyed	3 × 18"/45cm	24	Plus CR	VCP754T
	0	Violet	3 × 18"/45cm	24	Plus CR	VCP711T
	0	Undyed	1 × 27"/70cm	36	Plus	VCP534H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP535H
 <p>CP 40mm 1/2 Circle</p>	2-0	Violet	1 × 27"/70cm	36	Plus	VCP478H
	0	Undyed	1 × 27"/70cm	36	Plus	VCP197H
	0	Violet	1 × 27"/70cm	36		J479H
	0	Violet	1 × 36"/90cm	36		J485H
	1	Undyed	1 × 27"/70cm	36		J196H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP196H
	1	Violet	1 × 27"/70cm	36		J480H
	1	Violet	1 × 27"/70cm	36	Plus	VCP480H
	1	Violet	1 × 36"/90cm	36		J486H
	1	Violet	1 × 36"/90cm	36	Plus	VCP486H
	2	Undyed	1 × 27"/70cm	36		J195H
	2	Undyed	1 × 27"/70cm	36	Plus	VCP195H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture





Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>OS-8 40mm 1/2 Circle</p>	2-0	Violet	3 × 18"/45cm	24	Plus CR	VCP716T
	0	Undyed	3 × 18"/45cm	24	CR	J756T
	0	Undyed	3 × 18"/45cm	24	Plus CR	VCP756T
	0	Violet	3 × 18"/45cm	24	Plus CR	VCP717T
	0	Undyed	1 × 27"/70cm	36	Plus	VCP698H
	1	Undyed	3 × 18"/45cm	24	CR	J757T
	1	Undyed	3 × 18"/45cm	24	Plus CR	VCP757T
	1	Violet	3 × 18"/45cm	36	Plus CR	VCP718T
	1	Undyed	1 × 27"/70cm	36		J699H
	1	Undyed	1 × 27"/70cm	24	Plus CR	VCP699H
 <p>CPX 48mm 1/2 Circle</p>	2	Violet	3 × 18"/45cm		Plus CR	VCP719T
	0	Undyed	1 × 27"/70cm	36		J568H
	0	Undyed	1 × 27"/70cm	36	Plus	VCP568H
	1	Undyed	1 × 27"/70cm	36		J569H
	1	Undyed	1 × 27"/70cm	36	Plus	VCP569H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture



Taperpoint  , Reverse Cutting 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT 40mm/ FSL 30mm 1/2 Circle, 3/8 Circle	3-0	Undyed	1 × 54"/135cm	36	Plus	VCP940H
	3-0	Violet	1 × 54"/135cm	36		J531H
 CT-1 36mm/ CP-1 36mm 1/2 Circle/ 1/2 Circle	3-0	Undyed	1 × 54"/135cm	36	Plus	VCP931H
	2-0	Violet	1 × 54"/135cm	36		J532H

Straight Cutting 























Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Undyed	1 × 27"/70cm	36		J662H
	4-0	Undyed	1 × 27"/70cm	36	Plus	VCP662H
	3-0	Undyed	1 × 27"/70cm	36	Plus	VCP663H
	3-0	Violet	1 × 27"/70cm	36		J523H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




PRIME Reverse Cutting






















Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-6 8mm 1/4 Circle	6-0	Violet	1 × 18"/45cm	12		J551G
	7-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J488G
 P-1 11mm 3/8 Circle	6-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J489G
	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J490G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCPJ490G
 P-3 13mm 3/8 Circle	6-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J492G
	6-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J492H
	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J493G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP493G
	5-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J493H
	5-0	Undyed	1 × 18"/45cm	36	Plus  MultiPass [®]	VCP493H
	5-0	Violet	1 × 18"/45cm	12	 MultiPass [®]	J463G
	5-0	Violet	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP463G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J494G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP494G
	4-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J494H
	4-0	Undyed	1 × 18"/45cm	36	Plus  MultiPass [®]	VCP494H
	4-0	Violet	1 × 18"/45cm	12	 MultiPass [®]	J464G
	4-0	Violet	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP464G
	 PSL 30mm 3/8 Circle	0	Violet	1 × 27"/70cm	36	

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




PRIME Reverse Cutting

















Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-3 16mm 3/8 Circle	6-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J499G
	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J500G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP500G
	5-0	Undyed	1 × 18"/45cm	36	Plus  MultiPass [®]	VCP500H
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J495G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP495G
	5-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J495H
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP495G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J496G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP496ZG
	4-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J496H
	4-0	Undyed	1 × 18"/45cm	36	Plus  MultiPass [®]	VCP496ZH
	4-0	Undyed	1 × 27"/70cm	36	Plus  MultiPass [®]	VCP426H
	3-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J497G
	3-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J497H
	3-0	Undyed	1 × 18"/45cm	36	Plus  MultiPass [®]	VCP497H
	3-0	Undyed	1 × 27"/70cm	36	Plus  MultiPass [®]	VCP427H
	2-0	Undyed	1 × 27"/70cm	36	 MultiPass [®]	J428H
	2-0	Undyed	1 × 27"/70cm	36	Plus  MultiPass [®]	VCP428H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




PRIME Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 × 18"/45cm	12	Plus 	VCP682G
	4-0	Undyed	1 × 27"/70cm	36	Plus 	VCP935H
	3-0	Undyed	1 × 18"/45cm	12	Plus 	VCP683G
	3-0	Undyed	1 × 18"/45cm	36	Plus 	VCP683H
	3-0	Undyed	1 × 27"/70cm	36	Plus 	VCP936H
	4-0	Undyed	1 × 18"/45cm	36		J682H
	4-0	Undyed	1 × 18"/45cm	36	Plus 	VCP682J
	4-0	Undyed	1 × 27"/70cm	36		J935H
	4-0	Undyed	1 × 27"/70cm	36	Plus 	VCP935H
	3-0	Undyed	1 × 18"/45cm	12	Plus 	VCP683G
	3-0	Undyed	1 × 18"/45cm	36		J683H
	3-0	Undyed	1 × 18"/45cm	36	Plus 	VCP683H
	3-0	Undyed	1 × 27"/70cm	36	Plus 	VCP936H
 P-2 8mm 1/2 Circle	5-0	Undyed	1 × 18"/45cm	12	Plus	VCP503G
	4-0	Undyed	1 × 18"/45cm	12	Plus	VCP504G
 PS-6 10mm 1/2 Circle	6-0	Undyed	1 × 18"/45cm	12		J510G
	5-0	Undyed	1 × 18"/45cm	12	Plus	VCP511G

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture



PRIME Reverse Cutting ▾


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-5 13mm 1/2 Circle	5-0	Undyed	1 × 18"/45cm	12	Plus	VCP593G
	4-0	Undyed	1 × 18"/45cm	12	Plus	VCP594G
 PS-4 16mm 1/2 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass	J506G
	5-0	Undyed	1 × 18"/45cm	12	Plus 	VCP506G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass	J507G
	4-0	Undyed	1 × 18"/45cm	12	Plus 	VCP507G
 PS-4C 16mm Compound Curve	4-0	Undyed	1 × 18"/45cm	12		J656G
	4-0	Undyed	1 × 18"/45cm	12	Plus	VCP656G
	3-0	Undyed	1 × 18"/45cm	12		J657G

Cutting Edge Conventional ▲





















Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CCS-1 60mm 1/2 Circle	3	Undyed	1 × 27"/70cm	24		J808T

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




PRIME Conventional Cutting 



Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>PC-1 13mm 3/8 Circle</p>	6-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J833G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP834G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP835G
 <p>PC-3 16mm 3/8 Circle</p>	6-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J843G
	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J844G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP844G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J845G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP845G
 <p>PC-5 19mm 3/8 Circle</p>	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP822G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J823G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP823G
	4-0	Undyed	1 × 18"/45cm	36	 MultiPass [®]	J823H
	4-0	Undyed	1 × 18"/45cm	36	Plus  MultiPass [®]	VCP823H
	3-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J824G
	3-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	VCP824G
	3-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	J824G
 <p>PC-12 16mm 1/2 Circle</p>	5-0	Undyed	1 × 18"/45cm	12		J813G
	4-0	Undyed	1 × 18"/45cm	12		J814G

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture


Coated VICRYL™ (polyglactin 910) Suture





Spatula Conventional Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 S-29 7.6mm 1/4 Circle	6-0	Violet	1 × 12"/30cm	12		J556G
	6-0	Violet	1 × 18"/45cm	12		J555G
 S-28 7.6mm 1/2 Circle	6-0	Violet	1 × 18"/45cm	12		J562G

SABRELOC™ Centerpoint Spatula 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RD-1 8mm 1/4 Circle	6-0	Violet	1 × 12"/30cm	12		J446G

Spatula Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG100-8 6.5mm 1/4 Circle	6-0	Violet	1 × 18"/45cm	12		J544G
 S-14 8mm 1/4 Circle	6-0	Violet	1 × 8"/20cm	12		J578G
	6-0	Violet	1 × 8"/20cm	12		J590G
	6-0	Violet	1 × 12"/30cm	12		J552G
	6-0	Undyed	1 × 18"/45cm	12		J670G
	6-0	Violet	1 × 18"/45cm	12		J570G
	5-0	Violet	1 × 8"/20cm	12		J579G
	5-0	Violet	1 × 8"/20cm	12		J591G
	5-0	Violet	1 × 12"/30cm	12		J553G
	5-0	Undyed	1 × 18"/45cm	12		J671G
	5-0	Violet	1 × 18"/45cm	12		J571G


Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture



Spatula Reverse Cutting




Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG140-6 6.5mm 3/8 Circle	8-0	Violet	1 × 12"/30cm	12		J974G
 TG140-8 6.5mm 3/8 Circle	8-0	Violet	1 × 8"/20cm	12		J547G
	8-0	Violet	1 × 12"/30cm	12		J548G
	8-0	Violet	1 × 12"/30cm	12		J599G*
	7-0	Violet	1 × 12"/30cm	12		J566G
	7-0	Violet	1 × 18"/45cm	12		J546G
 TG160-6 5.5mm 1/2 Circle	8-0	Violet	1 × 12"/30cm	12		J975G
 TG160-8 5.5mm 1/2 Circle	8-0	Violet	1 × 8"/20cm	12		J574G
	7-0	Violet	1 × 12"/30cm	12		J575G
	7-0	Violet	1 × 18"/45cm	12		J576G

*Beaded retraction suture

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture




SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Undyed	3 × 18"/45cm	36		J643H
	4-0	Violet	3 × 18"/45cm	36		J633H
	4-0	Undyed	12 × 18"/45cm	12	Plus	VCP109G
	4-0	Violet	12 × 18"/45cm	12	Plus	VCP103G
	3-0	Undyed	3 × 18"/45cm	36		J644H
	3-0	Violet	3 × 18"/45cm	36		J634H
	3-0	Violet	3 × 18"/45cm	24	Plus	VCP634T
	3-0	Undyed	12 × 18"/45cm	12	Plus	VCP1110G
	3-0	Undyed	12 × 18"/45cm	24		J910T
	3-0	Violet	12 × 18"/45cm	12		J104G*
	3-0	Violet	12 × 18"/45cm	12	Plus	VCP104G
	3-0	Undyed	1 × 54"/135cm	36		J606H
	3-0	Violet	1 × 54"/135cm	36		J614H
	2-0	Undyed	3 × 18"/45cm	36		J645H
	2-0	Undyed	3 × 18"/45cm	24	Plus	VCP645T
	2-0	Violet	3 × 18"/45cm	36		J635H
	2-0	Violet	3 × 18"/45cm	24	Plus	VCP635T
	2-0	Undyed	12 × 18"/45cm	12	Plus	VCP111G
	2-0	Violet	12 × 18"/45cm	12	Plus	VCP105G
	2-0	Undyed	1 × 54"/135cm	36	Plus	VCP607H
	2-0	Violet	1 × 54"/135cm	36	Plus	VCP615H

Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture



SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	0	Undyed	3 × 18"/45cm	36		J646H
	0	Undyed	3 × 18"/45cm	24	Plus	VCP646T
	0	Violet	3 × 18"/45cm	36		J636H
	0	Violet	3 × 18"/45cm	24	Plus	VCP636T
	0	Undyed	6 × 18"/45cm	12	Plus	VCP112G
	0	Violet	6 × 18"/45cm	12	Plus	VCP106G
	0	Undyed	12 × 18"/45cm	12		J912G
	0	Violet	12 × 18"/45cm	12		J906G
	0	Undyed	1 × 54"/135cm	36	Plus	VCP608H
	0	Violet	1 × 54"/135cm	36	Plus	VCP616H


SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	1	Undyed	6 × 18"/45cm	24		J113T*
	1	Undyed	12 × 18"/45cm	12		J913G
	1	Violet	12 × 18"/45cm	12		J907G
	1	Violet	6 × 18"/45cm	12	Plus	VCP107G
	1	Violet	1 × 54"/135cm	36	Plus	VCP617H
	2	Violet	1 × 54"/135cm	36	Plus	VCP618H

*Single strand delivery


Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture

Coated VICRYL™ (polyglactin 910) Suture



LIGAPAK™ Dispensing Reel


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Undyed	1 × 54"/135cm	12		J283G
	4-0	Undyed	1 × 54"/135cm	12		J284G
	4-0	Undyed	1 × 54"/135cm	12	Plus	VCP284G
	4-0	Violet	1 × 54"/135cm	12		J204G
	3-0	Undyed	1 × 54"/135cm	12		J285G
	3-0	Undyed	1 × 54"/135cm	12	Plus	VCP285G
	3-0	Violet	1 × 54"/135cm	12	Plus	VCP205G
	3-0	Undyed	1 × 144"/366cm	12		J885G
	3-0	Undyed	1 × 144"/366cm	12	Plus	VCP285G
	2-0	Undyed	1 × 54"/135cm	12	Plus	VCP286G
	2-0	Violet	1 × 54"/135cm	12		J206G
	0	Undyed	1 × 54"/135cm	12	Plus	VCP287G
	0	Violet	1 × 54"/135cm	12	Plus	VCP207G
	1	Violet	1 × 54"/135cm	12	Plus	VCP208G




VICRYL™ (polyglactin 910) Suture



VICRYL™ (polyglactin 910) Suture








Taperpoint


Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 BV100-3 3.8mm 3/8 Circle	9-0	Violet/Monofilament	1 × 5"/13cm	12		V402G*

Spatula Reverse Cutting

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 TG140-8 6.5mm 3/8 Circle	9-0	Violet/Braided	1 × 12"/30cm	12		V5495G*

ULTIMA™ Reverse Cutting

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 CS90-6 7mm 1/4 Circle	10-0	Violet/Monofilament	1 × 4"/10cm	12		V980G*
 CS140-6 6.5mm 3/8 Circle	10-0	Violet/Monofilament	1 × 4"/10cm	12		V960G*
	9-0	Violet/Braided	1 × 4"/10cm	12		V950G*
 CS160-6 5.5mm 1/2 Circle	10-0	Violet/Monofilament	1 × 12"/30cm	12		V448G*
 CS160-8 5.5mm 1/2 Circle	10-0	Violet/Monofilament	1 × 12"/30cm	12		V450G*
	9-0	Violet/Monofilament	1 × 12"/30cm	12		V449G*
 CS-B-6 5mm Compound Curve	10-0	Violet/Monofilament	1 × 4"/10cm	12		V966G*



MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture



MONOCRYL™ (poliglecaprone 25) Suture

MONOCRYL Suture is a monofilament synthetic absorbable surgical suture prepared from a copolymer of glycolide and epsilon-caprolactone.¹

MONOCRYL Sutures are indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.¹


MONOCRYL Plus Sutures are ideal for subcuticular skin closure.

- High initial breaking strength to maintain wound approximation during the critical first few days of skin healing¹
- Shown in vitro to inhibit bacterial colonization of the suture for 7 days or more²






MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture


MONOCRYL™ Plus Suture




ETHIGUARD™ Blunt Point Needle

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SHB 26mm 1/2 Circle	3-0	Undyed	1 × 27"/70cm	36		YB416
	3-0	Undyed	1 × 27"/70cm	36	Plus	MCPB416H
	3-0	Violet	1 × 27"/70cm	36	Plus	MCPB316H
 CTB-1 36mm 1/2 Circle	3-0	Undyed	1 × 36"/90cm	36		YB944
	3-0	Undyed	1 × 36"/90cm	36	Plus	MCPB944H
	2-0	Undyed	1 × 36"/90cm	36	Plus	MCPB945H
	2-0	Violet	1 × 27"/70cm	36		YB257
	0	Violet	1 × 27"/70cm	36		YB258
	0	Undyed	1 × 36"/90cm	36		YB946
	0	Undyed	1 × 36"/90cm	36	Plus	MCPB946H
	0	Violet	1 × 36"/90cm	36	Plus	MCPB346H
	1	Undyed	1 × 36"/90cm	36		YB947
 CTXB 48mm 1/2 Circle	0	Undyed	1 × 36"/90cm	36		YB978H

Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	6-0	Undyed	1 × 27"/70cm	36		Y432H
	5-0	Undyed	1 × 27"/70cm	36		Y433H
	4-0	Undyed	1 × 27"/70cm	36		Y434H

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture




MONOCRYL™ Plus Suture

Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	5-0	Undyed	1 × 27"/70cm	36		Y213H
	5-0	Undyed	1 × 27"/70cm	36	Plus	MCP213H
	4-0	Undyed	1 × 27"/70cm	36		Y214H
	4-0	Undyed	1 × 27"/70cm	36	Plus	MCP214H
	3-0	Undyed	1 × 27"/70cm	36		Y215H
	5-0	Violet	1 × 27"/70cm	36		Y303H
	4-0	Violet	1 × 27"/70cm	36		Y304H
	4-0	Violet	1 × 27"/70cm	36	Plus	MCP304H
	3-0	Violet	1 × 27"/70cm	36		Y305H
 SH-1 22mm 1/2 Circle	3-0	Undyed	1 × 27"/70cm	36		Y219H
	3-0	Undyed	8 × 18"/45cm	12	CR	Y872D*
 SH 26mm 1/2 Circle	4-0	Undyed	1 × 27"/70cm	36		Y415H
	4-0	Violet	1 × 27"/70cm	36		Y315H
	4-0	Violet	1 × 27"/70cm	36	Plus	MCP415H
	4-0	Violet	1 × 27"/70cm	36	Plus	MCP315H
	3-0	Undyed	1 × 27"/70cm	36		Y416H
	3-0	Undyed	1 × 27"/70cm	36	Plus	MCP416H
	3-0	Violet	1 × 27"/70cm	36	Plus	MCP316H
	2-0	Undyed	1 × 27"/70cm	36		Y417H
	2-0	Undyed	1 × 27"/70cm	36	Plus	MCP417H
	2-0	Violet	1 × 27"/70cm	36		Y317H
	2-0	Violet	1 × 27"/70cm	36	Plus	MCP317H
	2.0	Violet	1 × 36"/90cm	36		Y732H
	0	Violet	1 × 27"/70cm	36		Y318H
 SH 26mm 1/2 Circle	3-0	Violet	1 × 36"/90cm	36		Y527H

*Safety organizer tray

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture




MONOCRYL™ Plus Suture

Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-2 26mm 1/2 Circle	2-0	Violet	1 × 27"/70cm	36		Y333H
	2-0	Violet	1 × 27"/70cm	36	Plus	MCP333H
	2-0	Violet	1 × 36"/90cm	36		Y762H
	0	Violet	1 × 27"/70cm	36		Y334H
	1	Violet	1 × 27"/70cm	36		Y335H
 MO-5 31mm 1/2 Circle	0	Violet	8 × 18"/45cm	12	CR	Y801D*
 MH 36mm 1/2 Circle	2-0	Undyed	1 × 27"/70cm	36		Y227H
 CT-1 36mm 1/2 Circle	3-0	Violet	8 × 18"/45cm	12	CR	Y738D*
	3-0	Violet	1 × 27"/70cm	36		Y338H
	3-0	Undyed	1 × 36"/90cm	36		Y944H
	3-0	Undyed	1 × 36"/90cm	36	Plus	MCP944H
	3-0	Violet	1 × 36"/90cm	36		Y344H
	3-0	Violet	1 × 36"/90cm	36	Plus	MCP344H
	2-0	Violet	8 × 18"/45cm	12	CR	Y739D*
	2-0	Undyed	1 × 36"/90cm	36		Y945H
	2-0	Undyed	1 × 36"/90cm	36	Plus	MCP945H
	2-0	Violet	1 × 36"/90cm	36		Y345H
	2-0	Violet	1 × 36"/90cm	36	Plus	MCP345H
	0	Violet	8 × 18"/45cm	12	CR	Y740D*
	0	Violet	1 × 27"/70cm	36		Y340H
	0	Undyed	1 × 36"/90cm	36		Y946H
	0	Violet	1 × 36"/90cm	36		Y346H

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture



MONOCRYL™ Plus Suture


Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1 36mm 1/2 Circle	0	Violet	4 × 27"/70cm	12	CR	YY31G
	0	Undyed	1 × 36"/90cm	36	Plus	MCP946H
	0	Violet	1 × 36"/90cm	36	Plus	MCP346H
	0	Violet	1 × 27"/70cm	36	Plus	MCP340H
	1	Violet	1 × 27"/70cm	36		Y341H
	1	Violet	4 × 27"/70cm	12	CR	YY30G
	1	Undyed	1 × 36"/90cm	36		Y947H
	1	Violet	1 × 36"/90cm	36		MCP347H
	1	Violet	1 × 36"/90cm	36	Plus	Y347H
 MO-4 36mm 1/2 Circle	0	Violet	8 × 18"/45cm	12	CR	Y701D*
	0	Violet	4 × 27"/70cm	12	CR	YY71G
	0	Violet	1 × 27"/70cm	24		Y436H
	1	Violet	1 × 27"/70cm	24		Y437H
	1	Violet	4 × 27"/70cm	12	CR	YY70G
 CT 40mm 1/2 Circle	3-0	Violet	1 × 27"/70cm	36		Y350H
	2-0	Violet	1 × 27"/70cm	36		Y351H
	2-0	Violet	1 × 36"/90cm	36	CR	Y357H
	0	Violet	1 × 27"/70cm	36		Y352H
	0	Violet	1 × 36"/90cm	36		Y358H
	0	Violet	1 × 36"/90cm	36	Plus	MCP358H
	0	Violet	8 × 18"/45cm	24	CR	Y752D
	1	Violet	1 × 27"/70cm	36		Y353H
	1	Violet	1 × 36"/90cm	36		Y359H
	 CTX 48mm 1/2 Circle	2-0	Violet	1 × 36"/90cm	36	
0		Violet	1 × 36"/90cm	36		Y398H
1		Violet	1 × 36"/90cm	36		Y399H

*Safety organizer tray



CR CONTROL RELEASE™ Needle Suture

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture




MONOCRYL™ Plus Suture


Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 UR-6 26mm 5/8 Circle	3-0	Violet	1 × 27"/70cm	36		Y604H
	3-0	Violet	1 × 27"/70cm	36	Plus	MCP604H
	2-0	Violet	1 × 27"/70cm	36		Y605H
	2-0	Violet	1 × 27"/70cm	36	Plus	MCP605H
	0	Violet	1 × 27"/70cm	36		Y606H
 UR-5 36mm 5/8 Circle	3-0	Violet	1 × 27"/70cm	36		Y377H
	2-0	Violet	1 × 27"/70cm	36	Plus	MCP377H
	0	Violet	1 × 27"/70cm	36		Y378H


TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-34 36mm 1/2 Circle	2-0	Undyed	1 × 27"/70cm	36		Y817H
	2-0	Undyed	1 × 36"/90cm	36		Y917H

Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-1 36mm 1/2 Circle	2-0	Undyed	1 × 27"/70cm	36		Y266H
	2-0	Undyed	1 × 27"/70cm	36	Plus	MCP266H
	0	Undyed	1 × 27"/70cm	36		Y267H

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture




MONOCRYL™ Plus Suture


Taperpoint Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1, CP-1 36mm 1/2 Circle	3-0	Undyed	1 x 54"/135cm	36		Y931H
	2-0	Undyed	1 x 54"/135cm	36		Y932H


Straight Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Undyed	1 x 27"/70cm	36		Y523H

PRIME Reverse Cutting





Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	Y489G
	6-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	MCP489G
 P-3 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	Y490G
	5-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	MCP490G
	6-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	Y492G
	6-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	MCP492G
	5-0	Undyed	1 x 18"/45cvm	12	MultiPass [®] NEEDLER	Y493G
	5-0	Undyed	1 x 18"/45cm	12	MultiPass [®] NEEDLER	MCP493G
	5-0	Violet	1 x 18"/45cm	12	MultiPass [®] NEEDLER	Y463G

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture




MONOCRYL™ Plus Suture

PRIME Reverse Cutting 









Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-3 13mm 3/8 Circle	5-0	Violet	1 × 18 [”] /45cm	12	Plus MultiPass [®] NEEDLER	MCP463G
	4-0	Undyed	1 × 18 [”] /45cm	12	MultiPass [®] NEEDLER	Y494G
	4-0	Undyed	1 × 18 [”] /45cm	12	Plus MultiPass [®] NEEDLER	MCP494G
	4-0	Violet	1 × 18 [”] /45cm	12	MultiPass [®] NEEDLER	Y464G
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 × 18 [”] /45cm	12	MultiPass [®] NEEDLER	Y495G
	5-0	Undyed	1 × 18 [”] /45cm	12	Plus MultiPass [®] NEEDLER	MCP495G
	4-0	Undyed	1 × 18 [”] /45cm	12	MultiPass [®] NEEDLER	Y496G
	4-0	Undyed	1 × 18 [”] /45cm	12	Plus MultiPass [®] NEEDLER	MCP496G
	4-0	Violet	1 × 18 [”] /45cm	12	MultiPass [®] NEEDLER	Y513G
	4-0	Violet	1 × 18 [”] /45cm	12	Plus MultiPass [®] NEEDLER	MCP513G
	4-0	Undyed	1 × 27 [”] /70cm	36	MultiPass [®] NEEDLER	Y426H
	4-0	Undyed	1 × 27 [”] /70cm	36	Plus MultiPass [®] NEEDLER	MCP426H
	3-0	Undyed	1 × 18 [”] /45cm	12	MultiPass [®] NEEDLER	Y497G
	3-0	Undyed	1 × 18 [”] /45cm	12	Plus MultiPass [®] NEEDLER	MCP497G
	3-0	Undyed	1 × 27 [”] /70cm	36	MultiPass [®] NEEDLER	Y427H
	3-0	Undyed	1 × 27 [”] /70cm	36	Plus MultiPass [®] NEEDLER	MCP427H
	 PS-1 24mm 3/8 Circle	4-0	Undyed	1 × 18 [”] /45cm	36	MultiPass [®] NEEDLER Plus
4-0		Undyed	1 × 18 [”] /45cm	36	MultiPass [®] NEEDLER	MCP682H
4-0		Undyed	1 × 27 [”] /70cm	36	MultiPass [®] NEEDLER	Y935H
4-0		Undyed	1 × 27 [”] /70cm	36	Plus MultiPass [®] NEEDLER	MCP935H
3-0		Undyed	1 × 27 [”] /70cm	36	MultiPass [®] NEEDLER	Y936H
3-0		Undyed	1 × 27 [”] /70cm	36	Plus MultiPass [®] NEEDLER	MCP936H
 PS-6 10mm 1/2 Circle	6-0	Undyed	1 × 18 [”] /45cm	12		Y510G

MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture




MONOCRYL™ Plus Suture

PRIME Conventional Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y833G
	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y834G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	MCP834G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y835G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	MCP835G
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y844G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	MCP844G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y845G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	MCP845G
 PC-5 19mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y822G
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	MCP822G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®]	Y823G
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass [®]	MCP823G
 PC-12 16mm 1/2 Circle	4-0	Undyed	1 × 18"/45cm	12		Y814G

SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Violet	6 × 18"/45cm	12		Y908G*
	2-0	Violet	6 × 18"/45cm	12		Y909G*
	0	Violet	6 × 18"/45cm	12		Y910G*
	1	Violet	6 × 18"/45cm	12		Y911G*



PDS™ Plus Antibacterial (polydioxanone) Suture

PDS™ II (polydioxanone) Suture


PDS II Suture is prepared from polyester, poly (p-dioxanone).¹

PDS II Sutures are indicated for use in soft tissue approximation, including use in pediatric cardiovascular tissue and ophthalmic surgery. PDS II Suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue.¹


These sutures are particularly useful where the combination of an absorbable suture and extended wound support (up to 6 weeks) is desirable.¹




PDS Plus Sutures combine the features of soft, pliable, monofilament construction with absorbability and extended wound support¹:




- Retains 60% of its original strength for 6 weeks, providing support to the fascia as it slowly heals²
- Shown in vitro to inhibit bacterial colonization of the suture³
- Monofilament design and polymer properties minimize drag force and elicit only a slight tissue reaction during absorption¹




PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture







Blunt Point 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>BP-2 48mm 1/2 Circle</p>	0	Violet	1 × 60"/150cm	12	Plus	PDPB990G
 <p>BP-1 65mm 1/2 Circle</p>	0	Violet	1 × 60" Looped 30"	12	Plus Looped	PDPB991G
	1	Violet	1 × 96" Looped 48"	12	Plus Looped	PDPB880G


ETHIGUARD™ Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SHB 26mm 1/2 Circle</p>	3-0	Violet	1 × 27"/70cm	36	Plus	PDPB316H
	2-0	Violet	1 × 27"/70cm	36	Plus	PDPB317H
 <p>CTB-2 26mm 1/2 Circle</p>	3-0	Violet	1 × 27"/70cm	36	Plus	PDPB332H
	2-0	Violet	1 × 27"/70cm	36	Plus	PDPB333H
	0	Violet	1 × 27"/70cm	36	Plus	PDPB334H










PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture




ETHIGUARD™ Blunt Point Needle						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTB-1 36mm 1/2 Circle</p>	2-0	Violet	1 × 36"/90cm	36		ZB345
	2-0	Violet	1 × 36"/90cm	36	Plus	PDPB345
	2-0	Violet	1 × 27"/70cm	36	Plus	PDPB339H
	0	Violet	8 × 18"/45cm	12	CR	ZB740
	0	Violet	1 × 27"/70cm	36		ZB340
	0	Violet	1 × 27"/70cm	36	Plus	PDPB340
	0	Violet	1 × 36"/90cm	36		ZB346
	0	Violet	1 × 36"/90cm	36	Plus	PDPB346
	1	Violet	1 × 36"/90cm	36		ZB347
	1	Violet	1 × 36"/90cm	36	Plus	PDPB347H
 <p>CTB 40mm 1/2 Circle</p>	0	Violet	1 × 36"/90cm	24	Plus	PDPB358T
	1	Violet	1 × 36"/90cm	24	Plus	PDPB359T
 <p>CTXB 48mm 1/2 Circle</p>	0	Violet	1 × 36"/90cm	24		ZB370
	0	Violet	1 × 36"/90cm	24	Plus	PDPB370T
	1	Violet	1 × 36"/90cm	24		ZB371
	1	Violet	1 × 36"/90cm	24	Plus	PDPB371T
 <p>TPB-1 65mm 1/2 Circle</p>	1	Violet	1 × 54"/135cm	12	Plus	PDPB879G





PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	7-0	Violet	1 × 24/60cm	36		Z155H
	7-0	Violet	1 × 30"/75cm	36		Z135H
	6-0	Violet	1 × 30"/75cm	36		Z117H
	6-0	Violet	1 × 30"/75cm	36	Plus	PDP117H
 BV 11mm 3/8 Circle	6-0	Violet	1 × 30"/75cm	36	Plus	PDP123H
	6-0	Violet	1 × 30"/75cm	36	Plus	PDP123H
 C-1 13mm 3/8 Circle	6-0	Violet	1 × 30"/75cm	36		Z127H
	6-0	Violet	1 × 30"/75cm	36	Plus	PDP149H
	5-0	Violet	1 × 30"/75cm	36		Z126H
	5-0	Violet	1 × 30"/75cm	36	Plus	PDP148H
 RB-2 13mm 1/2 Circle	6-0	Violet	1 × 30"/75cm	36		Z149H
	6-0	Violet	1 × 30"/75cm	36	Plus	PDP149H
	5-0	Violet	1 × 30"/75cm	36		Z148H
	5-0	Violet	1 × 30"/75cm	36	Plus	PDP148H
 TF 13mm 1/2 Circle	6-0	Violet	1 × 27"/70cm	36		Z432H
	6-0	Violet	1 × 27"/70cm	36	Plus	PDP432H
 UCL 13mm 1/2 Circle	2-0	Violet	1 × 27"/70cm	36		Z114H
	2-0	Violet	1 × 27"/70cm	36	Plus	PDP114H
 RB-1 17mm 1/2 Circle	5-0	Violet	1 × 27"/70cm	36		Z303H
	5-0	Violet	1 × 27"/70cm	36	Plus	PDP303H
	5-0	Violet	8 × 24"/60cm	12	Plus	PDP800G
	4-0	Violet	1 × 27"/70cm	36		Z304H
	4-0	Violet	1 × 27"/70cm	36	Plus	PDP304H
	3-0	Violet	1 × 27"/70cm	36		Z305H
	3-0	Violet	1 × 27"/70cm	36	Plus	PDP305H
	5-0	Violet	1 × 30"/75cm	36		Z320H
 RB-1 17mm 1/2 Circle	5-0	Violet	1 × 30"/75cm	36	Plus	PDP320H


PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture





Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH-1 22mm 1/2 Circle	4-0	Violet	1 × 27"/70cm	36		Z310H
	4-0	Violet	1 × 27"/70cm	36	Plus	PDP310H
	3-0	Violet	1 × 27"/70cm	36		Z311H
	3-0	Violet	1 × 27"/70cm	36	Plus	PDP311H
 CT-2 26mm 1/2 Circle	3-0	Violet	1 × 27"/70cm	36		Z332H
	3-0	Violet	1 × 27"/70cm	36	Plus	PDP332H
	2-0	Violet	1 × 27"/70cm	36		Z333H
	2-0	Violet	1 × 27"/70cm	36	Plus	PDP333H
	0	Violet	1 × 27"/70cm	36		Z334H
	0	Violet	1 × 27"/70cm	36	Plus	PDP334H
 CTX 48mm 1/2 Circle	0	Violet	8 × 18"/45cm	12	CR	Z776D*
	0	Violet	1 × 36"/90cm	24		Z370T
	0	Violet	1 × 36"/90cm	24	Plus	PDP370T
	0	Violet	1 × 60" Looped 30"	12	Looped	Z990G
	0	Violet	1 × 60" Looped 30"	12	Plus Looped	PDP990G
	1	Violet	8 × 18"/45cm	12	CR	Z765D*
	1	Violet	1 × 36"/90cm	24		Z371T
	1	Violet			Plus	PDP371T
	 TP-1 65mm 1/2 Circle	0	Violet	1 × 60" Looped 30"	12	Looped
0		Violet	1 × 60" Looped 30"	12	Plus	PDP9910G
1		Violet	2 × 27"/70cm	12		Z690G
1		Violet	1 × 54"/135cm	12		Z879G
1		Violet	1 × 96" Looped 48"	12	Looped	Z880G
1		Violet	1 × 96" Looped 48"	12	Plus Looped	PDP880G

*Safety organizer tray

PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture




Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 XLH 70mm 1/2 Circle	1	Violet	1 × 27"/70cm	12	Plus	PDP582G
	1	Violet	1 × 27"/70cm	12		Z583G
	1	Violet	1 × 54"/135cm	12		Z877G
	1	Violet	1 × 54"/135cm	12	Plus	PDP877G
	1	Violet	1 × 96" Looped 48"	12	Looped	Z881G
	1	Violet	1 × 96" Looped 48"	12	Plus Looped	PDP881G
 ST-4 19mm Straight	4-0	Violet	1 × 14"/35cm	12		Z420G
	4-0	Violet	1 × 14"/35cm	12	Plus	PDP420G




TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-5 17mm 1/2 Circle	4-0	Violet	1 × 36"/90cm	36		Z137H
	4-0	Violet	1 × 36"/90cm	36	Plus	PDP137H
 V-38 30mm 1/2 Circle	0	Violet	1 × 18"/45cm	12		Z00JG


PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture








Trocar						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
STP-10 254mm Straight	2-0	Violet	1 × 27"/70cm	12		Z997G

Reverse Cutting ▽						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	5-0	Undyed	1 × 27"/70cm	36		Z421H
	5-0	Undyed	1 × 27"/70cm	36	Plus	PDP421H
	4-0	Undyed	1 × 27"/70cm	36		Z422H
	4-0	Undyed	1 × 27"/70cm	36	Plus	PDP422H
	3-0	Undyed	1 × 27"/70cm	36		Z423H
	3-0	Undyed	1 × 27"/70cm	36	Plus	PDP423H
 FS-1 24mm 3/8 Circle	4-0	Undyed	1 × 27"/70cm	36		Z441H
	4-0	Undyed	1 × 27"/70cm	36	Plus	PDP441H
	3-0	Undyed	1 × 27"/70cm	36		Z442H
	3-0	Undyed	1 × 27"/70cm	36	Plus	PDP442H
	2-0	Undyed	1 × 27"/70cm	36		Z443H
	2-0	Undyed	1 × 27"/70cm	36	Plus	PDP443H
 OS-4 22mm 1/2 Circle	1	Violet	1 × 27"/70cm	36		Z695H
	1	Violet	1 × 27"/70cm	36	Plus	PDP695H

PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture




Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CP-2 26mm 1/2 Circle</p>	3-0	Violet	8 × 18"/45cm	12	GR	Z761D*
	3-0	Violet	1 × 27"/70cm	36		Z968H
	3-0	Violet	1 × 27"/70cm	36	Plus	PDP968H
	2-0	Violet	8 × 18"/45cm	12	GR	Z762D*
	2-0	Violet	1 × 27"/70cm	36		Z969H
	2-0	Violet	1 × 27"/70cm	36	Plus	PDP969H
	0	Violet	1 × 27"/70cm	36		Z970H
	0	Violet	1 × 27"/70cm	36	Plus	PDP970H
 <p>CP-1 36mm 1/2 Circle</p>	2-0	Violet	1 × 27"/70cm	36		Z466H
	2-0	Violet	1 × 27"/70cm	36	Plus	PDP466H
	0	Violet	1 × 27"/70cm	36		Z467H
	0	Violet	1 × 27"/70cm	36	Plus	PDP467H
	1	Violet	1 × 27"/70cm	36		Z468H
	1	Violet	1 × 27"/70cm	36	Plus	PDP468H
	 <p>OS-6 36mm 1/2 Circle</p>	0	Violet	1 × 27"/70cm	24	
0		Violet	1 × 27"/70cm	24	Plus	PDP534T
1		Violet	1 × 27"/70cm	36		Z535H
1		Violet	1 × 27"/70cm	36	Plus	PDP535H
 <p>CP 40mm 1/2 Circle</p>	2	Violet	1 × 27"/70cm	24		Z195T
	 <p>CPX 48mm 1/2 Circle</p>	1	Violet	1 × 27"/70cm	24	
1		Violet	1 × 27"/70cm	24	Plus	PDP569T


*Safety organizer tray

GR CONTROL RELEASE™ Needle Suture






















PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture




Straight Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Undyed	1 × 27"/70cm	36		Z662H
	4-0	Undyed	1 × 27"/70cm	36	Plus	PDP662H
	3-0	Undyed	1 × 27"/70cm	36		Z663H
	3-0	Undyed	1 × 27"/70cm	36	Plus	PDP663H


PRIME Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 P-1 11mm 3/8 Circle	6-0	Undyed	1 × 18"/45cm	12	 MultiPass	Z489G	
	6-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass	PDP489G	
	6-0	Violet	1 × 18"/45cm	12	 MultiPass	Z487G	
	6-0	Violet	1 × 18"/45cm	12	Plus  MultiPass	PDP487G	
 P-3 13mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass	Z493G	
	5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass	PDP493G	
	5-0	Violet	1 × 18"/45cm	12	 MultiPass	Z463G	
	5-0	Violet	1 × 18"/45cm	12	Plus  MultiPass	PDP463G	
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass	Z494G	
	4-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass	PDP494G	
	 PS-2 19mm 3/8 Circle	5-0	Undyed	1 × 18"/45cm	12	 MultiPass	Z495G
		5-0	Undyed	1 × 18"/45cm	12	Plus  MultiPass	PDP495G
4-0		Undyed	1 × 18"/45cm	12	 MultiPass	Z496G	
4-0		Undyed	1 × 18"/45cm	12	Plus  MultiPass	PDP496ZG	
4-0		Violet	1 × 18"/45cm	12	 MultiPass	Z513G	
4-0		Violet	1 × 18"/45cm	12	Plus  MultiPass	PDP513G	
3-0		Undyed	1 × 18"/45cm	12	 MultiPass	Z497G	
3-0		Undyed	1 × 18"/45cm	12	Plus  MultiPass	PDP497G	


PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture




PRIME Reverse Cutting 















Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 × 18”/45cm	12	 MultiPass	Z682G
	4-0	Undyed	1 × 18”/45cm	12	 Plus	PDP682G
	3-0	Undyed	1 × 18”/45cm	12	 MultiPass	Z683G
	3-0	Undyed	1 × 18”/45cm	12	 Plus	PDP683G
 P-2 8mm 1/2 Circle	5-0	Undyed	1 × 18”/45cm	12	 MultiPass	Z503G
	5-0	Undyed	1 × 18”/45cm	12	 Plus	PDP503G
	4-0	Undyed	1 × 18”/45cm	12	 MultiPass	Z504G
	4-0	Undyed	1 × 18”/45cm	12	 Plus	PDP504G
 PS-4 16mm 1/2 Circle	4-0	Undyed	1 × 18”/45cm	12	 MultiPass	Z507G
	4-0	Undyed	1 × 18”/45cm	12	 Plus	PDP507G


Conventional Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CCS-1 60mm 1/2 Circle	2	Violet	1 × 27”/70cm	24		Z807T

PDS™ Plus Antibacterial (polydioxanone) Suture PDS™ II (polydioxanone) Suture




PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Undyed	1 × 18 ^{3/8} /45cm	12	 MultiPass [®] NEEDLES	Z833G
	6-0	Undyed	1 × 18 ^{3/8} /45cm	12	 Plus MultiPass [®] NEEDLES	PDP833G
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 × 18 ^{3/8} /45cm	12	 MultiPass [®] NEEDLES	Z844G
	5-0	Undyed	1 × 18 ^{3/8} /45cm	12	 Plus MultiPass [®] NEEDLES	PDP844G
 PC-5 19mm 3/8 Circle	4-0	Undyed	1 × 18 ^{3/8} /45cm	12	 MultiPass [®] NEEDLES	Z845G
	4-0	Undyed	1 × 18 ^{3/8} /45cm	12	 Plus MultiPass [®] NEEDLES	PDP845G
	4-0	Undyed	1 × 18 ^{3/8} /45cm	12	 MultiPass [®] NEEDLES	Z823G
	4-0	Undyed	1 × 18 ^{3/8} /45cm	12	 Plus MultiPass [®] NEEDLES	PDP823G
	3-0	Undyed	1 × 18 ^{3/8} /45cm	12	 MultiPass [®] NEEDLES	Z824G
	3-0	Undyed	1 × 18 ^{3/8} /45cm	12	 Plus MultiPass [®] NEEDLES	PDP824G




Surgical Gut Suture – Chromic





Surgical Gut Suture is an absorbable, sterile surgical suture composed of purified connective tissue (mostly collagen) derived from either the serosal layer of beef (bovine) or the submucosal fibrous layer of sheep (ovine) intestines. Surgical Gut Sutures are available in plain or chromic.¹

Surgical Gut Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.¹




Surgical Gut Suture – Chromic




Blunt Point <input checked="" type="radio"/>							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 <p>BP-1 65mm 1/2 Circle</p>	0	Brown	1 × 27"/70cm	24		47T	
 <p>BP 80mm 1/2 Circle</p>	1	Brown	1 × 27"/70cm	12		48G	
Taperpoint <input type="radio"/>							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 <p>C-1 13mm 3/8 Circle</p>	5-0	Brown	1 × 27"/70cm	36		K895H	
 <p>RB-1 17mm 1/2 Circle</p>	5-0	Brown	1 × 27"/70cm	36		U202H	
	4-0	Brown	1 × 27"/70cm	36		U203H	
	3-0	Brown	1 × 27"/70cm	36		U204H	

Surgical Gut Suture – Chromic






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-3 22mm 1/2 Circle	3-0	Brown	1 × 27"/70cm	36		892H
	2-0	Brown	1 × 27"/70cm	36		893H
 SH-1 22mm 1/2 Circle	4-0	Brown	1 × 27"/70cm	36		G181H
	3-0	Brown	1 × 27"/70cm	36		G182H
 CT-2 26mm 1/2 Circle	3-0	Brown	1 × 27"/70cm	36		882H
	2-0	Brown	1 × 27"/70cm	36		883H
	0	Brown	1 × 27"/70cm	36		884H
	1	Brown	1 × 27"/70cm	36		885H
 SH 26mm 1/2 Circle	4-0	Brown	1 × 27"/70cm	36		G121H
	3-0	Brown	1 × 27"/70cm	36		G122H
	2-0	Brown	1 × 27"/70cm	36		G123H
	0	Brown	1 × 27"/70cm	36		G124H
 SH 26mm 1/2 Circle	3-0	Brown	1 × 36"/90cm	36		G172H
	2-0	Brown	1 × 36"/90cm	36		G173H
 MO-5 31mm 1/2 Circle	0	Brown	8 × 18"/45cm	12	CR	CC04G
	1	Brown	8 × 18"/45cm	12	CR	CC03G
	1	Brown	1 × 27"/70cm	36		4995H


Surgical Gut Suture – Chromic








Taperpoint


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	3-0	Brown	1 × 27"/70cm	36		810H
	3-0	Brown	1 × 36"/90cm	36		922H
	2-0	Brown	1 × 27"/70cm	36		811H
	2-0	Brown	3 × 27"/70cm	24		3811T
	2-0	Brown	1 × 36"/90cm	36		923H
	0	Brown	8 × 18"/45cm	12	CR	CC41G
	0	Brown	1 × 27"/70cm	36		812H
	0	Brown	3 × 27"/70cm	24		3812T
	O	Brown	1 × 36"/90cm	36		924H
	1	Brown	8 × 18"/45cm	12	CR	CC40G
	1	Brown	1 × 27"/70cm	36		813H
	1	Brown	3 × 27"/70cm	12		3813G
	1	Brown	1 × 36"/90cm	36		925H
	2	Brown	1 × 27"/70cm	36		814H
 <p>MH 36mm 1/2 Circle</p>	3-0	Brown	1 × 27"/70cm	36		G126H
	2-0	Brown	1 × 27"/70cm	36		G127H
	0	Brown	1 × 27"/70cm	36		G128H
 <p>MO-4 36mm 1/2 Circle</p>	0	Brown	8 × 18"/45cm	12	CR	CC02G
	O	Brown	1 × 27"/70cm	36		4924H
	1	Brown	8 × 18"/45cm	12	CR	CC01G
	1	Brown	1 × 27"/70cm	36		4925H
	1	Brown	3 × 27"/70cm	24		5589T
	1	Brown	8 × 27"/70cm	12	CR	CC38G
	2-0	Brown	1 × 27"/70cm	24		4923H







Surgical Gut Suture – Chromic




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT 40mm 1/2 Circle</p>	3-0	Brown	1 × 27"/70cm	36		800H
	3-0	Brown	1 × 36"/90cm	36		912H
	2-0	Brown	1 × 27"/70cm	36		801H
	2-0	Brown	1 × 36"/90cm	36		913H
	0	Brown	3 × 27"/70cm	24		3802T
	0	Brown	1 × 27"/70cm	36		802H
	0	Brown	1 × 36"/90cm	36		914H
	1	Brown	1 × 27"/70cm	36		803H
	1	Brown	3 × 27"/70cm	24		3803T
	1	Brown	1 × 36"/90cm	36		915H
 <p>CT 40mm 1/2 Circle</p>	3-0	Brown	1 × 54"/135cm	36		990H
	2-0	Brown	1 × 54"/135cm	36		991H
 <p>LH 40mm 1/2 Circle</p>	2-0	Brown	1 × 27"/70cm	36		G133H
 <p>MO-2 40mm 1/2 Circle</p>	1	Brown	1 × 27"/70cm	36		4915H



Surgical Gut Suture – Chromic



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	3-0	Brown	1 × 27"/70cm	36		862H
	2-0	Brown	1 × 27"/70cm	36		863H
	2-0	Brown	1 × 36"/90cm	36		903H
	0	Brown	1 × 27"/70cm	36		864H
	0	Brown	1 × 36"/90cm	36		904H
	1	Brown	1 × 27"/70cm	36		865H
	1	Brown	1 × 36"/90cm	36		905H
 <p>TP-1 65mm 1/2 Circle</p>	1	Brown	4 × 27"/70cm	12		M760G
	1	Brown	1 × 60" Looped 30"	12	Looped	GL30G
	2	Brown	1 × 60" Looped 30"	12	Looped	GL31G
 <p>XLH 70mm 1/2 Circle</p>	0	Brown	1 × 27"/70cm	24		44T
 <p>UR-6 26mm 5/8 Circle</p>	3-0	Brown	1 × 27"/70cm	36		N877H
	2-0	Brown	1 × 27"/70cm	36		N878H
	0	Brown	1 × 27"/70cm	36		N879H
 <p>UR-5 36mm 5/8 Circle</p>	2-0	Brown	1 × 27"/70cm	36		U245H
	1	Brown	1 X 27"/70cm	24		U247H
	0	Brown	1 × 27"/70cm	36		U246H


Surgical Gut Suture – Chromic









Taperpoint 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 <p>UR-4 40mm 5/8 Circle</p>	2-0	Brown	1 × 27"/70cm	36		U255H	
	0	Brown	1 × 27"/70cm	36		U256H	
 <p>SH 26mm/ ST-1 45mm 1/2 Circle/ Straight</p>	3-0	Brown	1 × 36"/90cm	36		G152H	


TAPERCUT™ 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 <p>V-5 17mm 1/2 Circle</p>	4-0	Brown	1 × 27"/70cm	36		U233H	
	3-0	Brown	1 × 36"/90cm	36		942H	
	2-0	Brown	1 × 27"/70cm	36		837H	
	2-0	Brown	1 × 36"/90cm	36		943H	
 <p>V-34 36mm 1/2 Circle</p>	0	Brown	1 × 27"/70cm	36		838H	
	0	Brown	1 × 36"/90cm	36		944H	
	1	Brown	1 × 27"/70cm	36		839H	
	1	Brown	1 × 36"/90cm	36		945H	



Surgical Gut Suture – Chromic






Reverse Cutting ▾							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 M-2 13mm 1/4 Circle	4-0	Brown	1 × 12"/30cm	12		752G	
 M-2 13mm 1/4 Circle	4-0	Brown	1 × 12"/30cm	12		751G	
 C-3 13mm 3/8 Circle	5-0	Brown	1 × 12"/30cm	12		723G	
 M-1 13mm 1/2 Circle	4-0	Brown	1 × 12"/30cm	12		744G	
 J-1 17mm 1/2 Circle	4-0	Brown	1 × 18"/45cm	12		724G	
 CP-2 26mm 1/2 Circle	4-0	Brown	1 × 27"/70cm	36		U213H	
	3-0	Brown	1 × 27"/70cm	36		886H	
	2-0	Brown	1 × 27"/70cm	36		887H	
	0	Brown	1 × 27"/70cm	24		888H	

Surgical Gut Suture – Chromic




Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-1 36mm 1/2 Circle	2-0	Brown	1 × 27"/70cm	36		816H
	0	Brown	1 × 27"/70cm	36		817H
	1	Brown	1 × 36"/90cm	36		975H
 CP 40mm 1/2 Circle	3-0	Brown	1 × 27"/70cm	36		832H
	2-0	Brown	1 × 27"/70cm	36		833H
	1	Brown	1 × 27"/70cm	36		805H

Taperpoint ⊙ , Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1, CP-1 36mm 1/2 Circle	3-0	Brown	1 × 54"/135cm	36		982H
	2-0	Brown	1 × 54"/135cm	36		983H
 CT, FSL 40mm 1/2 Circle	2-0	Brown	1 × 54"/135cm	36		977H

Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Brown	1 × 27"/70cm	36		654H


Surgical Gut Suture – Chromic




PRIME Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-6 8mm 1/4 Circle	6-0	Brown	1 × 18"/45cm	12		794G
 G-1 11mm 3/8 Circle	6-0	Blue	1 × 18"/45cm	12		795G
 G-1 11mm 3/8 Circle	6-0	Brown	1 × 18"/45cm	12		796G
 G-1 11mm 3/8 Circle	6-0	Brown	1 × 18"/45cm	12		790G
 G-3 13mm 3/8 Circle	5-0	Brown	1 × 18"/45cm	12		792G
 G-3 13mm 3/8 Circle	4-0	Brown	1 × 18"/45cm	12		793G


Surgical Gut Suture – Chromic








PRIME Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>G-2 8mm 1/2 Circle</p>	4-0	Brown	1 × 18"/45cm	12		798G


Surgical Gut Suture – Chromic




Spatula Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG100-8 6.5mm 1/4 Circle	7-0	Brown	1 × 18"/45cm	12		1744G
 S-14 8mm 1/4 Circle	6-0	Blue	1 × 18"/45cm	12		1791G
	5-0	Brown	1 × 12"/30cm	12		1766G
	5-0	Brown	1 × 18"/45cm	12		1792G
 S-4 8mm 1/4 Circle	4-0	Brown	1 × 18"/45cm	12		1793G
 TG140-8 65mm 3/8 Circle	7-0	Blue	1 × 18"/45cm	12		1797G
	6-0	Blue	1 × 18"/45cm	12		1731G
 S-2 8mm 1/2 Circle	4-0	Brown	1 × 18"/45cm	12		1798G

ULTIMA™ Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CS175-8 7mm 1/2 Circle	7-0	Brown	1 × 18"/45cm	12		1745G

Surgical Gut Suture – Chromic



LIGAPAK™ Dispensing Reel


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Brown	1 × 54"/135cm	12		L111G
	3-0	Brown	1 × 54"/135cm	12		L112G
	2-0	Brown	1 × 54"/135cm	12		L113G
	0	Brown	1 × 54"/135cm	12		L114G
	1	Brown	1 × 54"/135cm	12		L115G

Pre-Cut Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Brown	1 × 54"/135cm	36		S111H
	3-0	Brown	1 × 54"/135cm	36		S112H
	2-0	Brown	1 × 54"/135cm	36		S113H
	0	Brown	1 × 54"/135cm	36		S114H
	1	Brown	1 × 54"/135cm	36		S115H
	2	Brown	1 × 54"/135cm	36		S116H
	3	Brown	1 × 54"/135cm	36		S117H

SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	2-0	Brown	12 × 18"/45cm	24		SG13T
	0	Brown	12 × 18"/45cm	24		SG14T
	1	Brown	12 × 18"/45cm	24		SG15T




Surgical Gut Suture – Plain









Surgical Gut Suture is an absorbable, sterile surgical suture composed of purified connective tissue (mostly collagen) derived from either the serosal layer of beef (bovine) or the submucosal fibrous layer of sheep (ovine) intestines. Surgical Gut Sutures are available in plain or chromic.¹

Surgical Gut Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.¹










Surgical Gut Suture – Plain




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	4-0	Yellowish Tan	1 × 27”/70cm	36		U207H
 CT-3 22mm 1/2 Circle	3-0	Yellowish Tan	1 × 27”/70cm	36		N862H
	2-0	Yellowish Tan	1 × 27”/70cm	36		N863H
	0	Yellowish Tan	1 × 27”/70cm	36		N864H
 SH 26mm 1/2 Circle	4-0	Yellowish Tan	1 × 27”/70cm	36		G321H
	3-0	Yellowish Tan	1 × 27”/70cm	36		G322H
 CT-1 36mm 1/2 Circle	3-0	Yellowish Tan	3 × 18”/45cm	24		5503T
	3-0	Yellowish Tan	1 × 27”/70cm	36		842H
	2-0	Yellowish Tan	1 × 27”/70cm	36		843H
	0	Yellowish Tan	1 × 27”/70cm	36		844H
 CT 40mm 1/2 Circle	3-0	Yellowish Tan	8 × 18”/45cm	12	CR	CP53G
	3-0	Yellowish Tan	1 × 27”/70cm	36		852H
	2-0	Yellowish Tan	1 × 27”/70cm	36		853H
 CTX 48mm 1/2 Circle	3-0	Yellowish Tan	1 × 27”/70cm	36		873H
 XLH 70mm 1/2 Circle	2-0	Yellowish Tan	1 × 27”/70cm	36		872H
	3-0	Yellowish Tan	1 × 27”/70cm	24		52T
	2-0	Yellowish Tan	1 × 27”/70cm	24		53T

Surgical Gut Suture – Plain



TAPERCUT™ 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-7 26mm 1/2 Circle</p>	3-0	Yellowish Tan	1 × 18"/45cm	12		849G
 <p>V-34 36mm 1/2 Circle</p>	3-0	Yellowish Tan	1 × 27"/70cm	36		846H
Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>M-1 13mm 1/2 Circle</p>	4-0	Yellowish Tan	1 × 12"/30cm	12		742G
 <p>X-1 22mm 1/2 Circle</p>	3-0	Brown	1 × 18"/45cm	12		612G
 <p>CP-1 36mm 1/2 Circle</p>	2-0	Yellowish Tan	1 × 27"/70cm	36		898H




Surgical Gut Suture – Plain



Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-1 11mm 3/8 Circle	6-0	Yellowish Tan	1 × 18"/45cm	12		774G
 G-3 13mm 3/8 Circle	5-0	Yellowish Tan	1 X 18"/45cm	12		772G
	4-0	Yellowish Tan	1 X 18"/45cm	12		773G


Straight Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SC-1 13mm Straight	4-0	Yellowish Tan	1 × 18"/45cm	36		1824H
 SC-1 13mm Straight	4-0	Yellowish Tan	1 × 18"/45cm	36		1828H
 KS 60mm Straight	3-0	Yellowish Tan	1 × 27"/70cm	36		665H

PRIME Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-6 8mm 1/4 Circle	6-0	Blue	1 × 18"/45cm	12		775G
 G-1 11mm 3/8 Circle	6-0	Yellowish Tan	1 × 18"/45cm	12		770G


Surgical Gut Suture – Plain



PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-2 19mm 3/8 Circle	4-0	Yellowish Tan	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	1627G

PRIME Conventional Cutting  Fast-Absorbing Sutures						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	5-0	Yellowish Tan	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1915G

Surgical Gut Suture – Plain



Spatula Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 S-14 8mm 1/4 Circle	5-0	Yellowish Tan	1 × 12"/30cm	12		1775G
	5-0	Yellowish Tan	1 × 18"/45cm	12		1772G
 TG140-8 6.5mm 3/8 Circle	6-0	Yellowish Tan	1 × 18"/45cm	12		1735G


LIGAPAK™ Dispensing Reel

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Yellowish Tan	1 × 54"/135cm	12		L101G
	3-0	Yellowish Tan	1 × 54"/135cm	12		L102G
	2-0	Yellowish Tan	1 × 54"/135cm	12		L103G
	0	Yellowish Tan	1 × 54"/135cm	12		L104G

Pre-Cut Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Yellowish Tan	12 × 18"/45cm	24		SG2T
	3-0	Yellowish Tan	1 × 54"/135cm	36		S102H
	2-0	Yellowish Tan	1 × 54"/135cm	36		S103H
	0	Yellowish Tan	1 × 54"/135cm	36		S104H
	1	Yellowish Tan	1 × 54"/135cm	36		S105H

Sutures: Nonabsorbable




PRONOVA™ Poly (Hexafluoropropylene - VDF) Suture














PRONOVA Suture (pigmented) is a nonabsorbable, sterile surgical suture made from polymer blend of poly (vinylidene fluoride) and poly (vinylidene fluoride-co-hexafluoropropylene). The suture is pigmented blue to enhance visibility.¹

PRONOVA Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular and neurological procedures.¹




PRONOVA™ Poly (Hexafluoropropylene - VDF) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-5 6.5mm 3/8 Circle	8-0	Blue	2 × 24"/60cm	12	 MultiPass™ NEEDLES	M3739*
 BV175-6 8mm 3/8 Circle	8-0	Blue	2 × 24"/60cm	12	 MultiPass™ NEEDLES	M3743*
	7-0	Blue	4 × 24"/60cm	12	 MultiPass™ NEEDLES	M3735*
	7-0	Blue	2 × 24"/60cm	12	 MultiPass™ NEEDLES	M3737*
	7-0	Blue	1 × 24"/60cm	36	 MultiPass™ NEEDLES	PN3735H
 BV175-7 8mm 3/8 Circle	8-0	Blue	2 × 24"/60cm	12	 MultiPass™ NEEDLES	M3767*
	7-0	Blue	4 × 30"/75cm	12	 MultiPass™ NEEDLES	M3770*
	7-0	Blue	2 × 24"/60cm	12	 MultiPass™ NEEDLES	M3766*
	7-0	Blue	2 × 18"/45cm	12	 MultiPass™ NEEDLES	M3768*
 BV175-8 9.3mm 3/8 Circle	7-0	Blue	2 × 24"/60cm	12	 MultiPass™ NEEDLES	M3747*

*Straight pack


PRONOVA™ Poly (Hexafluoropropylene - VDF) Suture








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	7-0	Blue	4 × 30"/75cm	12		PNM3703*
	7-0	Blue	1 × 30"/75cm	36		PN3703H
	7-0	Blue	4 × 24"/60cm	12		PNM3702*
	7-0	Blue	1 × 24"/60cm	36		PN3702H
	7-0	Blue	1 × 24"/60cm	36	VB	3304H
	7-0	Blue	1 × 18"/45cm	36		PN3701H
	6-0	Blue	4 × 30"/75cm	12		M3709*
	6-0	Blue	1 × 30"/75cm	36		PN3709H
	6-0	Blue	4 × 24"/60cm	12		PNM3805*
	6-0	Blue	1 × 24"/60cm	36		PN3805H
	6-0	Blue	1 × 24"/60cm	36	VB	PN3305H
6-0	Blue	1 × 18"/45cm	36		PN3806H	
 BV 11mm 3/8 Circle	6-0	Blue	1 × 30"/75cm	36		3776H
	6-0	Blue	4 × 30"/75cm	12		PNM3706*
 C-1 13mm 3/8 Circle	6-0	Blue	1 × 30"/75cm	36		PN3706H
	6-0	Blue	1 × 30"/75cm	36	VB	3307H
	6-0	Blue	1 × 24"/60cm	36		PN3726H
	5-0	Blue	4 × 36"/90cm	12		M3720*
	5-0	Blue	1 × 36"/90cm	36		PN3720H
	5-0	Blue	1 × 36"/90cm	36	VB	3321H
	5-0	Blue	4 × 24"/60cm	12		PNM3725*
	5-0	Blue	1 × 24"/60cm	36		PN3725H
	5-0	Blue	1 × 18"/45cm	36		PN3717H


*Straight pack

PRONOVA™ Poly (Hexafluoropropylene - VDF) Suture





Taperpoint 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BB 17mm 3/8 Circle	5-0	Blue	1 × 36"/90cm	36		3580H
	4-0	Blue	1 × 36"/90cm	36		PN3581H
 RB-2 13mm 1/2 Circle	6-0	Blue	1 × 30"/75cm	36		PN3711H
	5-0	Blue	1 × 30"/75cm	36		PN3710H
 RB-1 17mm 1/2 Circle	5-0	Blue	1 × 36"/90cm	36		PN3556H
	5-0	Blue	1 × 24"/60cm	36		PN3555H
	4-0	Blue	1 × 36"/90cm	36		PN3557H
	3-0	Blue	1 × 36"/90cm	36		PN3558H
 SH-1 22mm 1/2 Circle	4-0	Blue	1 × 30"/75cm	36		3526H
	3-0	Blue	1 × 30"/75cm	36		3762H
 SH 26mm 1/2 Circle	4-0	Blue	1 × 36"/90cm	36		3521H
	3-0	Blue	1 × 36"/90cm	36		PN3522H
	2-0	Blue	1 × 36"/90cm	36		PN3523H

PRONOVA™ Poly (Hexafluoropropylene - VDF) Suture



TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 MH 36mm 1/2 Circle	3-0	Blue	1 × 36"/90cm	36		3842H
 CC-1 13mm 3/8 Circle	6-0	Blue	1 × 30"/75cm	36		3707H
 V-7 26mm 1/2 Circle	4-0	Blue	1 × 36"/90cm	36		3975H
	3-0	Blue	1 × 36"/90cm	36		PN3976



ETHIBOND EXCEL™

Polyester Suture

ETHIBOND EXCEL Suture is a nonabsorbable, braided, sterile, surgical suture composed of Poly (ethylene terephthalate).¹


The suture is uniformly coated with polybutylate or poly {oxy-1, 4 butanedioxy (1, 6-dioxo-1, 6 hexanedioyl)}. The highly adherent coating is a relatively nonreactive nonabsorbable compound which acts as a lubricant to mechanically improve the physical properties of the uncoated suture by improving handling qualities.¹







ETHIBOND EXCEL Suture is braided for optimal handling properties, and dyed green for enhanced visibility in the surgical field.¹

ETHIBOND EXCEL Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹



ETHIBOND EXCEL™ Polyester Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-1 13mm 3/8 Circle	6-0	Green	1 × 30"/75cm	36		X889H
	5-0	Green	1 × 30"/75cm	36		X890H
 C-1 13mm 3/8 Circle	6-0	Green	1 × 18"/45cm	36		X706H
	5-0	Green	1 × 18"/45cm	36		X717H
	5-0	Green	1 × 30"/75cm	36		X727H
	4-0	Green	1 × 30"/75cm	36	Pledget*	X726H
 BB 17mm 3/8 Circle	4-0	Green	1 × 30"/75cm	36		X881H
 BB 17mm 3/8 Circle	4-0	Green	1 × 30"/75cm	36		X581H
	3-0	Green	1 × 36"/90cm	36		X588H
 TF 13mm 1/2 Circle	5-0	Green	1 × 24"/60cm	36		X205H
	4-0	Green	1 × 24"/60cm	36		X204H
 UCL 13mm 1/2 Circle	2-0	Green	1 × 27"/70cm	36		X114H

*7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

ETHIBOND EXCEL™ Polyester Suture



Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	5-0	Green	1 × 30"/75cm	36		X550H
	5-0	Green	1 × 36"/90cm	36		X556H
	4-0	Green	1 × 30"/75cm	36		X551H
	4-0	Green	4 × 30"/75cm	12		MX551
	4-0	Green	1 × 36"/90cm	36		X557H
	3-0	Green	4 × 30"/75cm	12		MX552
	3-0	Green	1 × 36"/90cm	36		X558H
	2-0	Green	4 × 30"/75cm	12		MX553
	2-0	White	4 × 30"/75cm	12		MX523
	2-0	Grn/Wht	8 × 30"/75cm	6	Pledget*	PXX41N
	2-0	Grn/Wht	8 × 30"/75cm	6	Pledget†	PXX43N
	2-0	Grn/Wht	10 × 30"/75cm	6		10X42N
2-0	Grn/Wht	10 × 30"/75cm	6	Pledget‡	SXX43N	
 RB-1 17mm 1/2 Circle	5-0	Green	1 × 30"/75cm	36		X870H
	4-0	Green	1 × 30"/75cm	36		X871H
	3-0	Green	1 × 30"/75cm	36		X872H
	2-0	Green	1 × 30"/75cm	36		X873H
 SH-2 20mm 1/2 Circle	2-0	Green	1 × 30"/75cm	36	Pledget*	PX82H
	2-0	Green	1 × 30"/75cm	36		X582H
	2-0	White	1 × 30"/75cm	36	Pledget*	PX83H
	2-0	Grn/Wht	10 × 30"/75cm	6	Pledget*	PXX80N
	2-0	Grn/Wht	10 × 30"/75cm	6		PXX82N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget§	SXP82N
	2-0	Green	1 × 36"/90cm	36		X559H
 MO-7 22mm 1/2 Circle	2-0	Green	8 × 18"/45cm	12		CX42D
	0	Green	8 × 18"/45cm	12		CX41D

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget

†3mm x 3mm x 1.5mm Firm TFE Polymer Pledget

‡3mm x 3mm x 1.5mm Soft TFE Polymer Pledget

§7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

^{||}Safety organizer tray

ETHIBOND EXCEL™ Polyester Suture






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SH-1 22mm 1/2 Circle</p>	4-0	Green	1 × 36"/90cm	36		X761H
	3-0	Green	1 × 36"/90cm	36		X762H
	2-0	Green	4 × 30"/75cm	12		MX763
	2-0	Grn/Wht	10 × 30"/75cm	6	Pledget*	PXX92SCN
	2-0	Green	1 × 36"/90cm	36		X763H
 <p>CT-2 26mm 1/2 Circle</p>	2-0	Green	8 × 18"/45cm	12	CR	CX26D [†]
	2-0	Green	1 × 30"/75cm	36		X411H
	0	Green	8 × 18"/45cm	12	CR	CX27D [†]
	0	Green	1 × 30"/75cm	36		X412H
 <p>MO-6 26mm 1/2 Circle</p>	2-0	Green	8 × 18"/45cm	12	CR	CX46D [†]
	0	Green	8 × 18"/45cm	12	CR	CX45D [†]
 <p>SH 26mm 1/2 Circle</p>	4-0	Green	8 × 18"/45cm	12	CR	CX14D [†]
	4-0	Green	1 × 30"/75cm	36		X831H
	3-0	Green	8 × 18"/45cm	12	CR	CX13D [†]
	3-0	Green	1 × 30"/75cm	36		X832H
	2-0	Green	8 × 18"/45cm	12	CR	CX12D [†]
	2-0	Green	1 × 30"/75cm	36		X833H
	2-0	Green	4 × 30"/75cm	12		MX833
	2-0	White	4 × 30"/75cm	12		MX823
	0	Green	1 × 30"/75cm	36		X834H

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
[†]Safety organizer tray

ETHIBOND EXCEL™ Polyester Suture



Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SH 26mm 1/2 Circle</p>	4-0	Green	1 × 36"/90cm	36		X521H
	3-0	Green	1 × 30"/75cm	36		X562H
	3-0	Green	1 × 36"/90cm	36		X522H
	2-0	Green	1 × 30"/75cm	36		X563H
	2-0	Green	4 × 30"/75cm	12		MX563
	2-0	Grn/Wht	10 × 30"/75cm	6		10X82N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget*	PXX76N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget†	SXP76N
	2-0	Grn/Wht	10 × 30"/75cm	6	Pledget‡	SXX86N
	2-0	Green	1 × 36"/90cm	12	Pledget*	PX514G
	2-0	Green	1 × 36"/90cm	36	Pledget†	SXP02H
	2-0	Green	1 × 36"/90cm	36		X523H
	2-0	White	1 × 36"/90cm	36		X513H
	2-0	Grn/Wht	6 × 36"/90cm	6	Pledget*	PXX85N
	2-0	Grn/Wht	10 × 36"/90cm	6		10X86N
	2-0	Grn/Wht	10 × 36"/90cm	12	Pledget*	PXX86
	2-0	Grn/Wht	10 × 36"/90cm	12	Pledget†	SXP86N
0	Grn/Wht	10 × 30"/75cm	12	Pledget*	PXX87	
0	Green	1 × 36"/90cm	36		X524H	
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Green	8 × 18"/45cm	12	CR	CX22D§
	2-0	Green	1 × 30"/75cm	36		X423H
	0	Green	8 × 18"/45cm	12	CR	CX21D§
	0	Green	1 × 30"/75cm	36		X424H
	1	Green	1 × 30"/75cm	36		X425H
 <p>MH 36mm 1/2 Circle</p>	2-0	Green	1 × 36"/90cm	36		X843H
	0	Green	1 × 36"/90cm	36		X844H

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget

†7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

‡3mm x 3mm x 1.5mm Firm TFE Polymer Pledget

§Safety organizer tray

CR CONTROL RELEASE™ Needle Suture

ETHIBOND EXCEL™ Polyester Suture



Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	0	Green	8 × 18"/45cm	12	CR	CX31D*
	0	Green	1 × 30"/75cm	36		X864H
	1	Green	8 × 18"/45cm	12	CR	CX30D*
	1	Green	1 × 30"/75cm	36		X865H
 <p>XLH 70mm 1/2 Circle</p>	1	Green	1 × 30"/75cm	12		X190G

TAPERCUT™






Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CC-1 13mm 3/8 Circle</p>	5-0	Green	1 × 30"/75cm	36		X728H
 <p>V-4 17mm 3/8 Circle</p>	4-0	Green	1 × 30"/75cm	36		X942H

*Safety organizer tray

ETHIBOND EXCEL™ Polyester Suture



TAPERCUT™ 



Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-4 17mm 3/8 Circle	4-0	Green	1 × 30"/75cm	36		X925H
	3-0	Green	1 × 30"/75cm	36		X926H
	2-0	Green	1 × 30"/75cm	36		X927H
	2-0	Green	4 × 30"/75cm	12		MX927
	2-0	White	4 × 30"/75cm	12		MX944
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget*	PX22N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget†	SXP22N
 V-5 17mm 1/2 Circle	5-0	Green	1 X 30"/75cm	36		X951H
 V-7 26mm 1/2 Circle	3-0	Green	1 × 30"/75cm	36		X963H
	2-0	Green	1 × 30"/75cm	36		X964H
 V-34 36mm 1/2 Circle	0	Green	1 × 30"/75cm	36		X444H
 V-40 48mm 1/2 Circle	5	Green	4 × 30"/75cm	12		MB46G

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
 †7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

ETHIBOND EXCEL™ Polyester Suture




TAPERCUT™ 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-5 17mm 1/2 Circle</p>	5-0	Green	1 × 30"/75cm	36		X934H
	4-0	Green	1 × 30"/75cm	36		X935H
	3-0	Green	1 × 30"/75cm	36		X936H
	3-0	White	1 × 30"/75cm	36		X916H
	3-0	Grn/Wht	8 × 30"/75cm	12	Pledget [†]	SXX53N
	3-0	Grn/Wht	10 × 30"/75cm	12		10X53N
	3-0	Green	1 × 36"/90cm	36		X806H
	2-0	Green	1 × 30"/75cm	36	Pledget*	PX937H
	2-0	Green	1 × 30"/75cm	36		X937H
	2-0	White	1 × 30"/75cm	36		X917H
	2-0	Green	4 × 30"/75cm	12	Pledget*	MPX37
	2-0	White	4 × 30"/75cm	12	Pledget*	MPX17
	2-0	Grn/Wht	8 × 30"/75cm	6	Pledget*	PXX50N
	2-0	Grn/Wht	8 × 30"/75cm	6	Pledget [†]	SXX50N
	2-0	Grn/Wht	10 × 30"/75cm	6	Pledget*	PXX52N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget [†]	SXP52N
	2-0	Grn/Wht	10 × 30"/75cm	6		10X52N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget [†]	SXP54
	2-0	Grn/Wht	10 × 30"/75cm	6	Pledget [†]	SXX54N
	 <p>V-7 26mm 1/2 Circle</p>	2-0	Green	1 × 36"/90cm	36	Pledget [†]
2-0		Green	1 × 36"/90cm	36	Pledget*	PX54H
2-0		Green	1 × 36"/90cm	36		X807H
2-0		White	1 × 36"/90cm	36	Pledget [†]	SX17H
2-0		White	1 × 36"/90cm	36	Pledget*	PX17H
4-0		Green	1 × 30"/75cm	36		X975H
3-0		Green	1 × 30"/75cm	36		X976H
2-0		Green	1 × 30"/75cm	36	Pledget*	PX977H
2-0		Green	1 × 30"/75cm	36	Pledget [†]	SXP07H
2-0		Green	1 × 30"/75cm	36		X977H
2-0	White	1 × 30"/75cm	36		X987H	
2-0	Green	4 × 30"/75cm	12		MX977	

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
[†]3mm x 3mm x 1.5mm Firm TFE Polymer Pledget
[‡]7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

ETHIBOND EXCEL™ Polyester Suture




TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-7 26mm 1/2 Circle	2-0	Grn/Wht	10 × 30"/75cm	6	Pledget*	PXX77N
	2-0	Grn/Wht	10 × 30"/75cm	12	Pledget†	SXP77N
	2-0	Grn/Wht	10 × 30"/75cm	6		10X72N
	2-0	Green	1 × 36"/90cm	36	Pledget*	PX62H
	2-0	Green	1 × 36"/90cm	36		X966H
	2-0	White	1 × 36"/90cm	36	Pledget*	PX64H
	2-0	Grn/Wht	8 × 36"/90cm	6	Pledget*	PXX79N
	2-0	Grn/Wht	10 x 36"/90cm	12		10X66N
	0	Green	1 × 36"/90cm	36		X905H
	0	White	1 × 36"/90cm	36	Pledget*	PX65H
 V-37 40mm 1/2 Circle	2	Green	4 × 30"/75cm	12		MX69G
	5	Green	4 × 30"/75cm	12		MB66G

Trocar

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
STP-10 250mm Straight	2-0	Green	1 × 27"/70cm	12		X997G

Reverse Cutting







Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	4-0	Green	1 × 18"/45cm	12		X683G

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
 †7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

ETHIBOND EXCEL™ Polyester Suture



Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-1 24mm 3/8 Circle	3-0	Green	1 x 18"/45cm	36		X663H
 FS 26mm 3/8 Circle	2-0	Green	1 x 18"/45cm	36		X664H
 LR 75mm 3/8 Circle	2	Green	1 x 20"/50cm	24		X406T
	5	Green	1 x 20"/50cm	24		B409T
 LR 75mm 3/8 Circle	2	Green	1 x 30"/75cm	24		X496T
	5	Green	1 x 30"/75cm	24		B499T
 OS-2 17mm 1/2 Circle	0	Green	3 x 18"/45cm	24	CR	CX74T
	1	Green	3 x 18"/45cm	24	CR	CX75T
 OS-4 22mm 1/2 Circle	0	Green	1 x 30"/75cm	36		X517H
	1	Green	1 x 30"/75cm	36		X518H
	2	Green	1 x 30"/75cm	36		X519H


ETHIBOND EXCEL™ Polyester Suture



Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 OS-6 36mm 1/2 Circle	1	Green	1 × 30"/75cm	36		X538H
 OS-8 40mm 1/2 Circle	1	Green	1 × 30"/75cm	36		X548H

Straight Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Green	1 × 30"/75cm	36		X622H


PRIME Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-3 13mm 3/8 Circle	5-0	Green	1 × 18"/45cm	12		X698G
 PS-3 16mm 3/8 Circle	5-0	Green	1 × 18"/45cm	12		X681G
 PS-2 19mm 3/8 Circle	4-0	White	1 × 18"/45cm	12		X692G
 PS-4 16mm 1/2 Circle	4-0	White	1 × 18"/45cm	12		X695G

ETHIBOND EXCEL™ Polyester Suture




Cutting Edge Conventional

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CCS 48mm 1/2 Circle</p>	5	Green	4 × 30"/75cm	12		MB47G

SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Green	12 × 18"/45cm	36		X182H
	4-0	Green	12 × 18"/45cm	36		X183H
	4-0	Green	12 × 30"/75cm	36		X303H
	3-0	Green	12 × 18"/45cm	36		X184H
	3-0	Green	12 × 30"/75cm	36		X304H
	2-0	Green	12 × 18"/45cm	36		X185H
	2-0	Green	12 × 30"/75cm	36		X305H
	0	Green	6 × 18"/45cm	36		X186H
	0	Green	6 × 30"/75cm	36		X306H



ETHILON™


Nylon Suture


















ETHILON Suture is nonabsorbable, sterile surgical monofilament suture composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6,6. ETHILON Suture is dyed black or green to enhance visibility in tissue.¹

ETHILON Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹




ETHILON™ Nylon Suture




Taperpoint 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 BV130-3 3.8mm 3/8 Circle	10-O	Black	1 × 5"/13cm	12		2820G	
	9-O	Black	1 × 5"/13cm	12		2819G	
	8-O	Black	1 × 5"/13cm	12		2822G	
 BV50-3 3.8mm 3/8 Circle	11-O	Black	1 × 5"/13cm	12		2881G	
 BV75-3 3.8mm 3/8 Circle	10-O	Black	1 × 5"/13cm	12		2870G	
 BV130-4 5mm 3/8 Circle	9-O	Black	1 × 5"/13cm	12	 MultiPass™ NEEDLE	2813G	
	8-O	Black	1 × 5"/13cm	12	 MultiPass™ NEEDLE	2815G	
 BV100-4 5.1mm 3/8 Circle	10-O	Black	1 × 5"/13cm	12		2830G	
	9-O	Black	1 × 5"/13cm	12		2829G	
 BV75-4 5.1mm 3/8 Circle	10-O	Black	1 × 5"/13cm	12		2850G	
 BV130-5 6.5mm 3/8 Circle	10-O	Black	1 × 5"/13cm	12	 MultiPass™ NEEDLE	2810G	
	9-O	Black	1 × 5"/13cm	12	 MultiPass™ NEEDLE	2809G	
	8-O	Black	1 × 5"/13cm	12	 MultiPass™ NEEDLE	2808G	
 BV130-5 6.5mm 3/8 Circle	9-O	Black	1 × 10"/25cm	12	 MultiPass™ NEEDLE	2800G	
 BVH100-3 5mm 1/2 Circle	10-O	Black	1 × 5"/13cm	12		2860G	
 CT 40mm 1/2 Circle	0	Black	1 × 96" Looped 48"	24	Looped	L886T	




ETHILON™ Nylon Suture




Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTX 48mm 1/2 Circle	1	Black	1 × 30"/75cm	36		830H
 TP-1 65mm 1/2 Circle	0	Black	1 × 96" Looped 48"	12	Looped	L880G
	1	Black	1 × 60"/150cm	12		824G
	2	Black	1 × 60"/150cm	12		825G
 ST100-3 3mm Straight	10-0	Black	1 × 5"/13cm	12		2849G


TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V100-3 3.5mm 1/4 Circle	9-0	Black	1 × 5"/13cm	12		2894G
 V130-3 3.8mm 1/4 Circle	9-0	Black	1 × 5"/13cm	12		2898G
	8-0	Black	1 × 5"/13cm	12		2897G
 V75-3 3.5mm 1/4 Circle	10-0	Black	1 × 5"/13cm	12		2892G








Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-2 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		667G


ETHILON™ Nylon Suture




Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-3 13mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		668G
 FS-2 19mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		661G
	5-0	Black	1 × 18"/45cm	36		661H
	4-0	Black	1 × 18"/45cm	12		662G
	4-0	Black	1 × 18"/45cm	36		662H
 FS-1 24mm 3/8 Circle	4-0	Black	1 × 18"/45cm	36		1629H
	3-0	Black	1 × 18"/45cm	12		663G
	3-0	Black	1 × 30"/75cm	36		669H
 FS 26mm 3/8 Circle	2-0	Black	1 × 18"/45cm	12		664G
	2-0	Black	1 × 18"/45cm	36		664H
 FSL 30mm 3/8 Circle	4-0	Black	1 × 30"/75cm	36		1670BH
 PSL 30mm 3/8 Circle	4-0	Black	1 × 30"/75cm	36		562H
	3-0	Black	1 × 30"/75cm	36		1691H
 FSLX 36mm 3/8 Circle	2-0	Black	1 × 30"/75cm	36		1674BH


ETHILON™ Nylon Suture




Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 LR 75mm 3/8 Circle	2	Black	1 x 20"/50cm	24		460T
 LR 75mm 3/8 Circle	1	Black	1 x 30"/75cm	24		489T
 LR 75mm 3/8 Circle	2	Black	1 x 30"/75cm	24		490T
 X-1 22mm 1/2 Circle	3-0	Black	1 x 18"/45cm	12		642G
 PSLX 36mm 3/8 Circle	3-0	Black	1 x 30"/75cm	36		1683H
	2-0	Black	1 x 30"/75cm	36		1697H

Straight Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Black	1 x 30"/75cm	36		626H
	3-0	Black	1 x 30"/75cm	36		627H
	2-0	Black	1 x 30"/75cm	36		628H


ETHILON™ Nylon Suture







PRIME Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-6 8mm 1/4 Circle	7-0	Black	1 × 18"/45cm	12		1647G
 P-1 11mm 3/8 Circle	7-0	Black	1 × 18"/45cm	12		1696G
	6-0	Undyed	1 × 18"/45cm	12		689G
	6-0	Black	1 × 18"/45cm	12		697G
	6-0	Black	1 × 18"/45cm	36		697H
	6-0	Green	1 × 18"/45cm	12		G697G
	5-0	Green	1 × 18"/45cm	12		G695G
 P-3 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		1698G
	6-0	Green	1 × 18"/45cm	12		G696G
	5-0	Undyed	1 × 18"/45cm	12		690G
	5-0	Black	1 × 18"/45cm	12		698G
	5-0	Black	1 × 18"/45cm	36		698H
	5-0	Green	1 × 18"/45cm	12		G698G
	4-0	Undyed	1 × 18"/45cm	12		691G
	4-0	Black	1 × 18"/45cm	12		699G
	4-0	Black	1 × 18"/45cm	36		699H
	 PS-3 16mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12	
5-0		Black	1 × 18"/45cm	12		1668G
5-0		Black	1 × 18"/45cm	36		1668H
 PS-2 19mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		1666G
	5-0	Black	1 × 18"/45cm	36		1666H
	5-0	Green	1 × 18"/45cm	12		G666G
	4-0	Undyed	1 × 18"/45cm	12		1611G
	4-0	Black	1 × 18"/45cm	12		1667G
	4-0	Black	1 × 18"/45cm	36		1667ZH
	4-0	Green	1 × 18"/45cm	12		G667G
	3-0	Black	1 × 18"/45cm	36		1669H
	2-0	Black	1 × 18"/45cm	36		593H









ETHILON™ Nylon Suture




PRIME Reverse Cutting ▽

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-1 24mm 3/8 Circle	3-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1663G
	3-0	Black	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	1663H
 PS 26mm 3/8 Circle	2-0	Black	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	585H
 PS-4 16mm 1/2 Circle	6-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1660G
	5-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1661G
	4-0	Undyed	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1603G
	4-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1662G

PRIME Conventional Cutting ▲

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1856G
	5-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1855G
	4-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1854G
 PC-3 16mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1866G
	5-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1865G
	4-0	Black	1 × 18"/45cm	12	 MultiPass [®] NEEDLES	1864G


ETHILON™ Nylon Suture






PRIME Conventional Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-5 19mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	1895G
	5-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	
	4-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	1894G
	4-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	
	3-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	1893G
	3-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	
 PC-12 16mm 1/2 Circle	5-0	Black	1 × 18"/45cm	12		1845G


SABRELOC™ Centerpoint Spatula

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RD-1 8mm 1/4 Circle	5-0	Black	1 × 18"/45cm	12		749G








Spatula Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG100-8 6.5mm 1/4 Circle	8-0	Black	1 × 12"/30cm	12		1714G
 S-14 8mm 1/4 Circle	5-0	Black	1 × 18"/45cm	12		7731G
 VAS100-4 5mm 3/8 Circle	9-0	Black	1 × 6"/15cm	12		2890G



ETHILON™ Nylon Suture




Spatula Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG140-6 6.5mm 3/8 Circle	10-0	Black	1 × 12"/30cm	12		7770G
 TG140-8 6.5mm 3/8 Circle	10-0	Black	1 × 12"/30cm	12		7718G
	9-0	Black	1 × 12"/30cm	12		7717G
 TG160-4 5.3mm 1/2 Circle	10-0	Black	1 × 12"/30cm	12		7707G
 TG160-4-3M 4.3mm 1/2 Circle	10-0	Black	1 × 8"/20cm	12		7711G
 TG160-6 5.5mm 1/2 Circle	10-0	Black	1 × 12"/30cm	12		7756G
	10-0	Black	1 × 6"/15cm	12		7757G
	9-0	Black	1 × 12"/30cm	12		7760G
 TG175-8 7mm 1/2 Circle	8-0	Black	1 × 12"/30cm	12		1716G
 S-22 8mm 1/2 Circle	5-0	Black	1 × 18"/45cm	12		1740G









ULTIMA™ Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CS44-6 6.5mm 1/8 Circle	10-0	Black	1 × 4"/10cm	12		9065G
 CS140-6 6.5mm 3/8 Circle	10-0	Black	1 × 4"/10cm	12		9060G

ETHILON™ Nylon Suture



ULTIMA™ Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CSB-6 5mm Compound Curve	10-0	Black	1 × 12"/30cm	12		9007G
 CS140-6 6.5mm 3/8 Circle	9-0	Black	1 × 12"/30cm	12		9013G
 CS140-8 6.5mm 3/8 Circle	10-0	Black	1 × 12"/30cm	12		9003G
 CSM-6 4.3mm 1/2 Circle	10-0	Black	1 X 8"/20cm	12		9033G
 CS160-6 5.5mm 1/2 Circle	10-0	Black	1 × 12"/30cm	12		9006G
 CS160-6 5.5mm 1/2 Circle	10-0	Black	1 × 4"/10cm	12		9005G
 CS160-8 5.5mm 1/2 Circle	11-0	Black	1 X 8"/20cm	12		9061G
 CS175-6 7mm 1/2 Circle	10-0	Black	1 X 6"/15cm	12		9016G
 CS175-8 7mm 1/2 Circle	10-0	Black	1 X 6"/15cm	12		9001G
 CS-C-6 5mm Compound Curve	10-0	Black	1 × 12"/30cm	12		9000G
	9-0	Black	1 × 12"/30cm	12		9011G
	10-0	Black	1 × 12"/30cm	12		9031G
	10-0	Black	1 × 6"/15cm	12		9030G
	10-0	Black	1 × 12"/30cm	12		9002G
	10-0	Black	1 × 12"/30cm	12		9032G
	9-0	Black	1 × 12"/30cm	12		9036G
	8-0	Black	1 × 12"/30cm	12		9038G
	10-0	Black	1 × 6"/15cm	12		9010G



MERSILENE™ Polyester Fiber Suture

MERSILENE Suture is a nonabsorbable, braided, sterile surgical suture composed of Poly (ethylene terephthalate).¹



MERSILENE Sutures are indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹



MERSILENE™ Polyester Fiber Suture



Taperpoint

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 SH 26mm 1/2 Circle	3-0	Green/Braided	1 × 30"/75cm	36		R832H
	2-0	Green/Braided	1 × 30"/75cm	36		R833H
 CT-1 36mm 1/2 Circle	0	Green/Braided	1 × 30"/75cm	36		R424H
	0	Green/Braided	8 × 18"/45cm	24		MR21T


Reverse Cutting

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	4-0	White/Braided	1 × 18"/45cm	36		R633H
 FS-1 24mm 3/8 Circle	3-0	White/Braided	1 × 18"/45cm	36		R647H
 FS 26mm 3/8 Circle	2-0	White/Braided	1 × 18"/45cm	36		R665H




MERSILENE™ Polyester Fiber Suture



Reverse Cutting ▾

Needle	Suture USP	Color/Construction	Length	EA/BX	Product Code
 <p>LR 75mm 3/8 Circle</p>	1	Green/Braided	1 x 30"/75cm	24	R495T

PRIME Reverse Cutting ▾

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 <p>P-3 13mm 3/8 Circle</p>	5-0	White/Braided	1 x 18"/45cm	12	 MultiPass [®] NEEDLES	R690G
	4-0	White/Braided	1 x 18"/45cm	12	 MultiPass [®] NEEDLES	R691G







SABRELOC™ Centerpoint Spatula

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 <p>RD-1 8mm 1/4 Circle</p>	5-0	Green/Braided	1 x 18"/45cm	12		740G
	5-0	White/Braided	1 x 18"/45cm	12		746G



MERSILENE™ Polyester Fiber Suture



Spatula Reverse Cutting




Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 S-14 8mm 1/4 Circle	6-0	White/Braided	1 × 18"/45cm	12		1722G
	5-0	White/Braided	1 × 18"/45cm	12		1760G
	5-0	Green/Braided	1 × 18"/45cm	12		1761G
	5-0	White/Braided	1 × 18"/45cm	12		1764G
 S-24 8mm 1/4 Circle	4-0	Green/Braided	1 × 18"/45cm	12		1762G
 S-4 8mm 1/4 Circle	4-0	White/Braided	1 × 18"/45cm	12		1778G
 TG160-4 5.3mm 1/2 Circle	11-0	Green/ Monofilament	1 × 12"/30cm	12		R750G
 TG160-4-3M 4.3mm 1/2 Circle	11-0	Green/Monofilament	1 × 8"/20cm	12		R749G
 S-2 8mm 1/2 Circle	4-0	White/Braided	1 × 18"/45cm	12		1779G


ULTIMA™ Reverse Cutting

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 CS140-6 6.5mm 3/8 Circle	10-0	Green/Monofilament	1 × 12"/30cm	12		R770G
 CS160-6 5.5mm 1/2 Circle	10-0	Green/Monofilament	1 × 12"/30cm	12		R756G

MERSILENE™ Polyester Fiber Strip



Taperpoint 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	5mm Tape	White	1 × 12"/30cm	6		RS22
 <p>MO-4 36mm 1/2 Circle</p>	5mm Tape	White	1 × 12"/30cm	6		RS23




NUROLON™ Nylon Suture























NUROLON Suture is a nonabsorbable sterile surgical braided suture composed of the long-chain aliphatic polymer Nylon 6,6. It is dyed black to enhance visibility in tissue.¹

NUROLON Suture is indicated for use in general soft tissue approximating and/or ligation, including use in cardiovascular, ophthalmic procedures, and neurological procedures.¹



















NUROLON™ Nylon Suture




Taperpoint 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 TF 13mm 1/2 Circle	4-0	Black	8 × 18"/45cm	12		C584D*	
	4-0	Black	8 × 18"/45cm	24		N104T	
 RB-1 17mm 1/2 Circle	4-0	Black	8 × 18"/45cm	12		C554D*	
	4-0	Black	8 × 18"/45cm	24		N124T	
	3-0	Black	8 × 18"/45cm	12		C553D*	
 SH-1 22mm 1/2 Circle	3-0	Black	8 × 18"/45cm	12		C503D*	
	2-0	Black	8 × 18"/45cm	12		C502D*	
 MO-7 22mm 1/2 Circle	2-0	Black	8 × 18"/45cm	12		C542D*	
	0	Black	8 × 18"/45cm	12		C541D*	
 SH 26mm 1/2 Circle	4-0	Black	8 × 18"/45cm	12		C514D*	
	3-0	Black	8 × 18"/45cm	12		C513D*	
	2-0	Black	8 × 18"/45cm	12		C512D*	
 CT-2 26mm 1/2 Circle	2-0	Black	8 × 18"/45cm	12		C526D*	
	0	Black	8 × 18"/45cm	12		C527D*	
 MO-6 26mm 1/2 Circle	2-0	Black	8 × 18"/45cm	12		C546D*	
	0	Black	8 × 18"/45cm	12		C545D*	

NUROLON™ Nylon Suture






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Black	8 × 18"/45cm	12		C522D*
	0	Black	1 × 30"/75cm	36		5424H
	0	Black	8 × 18"/45cm	12		C521D*
	1	Black	1 × 30"/75cm	36		5425H
	1	Black	8 × 18"/45cm	12		C520D*
 <p>CTX 48mm 1/2 Circle</p>	0	Black	8 × 18"/45cm	12		C551D*
	1	Black	8 × 18"/45cm	12		C550D*
Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>FS 26mm 3/8 Circle</p>	4-0	Black	8 × 18"/45cm	12		C570D*
 <p>FSL 30mm 3/8 Circle</p>	3-0	Black	8 × 18"/45cm	12		C566D*
 <p>X-1 22mm 1/2 Circle</p>	3-0	Black	8 × 18"/45cm	12		C563D*




NUROLON™ Nylon Suture




Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-2 26mm 1/2 Circle	3-0	Black	8 × 18"/45cm	12	CR	C576D*
	2-0	Black	8 × 18"/45cm	12	CR	C575D*
 CP 40mm 1/2 Circle	2-0	Black	8 × 18"/45cm	12	CR	C572D*
 OS-8 40mm 1/2 Circle	1	Black	3 × 18"/45cm	24	CR	C595T

PRIME Conventional Cutting ▲

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	5666G
	5-0	Black	1 × 18"/45cm	12	 MultiPass™ NEEDLES	5665G



PERMA-HAND™ Silk Suture


PERMA-HAND Suture is a nonabsorbable, sterile, surgical suture composed of an organic protein called fibroin.¹


















The sutures are processed to remove the natural waxes and gums, and then are dyed black and coated with a special wax mixture. PERMA-HAND Suture is also available in its natural color, also available is PERMA-HAND Virgin Silk Suture in which the sericin gum is not removed and serves to hold the filaments together.¹

PERMA-HAND Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹












PERMA-HAND™ Silk Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/ BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	6-0	Black	1 × 18"/45cm	36		K802H
 BV-1 9.3mm 3/8 Circle	7-0	Black	1 × 18"/45cm	36		K809H
 C-1 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	36		K801H
 C-1 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		707G
 C-1 13mm 3/8 Circle	6-0	Black	1 × 30"/75cm	36		K889H
 C-1 13mm 3/8 Circle	5-0	Black	1 × 30"/75cm	36		K890H
 C-1 13mm 3/8 Circle	4-0	Black	1 × 30"/75cm	36		K891H
 C-1 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		706G
 C-1 13mm 3/8 Circle	4-0	Black	1 × 18"/45cm	36		709H
 BB 17mm 3/8 Circle	5-0	Black	1 × 30"/75cm	36		K880H
 BB 17mm 3/8 Circle	4-0	Black	1 × 30"/75cm	36		K881H
 BB 17mm 3/8 Circle	3-0	Black	1 × 30"/75cm	36		K882H
 BB 17mm 3/8 Circle	2-0	Black	1 × 30"/75cm	36		K883H
 BB 17mm 3/8 Circle	5-0	Black	1 × 24"/60cm	36		K580H
 BB 17mm 3/8 Circle	4-0	Black	1 × 24"/60cm	36		K581H
 RB-2 13mm 1/2 Circle	4-0	Black	1 × 18"/45cm	36		N267H

PERMA-HAND™ Silk Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	5-0	Black	1 × 18"/45cm	36		N266H
	4-0	Black	8 × 18"/45cm	12	CR	C084D*
	4-0	Black	8 × 18"/45cm	24		M104T
 RB-1 17mm 1/2 Circle	5-0	Black	1 × 30"/75cm	36		K870H
	4-0	Black	8 × 18"/45cm	12	CR	C054D*
	4-0	Black	1 × 30"/75cm	36		K871H
	3-0	Black	8 × 18"/45cm	12	CR	C053D*
	3-0	Black	1 × 30"/75cm	36		K872H
	2-0	Black	1 × 30"/75cm	36		K873H
 RB-1 17mm 1/2 Circle	5-0	Black	1 × 24"/60cm	36		K570H
	4-0	Black	1 × 24"/60cm	36		K571H
	3-0	Black	1 × 24"/60cm	36		K572H
 CT-3 22mm 1/2 Circle	2-0	Black	1 × 20"/50cm	36		413H
 MO-7 22mm 1/2 Circle	0	Black	8 × 18"/45cm	12	CR	C041D*
 SH-1 22mm 1/2 Circle	4-0	Black	8 × 18"/45cm	12	CR	C004D*
	3-0	Black	8 × 18"/45cm	12	CR	C003D*
	3-0	Black	8 × 18"/45cm	24		M3T
 CT-2 26mm 1/2 Circle	2-0	Black	8 × 18"/45cm	12	CR	CO26D*
	0	Black	8 × 18"/45cm	12	CR	CO27D*

*Safety organizer tray

CR CONTROL RELEASE™ Needle Suture


PERMA-HAND™ Silk Suture









Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>MO-6 26mm 1/2 Circle</p>	0	Black	8 × 18"/45cm	12		C045D*
	4-0	Black	8 × 18"/45cm	12		C014D*
	4-0	Black	1 × 30"/75cm	36		K831H
	4-0	Black	8 × 30"/75cm	12		C018D*
	3-0	Black	8 × 18"/45cm	12		C013D*
	3-0	Black	8 × 18"/45cm	24		M13T
	3-0	Black	1 × 30"/75cm	36		K832H
	3-0	Black	8 × 30"/75cm	12		C017D*
	2-0	Black	8 × 18"/45cm	12		C012D*
	2-0	Black	1 × 30"/75cm	36		K833H
	2-0	Black	8 × 30"/75cm	12		C016D*
	0	Black	1 × 30"/75cm	36		K834H
	1	Black	1 × 30"/75cm	36		K835H
	 <p>SH 26mm 1/2 Circle</p>	2-0	Black	1 × 24"/60cm	36	
0		Black	1 × 24"/60cm	36		K534H
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Black	8 × 18"/45cm	12		CO22D*
	2-0	Black	1 × 30"/75cm	36		423H
	0	Black	8 × 18"/45cm	12		CO21D*
	0	Black	1 × 30"/75cm	36		424H

*Safety organizer tray


PERMA-HAND™ Silk Suture










Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 MH 36mm 1/2 Circle	2-0	Black	8 × 24"/60cm	24		M72T
	2-0	Black	1 × 30"/75cm	36		K843H
	0	Black	1 × 30"/75cm	36		K844H
	1	Black	1 × 30"/75cm	36		K845H
 CT-1 36mm 1/2 Circle	2	Black	1 × 30"/75cm	36		K846H
 TF-5 8mm 5/8 Circle	5-0	Black	1 × 18"/45cm	36		N271H
TF-5 8mm 5/8 Circle	4-0	Black	1 × 18"/45cm	36		N272H
 ST-1 45mm Straight	3-0	Black	1 × 30"/75cm	36		K852H
 C-1 13mm 3/8 Circle	4-0	Black	1 x12"/30cm	12		7734G


TAPERCUT™ 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-4 17mm 3/8 Circle	4-0	Black	1 × 18"/45cm	12		K945G
 V-5 17mm 1/2 Circle	4-0	Black	1 × 30"/75cm	36		K952H
 V-7 26mm 1/2 Circle	2-0	Black	1 × 30"/75cm	36		K964H
 CC-1 13mm 3/8 Circle	6-0	Black	1 × 12"/30cm	12		7729G












PERMA-HAND™ Silk Suture




Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-2 13mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		714G
 C-3 13mm 3/8 Circle	4-0	Black	1 × 18"/45cm	12		736G
 C-3 13mm 3/8 Circle	4-0	Black	1 × 12"/30cm	12		735G
 FS-2 19mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		682G
	4-0	Black	1 × 18"/45cm	12		683G
 FS-1 24mm 3/8 Circle	4-0	Black	1 × 18"/45cm	12		629G
	3-0	Black	1 × 18"/45cm	12		684G
 FS 26mm 3/8 Circle	2-0	Black	1 × 18"/45cm	12		685G
 FSL 30mm 3/8 Circle	3-0	Black	1 × 30"/75cm	36		675H
	2-0	Black	1 × 18"/45cm	12		677G
	2-0	Black	1 × 30"/75cm	36		679H
	0	Black	1 × 18"/45cm	12		678G
	0	Black	1 × 30"/75cm	36		680H

PERMA-HAND™ Silk Suture







Reverse Cutting ▾							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 <p>PSLX 36mm 3/8 Circle</p>	2-0	Black	1 x 30"/75cm	36		591H	
 <p>LR 75mm 3/8 Circle</p>	2	Black	1 x 30"/75cm	24		482T	
 <p>CP-1 36mm 1/2 Circle</p>	5	Black	1 x 30"/75cm	24		485T	
 <p>J-1 17mm 1/2 Circle</p>	4-0	Black	1 x 18"/45cm	12		734G	
 <p>X-1 22mm 1/2 Circle</p>	3-0	Black	1 x 18"/45cm	12		632G	
 <p>CP-1 36mm 1/2 Circle</p>	2-0	Black	1 x 18"/45cm	12		737G	
 <p>CP-1 36mm 1/2 Circle</p>	2-0	Black	1 x 30"/75cm	36		443H	
 <p>G-7 8mm 1/2 Circle</p>	7-0	Black	1 x 18"/45cm	12		766G	
 <p>G-7 8mm 1/2 Circle</p>	6-0	Black	1 x 18"/45cm	12		765G	
 <p>G-6 8mm 1/4 Circle</p>	7-0	Black	1 x 18"/45cm	12		768G	
 <p>G-6 8mm 1/4 Circle</p>	6-0	Black	1 x 18"/45cm	12		769G	


PERMA-HAND™ Silk Suture










Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-1 11mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		786G
 G-1 11mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		780G
 G-3 13mm 3/8 Circle	4-0	Black	1 × 18"/45cm	12		789G
 G-3 13mm 3/8 Circle	4-0	Black	1 × 18"/45cm	12		783G


Straight Cutting ▾







Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Black	1 × 30"/75cm	36		621H
	3-0	Black	1 × 30"/75cm	36		622H
	2-0	Black	1 × 30"/75cm	36		623H
	0	Black	1 × 30"/75cm	36		624H






PRIME Reverse Cutting ▾




Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Black	1 × 18"/45cm	12		639G
	6-0	Black	1 × 18"/45cm	12		1639G
 P-3 13mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		640G
	4-0	Black	1 × 18"/45cm	12		641G
 PS-3 16mm 3/8 Circle	6-0	Black	1 × 18"/45cm	36		1675H
	5-0	Black	1 × 18"/45cm	36		1678H
 PS-2 19mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		1676G
	4-0	Black	1 × 18"/45cm	12		1677G
	3-0	Black	1 × 18"/45cm	36		1679H
 PS-1 24mm 3/8 Circle	2-0	Black	1 × 18"/45cm	36		583H
	3-0	Black	1 × 18"/45cm	12		1684G

PERMA-HAND™ Silk Suture




PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PSL 30mm 3/8 Circle	2-0	Black	1 × 18"/45cm	36		673H
	2-0	Black	1 × 30"/75cm	36		486H
	0	Black	1 × 18"/45cm	36		580H
	0	Black	1 × 30"/75cm	36		590H
 PS 26mm 3/8 Circle	2-0	Black	1 × 18"/45cm	36		1588H
 PS-4 16mm 1/2 Circle	4-0	Black	1 × 18"/45cm	12		1687G

PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-3 16mm 3/8 Circle	5-0	Black	1 × 18"/45cm	12		1925G
 PC-5 19mm 3/8 Circle	4-0	Black	1 × 18"/45cm	12		1984G

Spatula Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 S-14 8mm 1/4 Circle	6-0	Black	1 × 18"/45cm	12		1780G
 TG140-8 6.5mm 3/8 Circle	7-0	Black	1 × 18"/45cm	12		7733G
	6-0	Black	1 × 18"/45cm	12		1732G

PERMA-HAND™ Silk Suture




LIGAPAK™ Dispensing Reel

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Black	1 x 144"/365cm	12		LA53G
	3-0	Black	1 x 144"/365cm	12		LA54G
	2-0	Black	1 x 144"/365cm	12		LA55G

SUTUPAK™ Pre-Cut Sterile Sutures


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Black	17 x 18"/45cm	36		SA63H
	4-0	Black	13 x 24"/60cm	36		SA73H
	4-0	Black	10 x 30"/75cm	36		SA83H
	3-0	Black	17 x 18"/45cm	24		SA64T
	3-0	Black	13 x 24"/60cm	36		SA74H
	3-0	Black	10 x 30"/75cm	36		SA84H
	2-0	Black	17 x 18"/45cm	24		SA65T
	2-0	Black	13 x 24"/60cm	36		SA75H
	2-0	Black	10 x 30"/75cm	36		SA85H
	2-0	Black	2 x 60"/150cm	36		SA5H
	0	Black	17x18"/45cm	12		SA66G
	0	Black	13 x 24"/60cm	12		SA76G
	0	Black	10 x 30"/75cm	12		SA86G
	0	Black	2 x 60"/150cm	36		SA6H
	1	Black	13 x 24"/60cm	12		SA77G
	1	Black	10 x 30"/75cm	12		SA87G
	1	Black	2 x 60"/150cm	36		SA7H
	2	Black	2 x 60"/150cm	36		SA8H
	3	Black	2 x 60"/150cm	12		SA9G
	4	Black	2 x 60"/150cm	12		SA10G
	5	Black	2 x 60"/150cm	12		SA11G

PERMA-HAND™ Silk Suture



SUTUPAK Pre-Cut Sutures in Labyrinth Package

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Black	12 × 18"/45cm	36		A182H
	5-0	Black	12 × 30"/75cm	36		A302H
	4-0	Black	12 × 18"/45cm	36		A183H
	4-0	Black	12 × 30"/75cm	36		A303H
	3-0	Black	12 × 18"/45cm	36		A184H
	3-0	Black	12 × 30"/75cm	36		A304H
	2-0	Black	12 × 18"/45cm	36		A185H
	2-0	Black	12 × 30"/75cm	36		A305H
	1	Black	6 × 18"/45cm	36		A187H
	1	Black	6 × 30"/75cm	36		A307H
	0	Black	6 × 18"/45cm	36		A186H
	0	Black	6 × 30"/75cm	36		A306H



PROLENE™


Polypropylene Suture

Widely used in general, cardiovascular, plastic, and orthopedic surgery, PROLENE Sutures do not adhere to tissue. PROLENE Sutures are relatively biologically inert, offering proven strength, reliability, and versatility.^{1,2}



PROLENE suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.^{1,3}














PROLENE™ Polypropylene Suture




ETHIGUARD™ Blunt Point Needle





































Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTB-1 36mm 1/2 Circle	0	Blue	1 × 30"/75cm	36		PB21H
	1	Blue	1 × 30"/75cm	36		PB20H
 CTXB 48mm 1/2 Circle	1	Blue	1 × 30"/75cm	36		PB30H

Taperpoint


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV75-3 3.8mm 3/8 Circle	10-O	Blue	1 × 5"/13cm	12		2790G
 BV100-4 5.1mm 3/8 Circle	10-O	Blue	1 × 5"/13cm	12		2794G
	9-O	Blue	1 × 5"/13cm	12		2793G
 BV130-4 5.1mm 3/8 Circle	10-O	Blue	1 × 2"/5cm	12		2795G
 BV130-5 6.5mm 3/8 Circle	8-O	Blue	1 × 5"/13cm	12		2775G
 BV130-5 6.5mm 3/8 Circle	8-O	Blue	1 × 18"/45cm	36		8730H
	8-O	Blue	1 × 18"/45cm	36		EP8730H
	8-O	Blue	2 × 18"/45cm	12		M8738
	8-O	Blue	2 × 18"/45cm	12		EPM8738
	8-O	Blue	4 × 18"/45cm	12		M8730
	8-O	Blue	4 × 18"/45cm	12		EPM8730





































PROLENE™ Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-5 6.5mm 3/8 Circle	8-0	Blue	1 × 24"/60cm	36		8732H
	8-0	Blue	1 × 24"/60cm	36	 	EP8732H
	8-0	Blue	2 × 24"/60cm	12		M8739
	8-0	Blue	2 × 24"/60cm	12	 	EPM8739
	8-0	Blue	4 × 24"/60cm	12		M8732
	8-0	Blue	4 × 24"/60cm	12	 	EPM8732
 BV175-6 8mm 3/8 Circle	8-0	Blue	1 × 18"/45cm	36		8740H
	8-0	Blue	1 × 18"/45cm	36	 	EP8740H
	8-0	Blue	2 × 18"/45cm	12		M8742
	8-0	Blue	2 × 18"/45cm	12	 	EPM8742
	8-0	Blue	1 × 24"/60cm	36		8741H
	8-0	Blue	1 × 24"/60cm	36	 	EP8741H
	8-0	Blue	2 × 24"/60cm	12		M8743
	8-0	Blue	2 × 24"/60cm	12	 	EPM8743
	8-0	Blue	4 × 24"/60cm	12		M8741
	8-0	Blue	4 × 24"/60cm	12	 	EPM8741
	7-0	Blue	1 × 18"/45cm	36		8734H
	7-0	Blue	1 × 18"/45cm	36	 	EP8734H
	7-0	Blue	2 × 18"/45cm	12		M8736
	7-0	Blue	2 × 18"/45cm	12	 	EPM8736
7-0	Blue	4 × 18"/45cm	12		M8734	
7-0	Blue	4 × 18"/45cm	12	 	EPM8734	





















PROLENE™ Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV175-6 8mm 3/8 Circle	7-0	Blue	1 × 24"/60cm	36	 MultiPass [®] NEEDLES	8735H
	7-0	Blue	1 × 24"/60cm	36	 EVP  MultiPass [®] NEEDLES	EP8735H
	7-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8737
	7-0	Blue	2 × 24"/60cm	12	 EVP  MultiPass [®] NEEDLES	EPM8737
	7-0	Blue	4 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8735
	7-0	Blue	4 × 24"/60cm	12	 EVP  MultiPass [®] NEEDLES	EPM8735
	7-0	Blue	2 × 30"/75cm	12	 MultiPass [®] NEEDLES	M8733
	7-0	Blue	2 × 30"/75cm	12	 EVP  MultiPass [®] NEEDLES	EPM8733
	7-0	Blue	4 × 30"/75cm	12	 MultiPass [®] NEEDLES	M8731
	7-0	Blue	4 × 30"/75cm	12	 EVP  MultiPass [®] NEEDLES	EPM8731
 BV175-7 8mm 3/8 Circle	8-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8767
	8-0	Blue	2 × 24"/60cm	12	 EVP  MultiPass [®] NEEDLES	EPM8767
	7-0	Blue	2 × 18"/45cm	12	 MultiPass [®] NEEDLES	M8765
	7-0	Blue	2 × 18"/45cm	12	 EVP  MultiPass [®] NEEDLES	EPM8765
	7-0	Blue	4 × 18"/45cm	12	 MultiPass [®] NEEDLES	M8764
	7-0	Blue	1 × 24"/60cm	36	 MultiPass [®] NEEDLES	8766H
	7-0	Blue	1 × 24"/60cm	36	 EVP  MultiPass [®] NEEDLES	EP8766H
	7-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8766
	7-0	Blue	2 × 24"/60cm	12	 EVP  MultiPass [®] NEEDLES	EPM8766
	7-0	Blue	4 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8775
	7-0	Blue	4 × 24"/60cm	12	 EVP  MultiPass [®] NEEDLES	EPM8775
	7-0	Blue	2 × 30"/75cm	12	 MultiPass [®] NEEDLES	M8768
	7-0	Blue	2 × 30"/75cm	12	 EVP  MultiPass [®] NEEDLES	EPM8768
	7-0	Blue	4 × 30"/75cm	12	 MultiPass [®] NEEDLES	M8785
	7-0	Blue	4 × 30"/75cm	12	 EVP  MultiPass [®] NEEDLES	EPM8785
	8-0	Blue	1x24"/60cm	12		EP8767LG












PROLENE™ Polypropylene Suture




Taperpoint 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 BV175-8 9.3mm 3/8 Circle	8-0	Blue	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	8752H	
	8-0	Blue	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	EP8752H	
	8-0	Blue	2 × 18"/45cm	12	 MultiPass [®] NEEDLES	M8752	
	8-0	Blue	2 × 18"/45cm	12	 MultiPass [®] NEEDLES	EPM8752	
	8-0	Blue	1 × 24"/60cm	36	 MultiPass [®] NEEDLES	8753H	
	8-0	Blue	1 × 24"/60cm	36	 MultiPass [®] NEEDLES	EP8753H	
	8-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8753	
	8-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	EPM8753	
	8-0	Blue	1 × 30"/75cm	36	 MultiPass [®] NEEDLES	8754H	
	8-0	Blue	1 × 30"/75cm	36	 MultiPass [®] NEEDLES	EP8754H	
	8-0	Blue	2 × 30"/75cm	12	 MultiPass [®] NEEDLES	M8754	
	8-0	Blue	2 × 30"/75cm	12	 MultiPass [®] NEEDLES	EPM8754	
	7-0	Blue	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	8746H	
	7-0	Blue	1 × 18"/45cm	36	 MultiPass [®] NEEDLES	EP8746H	
	7-0	Blue	2 × 18"/45cm	12	 MultiPass [®] NEEDLES	M8746	
	7-0	Blue	2 × 18"/45cm	12	 MultiPass [®] NEEDLES	EPM8746	
	7-0	Blue	1 × 24"/60cm	36	 MultiPass [®] NEEDLES	8747H	
	7-0	Blue	1 × 24"/60cm	36	 MultiPass [®] NEEDLES	EP8747H	



PROLENE™ Polypropylene Suture




Taperpoint 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 BV175-8 9.3mm 3/8 Circle	7-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8747	
	7-0	Blue	2 × 24"/60cm	12	 MultiPass [®] NEEDLES	EPM8747	
	7-0	Blue	2 × 30"/75cm	12	 MultiPass [®] NEEDLES	M8755	
	7-0	Blue	2 × 30"/75cm	12	 MultiPass [®] NEEDLES	EPM8755	
	7-0	Blue	4 × 24"/60cm	12	 MultiPass [®] NEEDLES	M8745	
	7-0	Blue	4 × 24"/60cm	12	 MultiPass [®] NEEDLES	EPM8745	
	7-0	Blue	4 × 30"/75cm	12	 MultiPass [®] NEEDLES	EPM8751	
	7-0	Blue	1 × 30"/75cm	36	 MultiPass [®] NEEDLES	8755H	
	7-0	Blue	1 × 30"/75cm	36	 MultiPass [®] NEEDLES	EP8755H	




PROLENE™ Polypropylene Suture


Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	7-0	Blue	1 x 18"/45cm	36		8701H
	7-0	Blue	1 x 18"/45cm	36	VB	8301H
	7-0	Blue	4 x 18"/45cm	12		M8701
	7-0	Blue	4 x 18"/45cm	12	VB	M8301
	7-0	Blue	1 x 24"/60cm	36		8702H
	7-0	Blue	1 x 24"/60cm	36	VB	8304H
	7-0	Blue	4 x 24"/60cm	12		M8702
	7-0	Blue	4 x 24"/60cm	12	VB	M8304
	7-0	Blue	4 x 24"/60cm	12	EVP	EPM8702
	7-0	Blue	1 x 30"/75cm	36		8703H
	7-0	Blue	2 x 30"/75cm	12		M8203
	7-0	Blue	4 x 30"/75cm	12		M8703
	7-0	Blue	4 x 30"/75cm	12	VB	M8303
	7-0	Blue	4 x 30"/75cm	12	EVP	EPM8703
	7-0	Blue	1 x 18"/45cm	36		EP8701H
	7-0	Blue	4 x 24"	36		EP8702H
	7-0	Blue	2X30"/75cm	12		EPM8203
	7-0	Blue	4X18"/45cm	12		EPM8701
	7-0	Blue	1X30"/75cm	36		EP8703H
	6-0	Blue	1 x 18"/45cm	36		8806H
	6-0	Blue	1 x 18"/45cm	36	EVP	EP8806H
	6-0	Blue	4 x 18"/45cm	12		M8806
	6-0	Blue	4 x 18"/45cm	12	EVP	EPM8806
	6-0	Blue	1 x 24"/60cm	36		8805H
	6-0	Blue	1 x 24"/60cm	36	EVP	EP8805H
	6-0	Blue	1 x 24"/60cm	36	VB	8305H
	6-0	Blue	2 x 24"/60cm	12		M8205
	6-0	Blue	2 x 24"/60cm	12	EVP	EPM8205
	6-0	Blue	4 x 24"/60cm	12		M8805
	6-0	Blue	4 x 24"/60cm	12	VB	M8305
	6-0	Blue	4 x 24"/60cm	12	EVP	EPM8805
	6-0	Blue	1 x 30"/75cm	36		8709H
	6-0	Blue	1 x 30"/75cm	36	EVP	EP8709H
	6-0	Blue	4 x 30"/75cm	12		M8709
	6-0	Blue	4 x 30"/75cm	12	EVP	EPM8709
5-0	Blue	1 x 24"/60cm	36	HS	9702H	
5-0	Blue	1 x 24"/60cm	36	HS EVP	EPH9702H	


PROLENE™ Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV 11mm 3/8 Circle	7-0	Blue	4 × 24"/60cm	12		M8600
	6-0	Blue	1 × 24"/60cm	36		8610H
	6-0	Blue	1 × 24"/60cm	36	EVP	EP8610H
	6-0	Blue	4 × 24"/60cm	12		M8610
	6-0	Blue	4 × 24"/60cm	12	EVP	EPM8610
	6-0	Blue	1 × 30"/75cm	36		8776H
	6-0	Blue	1 × 30"/75cm	36	EVP	EP8776H
	6-0	Blue	4 × 30"/75cm	12		M8776
	6-0	Blue	4 × 30"/75cm	12	EVP	EPM8776
	6-0	Blue	4X24"/60cm	12		EPM8610
 C-1 13mm 3/8 Circle	6-0	Blue	1 × 30"/75cm	36		8889H
	6-0	Blue	1 × 30"/75cm	36	EVP	EP8889H
	5-0	Blue	1 × 30"/75cm	36		8890H
	5-0	Blue	1 × 30"/75cm	36	HS EVP	EPH8890H
 C-1 13mm 3/8 Circle	7-0	Blue	1 × 18"/45cm	36		8800H
	7-0	Blue	1 × 18"/45cm	36	EVP	EP8800H
	7-0	Blue	4 × 18"/45cm	12		M8800
	7-0	Blue	4 × 18"/45cm	12	EVP	EPM8800
	7-0	Blue	4 × 24"/60cm	12		M8727
	7-0	Blue	4 × 24"/60cm	12	EVP	EPM8727
	7-0	Blue	1 × 30"/75cm	36		8708H
	7-0	Blue	1 × 30"/75cm	36	EVP	EP8708H

PROLENE™ Polypropylene Suture

Taperpoint 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>C-1 13mm 3/8 Circle</p>	6-0	Blue	1 × 18"/45cm	36		8718H
	6-0	Blue	1 × 18"/45cm	36	VB	8318H
	6-0	Blue	1 × 18"/45cm	36	EVP	EP8718H
	6-0	Blue	4 × 18"/45cm	12		M8718
	6-0	Blue	4 × 18"/45cm	12	VB	M8318
	6-0	Blue	4 × 18"/45cm	12	EVP	EPM8718
	6-0	Blue	1 × 24"/60cm	36		8726H
	6-0	Blue	1 × 24"/60cm	36	EVP	EP8726H
	6-0	Blue	4 × 24"/60cm	12		M8726
	6-0	Blue	4 × 24"/60cm	12	EVP	EPM8726
	6-0	Blue	1 × 30"/75cm	36		8706H
	6-0	Blue	1 × 30"/75cm	36	VB	8307H
	6-0	Blue	1 × 30"/75cm	36	EVP	EP8706H
	6-0	Blue	2 × 30"/75cm	12		M8206
	6-0	Blue	2 × 30"/75cm	12	EVP	EPM8206
	6-0	Blue	4 × 30"/75cm	12		M8706
	6-0	Blue	4 × 30"/75cm	12	VB	M8307
	6-0	Blue	4 × 30"/75cm	12	EVP	EPM8706
	5-0	Blue	1 × 18"/45cm	36		8717H
	5-0	Blue	4 × 18"/45cm	12		M8717
	5-0	Blue	1 × 24"/60cm	36		8725H
	5-0	Blue	1 × 24"/60cm	36	VB	8325H
	5-0	Blue	1 × 24"/60cm	36	HS EVP	EPH8725H
	5-0	Blue	4 × 24"/60cm	12		M8725
	5-0	Blue	4 × 24"/60cm	12	VB	M8325
	5-0	Blue	1 × 30"/75cm	36	HS EVP	EPH7477H
	5-0	Blue	1 × 36"/90cm	36		8720ZH
	5-0	Blue	1 × 36"/90cm	36	VB	8321H
	5-0	Blue	1 × 36"/90cm	36	HS EVP	EPH8720H
	5-0	Blue	2 × 36"/90cm	12		M8220
5-0	Blue	4 × 36"/90cm	12		M8720	
5-0	Blue	4 × 36"/90cm	12	VB	M8321	







PROLENE™ Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Blue	1 × 30"/75cm	36	HS	9706H
C-1 13mm 3/8 Circle	4-0	Blue	1 × 36"/90cm	36		8519H
	4-0	Blue	1 × 30"/75cm	36		8881H
BB 17mm 3/8 Circle	5-0	Blue	1 × 24"/60cm	36		8860H
	5-0	Blue	1 × 36"/90cm	36		8580H
BB 17mm 3/8 Circle	4-0	Blue	1 × 24"/60cm	36		8861H
	4-0	Blue	1 × 36"/90cm	36		8581H
	4-0	Blue	1 × 36"/90cm	36	VB	8381H
	3-0	Blue	1 × 24"/60cm	36		8862H
	7-0	Blue	1 × 24"/60cm	36		8206H
RB-3 11mm 1/2 Circle	6-0	Blue	1 × 18"/45cm	36		8714H
	6-0	Blue	1 × 18"/45cm	36	EVP	EP8714H
	6-0	Blue	1 × 30"/75cm	36		8711H
RB-2 13mm 1/2 Circle	6-0	Blue	1 × 30"/75cm	36	EVP	EP8711H
	6-0	Blue	4 × 30"/75cm	12		M8711
	6-0	Blue	4 × 30"/75cm	12	EVP	EPM8711
	5-0	Blue	1 × 18"/45cm	36		8713H
	5-0	Blue	1 × 24"/60cm	36		8554H
	5-0	Blue	1 × 30"/75cm	36		8710H
	5-0	Blue	1 × 30"/75cm	36	HS EVP	EPH8710H
	5-0	Blue	4 × 30"/75cm	12		M8710
	5-0	Blue	1 × 36"/90cm	36	HS EVP	EPH8716H







PROLENE™ Polypropylene Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Blue	1 × 24"/60cm	36		8205ZH
TF 13mm 1/2 Circle	4-0	Blue	1 × 24"/60cm	36		8204H
	5-0	Blue	1 × 30"/75cm	36		8870H
RB-1 17mm 1/2 Circle	4-0	Blue	1 × 30"/75cm	36		8871H
	3-0	Blue	1 × 30"/75cm	36		8872H
	5-0	Blue	1 × 18"/45cm	36		8756H
RB-1 17mm 1/2 Circle	5-0	Blue	1 × 24"/60cm	36		85550H
	5-0	Blue	1 × 30"/75cm	36	HS	HS6856H
	5-0	Blue	1 × 36"/90cm	36		8556H
	5-0	Blue	1 × 36"/90cm	36	VB	8356H
	5-0	Blue	1 × 36"/90cm	36	HS	HS8556H
	5-0	Blue	4 × 36"/90cm	12		M8556
	4-0	Blue	1 × 18"/45cm	36		8757H
	4-0	Blue	1 × 36"/90cm	36		8557H
	4-0	Blue	1 × 36"/90cm	36	VB	8357H
	4-0	Blue	1 × 36"/90cm	36	HS	HS6857H
	4-0	Blue	4 × 36"/90cm	12		M8557
	3-0	Blue	1 × 36"/90cm	36		8558H
	3-0	Blue	1 × 36"/90cm	36	HS	9556H
	2-0	Blue	1 × 36"/90cm	36		8559H
	4-0	Blue	1 × 36"/90cm	36	HS	HS6861H
SH-2 20mm 1/2 Circle						
	0	Blue	8 × 18"/45cm	12	CR	C841G
MO-7 22mm 1/2 Circle						


PROLENE™ Polypropylene Suture









Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SH-1 22mm 1/2 Circle</p>	4-0	Blue	1 × 30"/75cm	36		8526H
	3-0	Blue	1 × 30"/75cm	36		8762H
 <p>CT-2 26mm 1/2 Circle</p>	3-0	Blue	1 × 30"/75cm	36		8410H
	2-0	Blue	8 × 18"/45cm	12	CR	C826G
	2-0	Blue	1 × 30"/75cm	36		8411H
	0	Blue	8 × 18"/45cm	12	CR	C827G
	0	Blue	1 × 30"/75cm	36		8412H
 <p>MO-6 26mm 1/2 Circle</p>	2-0	Blue	1 × 30"/75cm	36		8417H
	0	Blue	8 × 18"/45cm	12	CR	C845G
	0	Blue	1 × 30"/75cm	36		8418H
 <p>SH 26mm 1/2 Circle</p>	4-0	Blue	1 × 30"/75cm	36		8831H
	3-0	Blue	1 × 30"/75cm	36		8832H
	2-0	Blue	1 × 30"/75cm	36		8833H
	0	Blue	1 × 30"/75cm	36		8834H
 <p>SH 26mm 1/2 Circle</p>	3-0	Blue	1 × 36"/90cm	36	HS	HS6822H
	4-0	Blue	1 × 36"/90cm	36		8521H
	4-0	Blue	1 × 36"/90cm	36	HS	HS6855H
	3-0	Blue	1 × 36"/90cm	36		8522H
	3-0	Blue	1 × 36"/90cm	36	Pledget*	P8522H
	3-0	Blue	1 × 48"/122cm	36		8534H
	2-0	Blue	1 × 36"/90cm	36		8523H
	2-0	Blue	1 × 48"/122cm	36		8533H


*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget

PROLENE™ Polypropylene Suture





Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	3-0	Blue	1 × 30"/75cm	36		8422H
	2-0	Blue	8 × 18"/45cm	12	CR	C822G
	2-0	Blue	1 × 30"/75cm	36		8423H
	0	Blue	8 × 18"/45cm	12	CR	C821G
	0	Blue	1 × 30"/75cm	36		8424H
	1	Blue	1 × 30"/75cm	36		8425H
 <p>MH 36mm 1/2 Circle</p>	3-0	Blue	1 × 36"/90cm	36		8842H
	2-0	Blue	1 × 36"/90cm	36		8843H
	2-0	Blue	1 × 48"/122cm	36		8853H
 <p>CT 40mm 1/2 Circle</p>	2-0	Blue	1 × 30"/75cm	36		8433H
	0	Blue	1 × 30"/75cm	36		8434H
	1	Blue	1 × 30"/75cm	36		8435H
 <p>CTX 48mm 1/2 Circle</p>	0	Blue	8 × 18"/45cm	12	CR	C811G
	0	Blue	1 × 30"/75cm	36		8454H
	1	Blue	1 × 30"/75cm	36		8455H
 <p>TP-1 65mm 1/2 Circle</p>	1	Blue	1 × 60"/150cm	12		8824G
	2	Blue	1 × 60"/150cm	12		8825G






PROLENE™ Polypropylene Suture




Taperpoint




























Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 XLH 70mm 1/2 Circle	0	Blue	3 × 24"/60cm	12		3844G
	1	Blue	1 × 60"/150cm	12		8845G
 XLH 70mm 1/2 Circle	1	Blue	1 × 30"/75cm	12		8190G

TAPERCUT™


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CIF-4 13mm 1/4 Circle	10-0	Blue	1 × 8"/20cm	12		788G
 CC 9.3mm 3/8 Circle	7-0	Blue	1 × 18"/45cm	36		8705H
	7-0	Blue	4 × 18"/45cm	12		M8705
	7-0	Blue	4 × 18"/45cm	12	 	EPM8705


PROLENE™ Polypropylene Suture




TAPERCUT™ 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CC 9.3mm 3/8 Circle</p>	7-0	Blue	1 × 24"/60cm	36		8704H
	7-0	Blue	4 × 24"/60cm	12		M8704
	7-0	Blue	4 × 24"/60cm	12	 	EPM8704
	7-0	Blue	1 × 24"	36		EP8704H
	7-0	Blue	1 × 18"	36		EP8705H
	7-0	Blue	4×18"/45cm	12		EPM8705
	6-0	Blue	1 × 24"/60cm	36		8807H
	6-0	Blue	1 × 24"/60cm	36	 	EP8807H
	6-0	Blue	4 × 24"/60cm	12		M8807
	6-0	Blue	4 × 24"/60cm	12	 	EPM8807
 <p>CC-1 13mm 3/8 Circle</p>	6-0	Blue	1 × 18"/45cm	36		8719H
	6-0	Blue	1 × 18"/45cm	36		EP8719H
	6-0	Blue	4 × 18"/45cm	12		M8719
	6-0	Blue	4 × 18"/45cm	12	 	EPM8719
	6-0	Blue	1 × 30"/75cm	36		8707H
	6-0	Blue	1 × 30"/75cm	36	 	EP8707H
	6-0	Blue	4 × 30"/75cm	12		M8707
	6-0	Blue	4 × 30"/75cm	12	 	EPM8707
	6-0	Blue	2×30"/75cm	12		EPM8723
	5-0	Blue	1 × 36"/90cm	36		8721H
 <p>V-4 17mm 3/8 Circle</p>	5-0	Blue	4 × 36"/90cm	12		M8721
	4-0	Blue	1 × 36"/90cm	36		8925H
	3-0	Blue	1 × 36"/90cm	36		8926H

PROLENE™ Polypropylene Suture







TAPERCUT™ 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
CC175-8 9.3mm 3/8 Circle  V-5 17mm 1/2 Circle	7-0	Blue	4 × 24"/60cm	12		EPM8704SL
	7-0	Blue	4 × 30"/75cm	12		EPM8812SL
 V-7 26mm 1/2 Circle  V-34 36mm 1/2 Circle	5-0	Blue	1 × 36"/90cm	36		8934H
	4-0	Blue	1 × 36"/90cm	36		8935H
	4-0	Blue	1 × 36"/90cm	36	Pledget*	P8935H
	3-0	Blue	1 × 36"/90cm	36		8936H
	2-0	Blue	1 × 36"/90cm	36		8937H
 V-47 36mm Straight	4-0	Blue	1 × 36"/90cm	36		8975H
	3-0	Blue	1 × 36"/90cm	36		8976H
	2-0	Blue	1 × 36"/90cm	36		8977H
	0	Blue	1 × 36"/90cm	36		8978H
	0	Blue	1 × 30"/75cm	36		8444H
	4-0	Blue	1 × 14"/35cm	12		8592G



Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	5-0	Blue	1 × 18"/45cm	12		8661G
	4-0	Blue	1 × 18"/45cm	12		8683G
	3-0	Blue	1 × 18"/45cm	12		8665G

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget


PROLENE™ Polypropylene Suture



























Reverse Cutting ▾							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 FS-1 24mm 3/8 Circle	4-0	Blue	1 × 18"/45cm	36		8629H	
	3-0	Blue	1 × 18"/45cm	12		8684G	
 FS 26mm 3/8 Circle	2-0	Blue	1 × 18"/45cm	36		8685H	
 FSL 30mm 3/8 Circle	3-0	Blue	1 × 30"/75cm	36		8675H	
 FSLX 36mm 3/8 Circle	3-0	Blue	1 × 30"/75cm	36		8649H	
	2-0	Blue	1 × 30"/75cm	36		8689H	
	0	Blue	1 × 30"/75cm	36		8690H	
 LR 75mm 3/8 Circle	2	Blue	3 × 20"/50cm	24		3846T	
 STC-6 16mm Straight	10-0	Blue	1 × 8"/20cm	12		1713G	


Straight Cutting ▾							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 TS 48mm Straight	2-0	Blue	1 × 14"/35cm	12		8571G	
 KS 60mm Straight	4-0	Blue	1 × 30"/75cm	36		8621H	
	3-0	Blue	1 × 30"/75cm	36		8622H	
	2-0	Blue	1 × 30"/75cm	36		8623H	


PROLENE™ Polypropylene Suture








PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-6 8mm 1/4 Circle	7-0	Blue	1 × 18"/45cm	12		8648G
 P-1 11mm 3/8 Circle	7-0	Blue	1 × 18"/45cm	12		8696G
	6-0	Undyed	1 × 18"/45cm	12		8606G
	6-0	Blue	1 × 18"/45cm	12		8697G
 P-3 13mm 3/8 Circle	6-0	Blue	1 × 18"/45cm	12		8695G
	5-0	Undyed	1 × 18"/45cm	12		8605G
	5-0	Blue	1 × 18"/45cm	12		8698G
	4-0	Undyed	1 × 18"/45cm	12		8604G
	4-0	Blue	1 × 18"/45cm	36		8699G
 PS-3 16mm 3/8 Circle	6-0	Blue	1 × 18"/45cm	12		8680G
	5-0	Blue	1 × 18"/45cm	12		8681G
 PS-2 19mm 3/8 Circle	5-0	Blue	1 × 18"/45cm	12		8686G
	4-0	Blue	1 × 18"/45cm	12		8682G
	4-0	Blue	1 × 18"/45cm	36		8682H
	3-0	Blue	1 × 18"/45cm	36		8687H
 PS-1 24mm 3/8 Circle	3-0	Blue	1 × 18"/45cm	12		8663G
 PS-5 13mm 1/2 Circle	5-0	Blue	1 × 18"/45cm	12		8655G
	4-0	Blue	1 × 18"/45cm	12		8656G
 PS-4 16mm 1/2 Circle	4-0	Undyed	1 × 18"/45cm	12		8603G

PROLENE™ Polypropylene Suture




PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>PC-1 13mm 3/8 Circle</p>	6-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8617G
	5-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8618G
	4-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8619G
 <p>PC-3 16mm 3/8 Circle</p>	6-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8636G
	5-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8635G
	4-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8634G
 <p>PC-5 19mm 3/8 Circle</p>	5-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8630G
	4-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8631G
	3-0	Blue	1 × 18"/45cm	12	 MultiPass [®]	8632G






Spatula Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SM-1 6mm 1/8 Circle</p>	5-0	Blue	1 × 12"/30cm	12		7740G
 <p>TG140-8 6.5mm 3/8 Circle</p>	10-0	Blue	1 × 6"/15cm	12		1757G
	10-0	Blue	1 × 12"/30cm	12		1771G
	9-0	Blue	1 × 6"/15cm	12		1754G

ULTIMA™ Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTC-6L 13mm 1/4 Circle</p>	10-0	Blue	1 × 4"/10cm	12		9091G

PROLENE™ Polypropylene Suture



ULTIMA™ Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTC-6 13mm/ CS160-6 5.5mm 1/4 Circle/1/2 Circle	10-0	Blue	1 × 4"/10cm	12		9090G
 CS140-6 6.5mm 3/8 Circle	10-0	Blue	1 × 12"/30cm	12		1710G
 CS160-6 5.5mm 1/2 Circle	10-0	Blue	1 × 12"/30cm	12		7761G
 CS160-8 5.5mm 1/2 Circle	10-0	Blue	1 × 6"/15cm	12		1755G
 CSC-6 5mm Compound Curve	10-0	Blue	1 × 12"/30cm	12		7796G


Surgical Stainless Steel Suture

Surgical stainless steel suture is a nonabsorbable, surgical suture composed of 316L stainless steel. The suture is available as a monofilament and multifilament suture.^{1,2}

Surgical stainless steel suture is indicated for use in abdominal wound closure, hernia repair, sternal closure and orthopaedic procedures including cerclage and tendon repair.^{1,2}






Surgical Stainless Steel Suture




Taperpoint

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Silver	8 × 18"/45cm	12		M682G
 <p>CTX 48mm 1/2 Circle</p>	2-0	Silver	8 × 24"/60cm	12		M414G







TAPERCUT™

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-40 48mm 1/2 Circle</p>	4	Silver	4 × 18"/45cm	12		M651G
	5	Silver	4 × 18"/45cm	12		M650G
	5	Silver	2 × 18"/45cm	12		M656G
	6	Silver	4 × 18"/45cm	12		M649G
 <p>V-56 60mm 1/2 Circle</p>	5	Silver	4 × 30"/75cm	12		M435G
 <p>V-37 40mm 1/2 Circle</p>	1	Silver	4 × 18"/45cm	12		M660G


Surgical Stainless Steel Suture



Reverse Cutting ▾

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	4-0	Silver	1 × 18"/45cm	12		603G
	3-0	Silver	1 × 18"/45cm	12		607G
 FS-1 24mm 3/8 Circle	3-0	Silver	1 × 18"/45cm	12		631G
 FSLX 36mm 3/8 Circle	2-0	Silver	1 × 18"/45cm	12		618G
 LS-1 45mm 3/8 Circle	2	Silver	4 × 20"/50cm	12		M400G
 LR 75mm 3/8 Circle	2	Silver	1 × 30"/75cm	12		467G
 CP-1 36mm 1/2 Circle	3-0	Silver	1 × 20"/50cm	12		497G





Surgical Stainless Steel Suture



Cutting Edge Conventional


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CCS 48mm 1/2 Circle	4	Silver	4 × 18"/45cm	12		M652G
	5	Silver	2 × 18"/45cm	12		M657G
	5	Silver	4 × 18"/45cm	12		M653G
	6	Silver	4 × 18"/45cm	12		M654G
	7	Silver	4 × 18"/45cm	12		M655G
 CCS-1 60mm 1/2 Circle	5	Silver	2 × 30"/75cm	12		M436G

Taperpoint , Straight Cutting Temporary Pacing Wire



Needle	Suture USP	Color	Length	EA/BX	Product Code
 BB 17mm/ SKS 60mm 3/8 Circle/Straight	2-0	Dark Blue	1 × 24"/60cm	12	TPW20*
 BB-1 22mm/SKS-3 90mm 3/8 Circle/Straight	2-0	Dark Blue	1 × 24"/60cm	12	TPW40**
 BB-1 22mm/SKS-3 90mm 3/8 Circle/Straight	2-0	Light Blue/Dark Blue	2 × 24"/60cm	12	TPW42**†
 RB-1 17mm/ SKS 60mm 1/2 Circle/Straight	2-0	Dark Blue	1 × 24"/60cm	12	TPW10*

*One dozen packets per box †1.5" long score marks; others .875" ‡Multistrand 2"- 60 cm


Surgical Stainless Steel Suture



Taperpoint  , Straight Cutting Temporary Pacing Wire 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>SH 26mm/ SKS-3 90mm 1/2 Circle/Straight</p>  <p>TE 32mm/ SKS-3 90mm 3/8 Circle/Straight</p>	2-0	Dark Blue	1 × 24”/ 60cm	12	TPW30* [†]
	2-0	Dark Blue	1 × 24”/ 60cm	12	TPW90* ^{†§}
	2-0	Light Blue/Dark Blue	2 × 24”/ 60cm	12	TPW32* ^{†‡}
	2-0	Light Blue/Dark Blue	2 × 24”/ 60cm	12	TPW72* ^{†‡}
	2-0	Light Blue/Dark Blue	2 × 24”/ 60cm	12	TPW92* ^{†‡§}
	2-0	Dark Blue	1 × 24”/ 60cm	12	TPW80* ^{†§}

Straight Cutting 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>KS 60mm Straight</p>	5-0	Silver	8 × 16”/40cm	12		M645G

SUTUPAK™ Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Silver	17 × 18”/45cm	12		SS30G
	2-0	Silver	17 × 18”/45cm	12		SS28G
	0	Silver	17 × 18”/45cm	12		SS26G

*One dozen packets per box †1.5” long score marks; others .875” ‡Multistrand 2”- 60 cm §TP wave codes

Surgical Stainless Steel Suture




ETHI-PACK™ Pre-cut Sutures (Non-sterile)

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Silver	50 × 18"/45cm	12		DS35
	4-0	Silver	50 × 18"/45cm	12		DS32
	3-0	Silver	50 × 18"/45cm	12		DS30
	2-0	Silver	50 × 18"/45cm	12		DS28
	0	Silver	50 × 18"/45cm	12		DS26
	1	Silver	50 × 18"/45cm	12		DS25
	2	Silver	12 × 18"/45cm	12		DS24
	3	Silver	12 × 18"/45cm	12		DS23
	4	Silver	12 × 18"/45cm	12		DS22
	5	Silver	12 × 18"/45cm	12		DS20
	7	Silver	12 × 18"/45cm	12		DS18

Primary Suture Storage Carts, Carrying Racks, and Other Wound Closure Devices


Other Wound Closure Devices

ENDOSUTURE ENDO-HOLDER™



Product Code	Description	EA/BX
EC11	EN-3 Needle - O, 42" ETHIBOND Suture	12
EC12	EN-3 Needle - 2-0, 42" ETHIBOND Suture	12
ESS15	ENDO-HOLDER Cannula (for use with 5mm diameter trocar)	1

ENDOKNOT™ Pre-Tied Suture



Product Code	Description	EA/BX
EX10G	ETHIBOND™ Extra Polyester Suture, O, 42"	12
JK10G	Coated VICRYL™ (polyglactin 910) Suture, ST-3 needle, O, 42"	12
ZK10G	PDS™ II (polydioxanone) Suture, ST-3 needle, O, 42"	12

ENDOLOOP™ Ligature made from PDS™ II (polydioxanone) Suture Material



Product Code	Description	EA/BX
EZ10G	ENDOLOOP Ligature PDS II Suture, O, 18"	12

ENDOLOOP™ Ligature made from VICRYL™ (polyglactin 910) Suture Material



Product Code	Description	EA/BX
EJ10G	ENDOLOOP Ligature VICRYL Suture, O, 18"	12

ENDOPATH™ Needle Holder



Product Code	Description	EA/BX
E705R	ENDOPATH Needle Holder, reusable, 5mm, 30cm length	1
E705RL	ENDOPATH Needle Holder, reusable, 5mm, 45cm length	1
SRNHI	Self Righting Needle Holder, reusable, 5mm, 30cm length	1
SRNHIL	Self Righting Needle Holder, reusable, 5mm, 45cm length	1

Other Wound Closure Devices

LAPRA-TY™ Suture Clip



Product Code	Description	EA/BX
KA200	Lapra-TY Suture Clip Applier.	1
XC200	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	6

PROXIMATE™ Rotating Head Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PRR35	35	3.9	0.53	5.7	Regular	6
PRW35	35	3.9	0.58	6.9	White	6

PROXIMATE™ Plus MD Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PMR35	35	3.9	0.53	5.7	Regular	6
PMW35	35	3.9	0.58	6.9	Wide	6

Other Wound Closure Devices

PROXIMATE™ PX Fixed Head
Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PXR35	35	3.9	0.53	5.7	Regular	6
PXW35	35	3.9	0.58	6.9	Wide	6

Staple Extractor



Product Code	Number of Staples	EA/BX
PSX	Squeeze handle staple extractor	12

Bone Wax



Product Code	Description
W31G	2.5 Grams, 1 Dozen per Box

TFE Polymer Pledgets



Product Code	Description
PCP20	TFE Pledget, White, Firm. 7.0mm x 3.5mm x 1.5mm. 36 envelopes/box (216 pledgets).
PCP30	TFE Pledget, White, Firm. 9.5mm x 4.8mm x 1.5mm. 36 envelopes/box (216 pledgets).
PCP40	TFE Pledget, White, Soft. 9.5mm x 4.8mm x 1.5mm. 36 envelopes/box (216 pledgets).
PCP50	TFE Pledget, White, Soft. 7.0mm x 3.5mm x 1.5mm. 36 envelopes/box (216 pledgets).


Other Wound Closure Devices

Suture Cart




Product Code	Description
CART	38"L x 18"W x 38"H. No suture racks included. Brackets to attach racks included. Some assembly required.

3-Wide Lazy Susan




Product Code	Description
LS3	15 3/8"L X 11 1/2"W X 20 3/4"H. Includes 6 modular racks. Can be expanded with additional modular racks.

4-Wide Lazy Susan




Product Code	Description
LS4	20 1/2"L X 11 1/2"W X 20 3/4"H. Includes 8 modular racks. Can be expanded with additional modular racks.

Dispenser Rack



Product Code	Description
MR20	25 3/4"L X 5 3/4"W X 27"H. Includes 20 modules, 10 back plates (each with 4 fasteners), 20 module clips, 30 divider panels.


Dispenser Rack



Product Code	Description
MR4	10 1/4"L X 5 3/4"W X 13 1/2"H. Includes 4 modules, 3 back plates (each with 4 fasteners), 4 module clips, 8 divider panels.

Other Wound Closure Devices

Portable Carrying Rack




Product Code	Description
PCR1	15 3/8"L x 6 1/8"W x 17 5/8"H. Includes 6 modular racks.

2-Sided Portable Carrying Rack



Product Code	Description
PCR2	15 3/8"L x 12 1/4"W x 16 3/4"H. Includes 12 modular racks.


I.V. Pole with Rack











Product Code	Description
SPR1	18 1/4"L x 16 1/2"W x 56 1/2"H. Includes 18 modules racks, 15 module clips, 30 divider panels. Includes 18 modules racks, 15 module clips, 30 divider panels.

Suture Chart

Absorbable Sutures: The Strength You Need for the Choices You Make



Nonabsorbable Sutures: The Strength You Need for the Choices You Make

	SUTURE BRAND	MATERIAL	SIZES	FREQUENT USES	BSR† PROFILE	
NATURAL	PERMA-HAND® Silk Suture	Braided Black	5, 2 through 7/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Gradual loss of all tensile strength over time	
	SURGICAL STAINLESS STEEL Suture	Monofilament Silver Metallic	7 through 5/0	Soft tissue approximation: n Abdominal n Sternal Closure n Orthopedic procedures including cerclage and tendon repair	Indefinite	
SYNTHETIC	NUROLON® Nylon Suture	Braided Black	1 through 6/0	General soft tissue approximation and/or ligation: n Neurological n Cardiovascular n Ophthalmic	Gradual loss of tensile strength over time	
	ETHILON® Nylon Suture	Monofilament Black/Green/ Undyed (Clear)	2 through 11/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological n Skin Closure	Gradual loss of tensile strength over time	
	MERSILENE® Polyester Fiber Suture	Braided Green / Undyed (white) Monofilament Green	1 through 6/0 10/0 through 11/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	
	ETHIBOND EXCEL® Polyester Suture	Braided Green / Undyed (white)	5, 2 through 7/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	
	PROLENE® Polypropylene Suture	Monofilament Blue / Clear	2 through 10/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	
	PRONOVA® Poly (Hexafluoropropylene-VDF) Suture	Monofilament Blue	2/0 through 8/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	

† BSR=breaking strength retention (tensile strength of suture in vivo)

Ethicon Wound Closure Catalog – Index

Absorbable Sutures Chart	215	PCR2.....	213	SXPP2B405.....	45
Contact and More Information	231	PMR35.....	210	SXPP2B407.....	44
Explanation of a Typical Selection and Page	5	PMW35.....	210	SXPP2B408.....	44
General Information.....	4	PRR35.....	210	SXPP2B409.....	44
Nonabsorbable Sutures Chart	216	PRW35.....	210	SXPP2B410.....	44
Selection and Use of Needleholders	7	PSX.....	211	SXPP2B411.....	44
Shapes	8	PXR35.....	211	SXPP2B412.....	44
Steps to Create Performance	8	PXW35.....	211	SXPP2B413.....	43
Table of Contents	2	SPR1.....	213	SXPP2B414.....	43
Trademarks.....	10	SRNH1.....	209	SXPP2B415.....	46
Types	9	SRNHIL.....	209	SXPP2B417.....	46
DERMABOND™ Portfolio.....	34-35	W31G.....	211	SXPP2B418.....	46
DERMABOND ADVANCED™ Topical Skin Adhesive.....	35	XC200.....	210	SXPP2B419.....	46
DNX6.....	35	ZK10G.....	209	SXPP2B420.....	46
DNX12.....	35	STRATAFIX™ Knotless Tissue Control Devices	36-55	SXPP2B423.....	46
DERMABOND™ Mini Topical Skin Adhesive.....	35	STRATAFIX™ Spiral Polypropylene Knotless Tissue Control Device.....	38	SXPP2B424.....	47
DHVM12.....	35	SXPL1B400.....	38	SXPP2B425.....	47
DERMABOND™ PRINEO™ Skin Closure System.....	35	SXPL2B400.....	38	SXPP2B426.....	47
CLR222US.....	35	SXPL2B401.....	38	SXPP2B427.....	43
CLR422US.....	35	Bidirectional STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device	38-40	SXPP2B428.....	43
CLR602US.....	35	SXMP2B150.....	40	SXPP2B430.....	43
Needle Index.....	12-28	SXMP2B151.....	40	SXPP2B431.....	43
1/4 Circle.....	12-13	SXMP2B152.....	40	SXPP2B432.....	45
3/8 Circle.....	13-17	SXMP2B400.....	39	SXPP2B433.....	45
1/2 Circle.....	18-24	SXMP2B401.....	39	SXPP2B434.....	44
5/8 Circle.....	24	SXMP2B402.....	38	SXPP2B435.....	45
Bi-Curve.....	27	SXMP2B403.....	38	SXPP2B436.....	45
Center Point Spatula.....	28	SXMP2B404.....	40	Unidirectional STRATAFIX™ PDS™ Plus Knotless Tissue Control Device	48-51
Compound Curve.....	27	SXMP2B405.....	40	SXPP1B101.....	50
Endoscopy.....	28	SXMP2B406.....	40	SXPP1B102.....	50
Sternum.....	27	SXMP2B407.....	40	SXPP1B103.....	51
Straight.....	25	SXMP2B408.....	40	SXPP1B104.....	51
Trocar.....	26	SXMP2B409.....	40	SXPP1B105.....	51
Primary Suture Storage Carts, Carrying Racks, and Other Wound Closure Devices	208-213	SXMP2B410.....	40	SXPP1B106.....	51
CART.....	212	SXMP2B411.....	40	SXPP1B107.....	51
E705R.....	209	SXMP2B412.....	40	SXPP1B108.....	51
E705RL.....	209	SXMP2B413.....	40	SXPP1B109.....	51
EC11.....	209	SXMP2B414.....	40	SXPP1B110.....	51
EC12.....	209	SXMP2B415.....	38	SXPP1B111.....	51
EJ1OG.....	209	SXMP2B416.....	38	SXPP1B112.....	51
ESS15.....	209	SXMP2B417.....	39	SXPP1B113.....	51
EX1OG.....	209	SXMP2B418.....	39	SXPP1B114.....	51
EZ1OG.....	209	SXMP2B419.....	39	SXPP1B116.....	51
JK10G.....	209	SXMP2B420.....	39	SXPP1B117.....	51
KA200.....	210	SXMP2B421.....	39	SXPP1B118.....	51
LS3.....	212	SXMP2B422.....	38	SXPP1B119.....	51
LS4.....	212	SXMP2B423.....	39	SXPP1B120.....	50
MR20.....	212	SXMP2B424.....	39	SXPP1B201.....	50
MR4.....	212	SXMP2B427.....	39	SXPP1B202.....	50
PCP20.....	211	Bidirectional STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device	43-47	SXPP1B203.....	50
PCP30.....	211	SXPP2B100.....	46	SXPP1B204.....	50
PCP40.....	211	SXPP2B102.....	46	SXPP1B205.....	49
PCP50.....	211	SXPP2B200.....	46	SXPP1B401.....	49
PCR1.....	213	SXPP2B201.....	46	SXPP1B402.....	49
		SXPP2B202.....	47	SXPP1B403.....	50
		SXPP2B203.....	47	SXPP1B404.....	49
		SXPP2B300.....	47	SXPP1B405.....	48
		SXPP2B301.....	47	SXPP1B406.....	48
		SXPP2B400.....	44	SXPP1B407.....	48
		SXPP2B401.....	44	SXPP1B408.....	48
		SXPP2B402.....	44	SXPP1B409.....	48
		SXPP2B403.....	44	SXPP1B410.....	48
		SXPP2B404.....	44	SXPP1B411.....	48
				SXPP1B412.....	48
				SXPP1B413.....	48
				SXPP1B414.....	48
				SXPP1B415.....	48

SXPP1B416.....	48	SXMP1B438.....	41	SXPP1A437.....	54
SXPP1B417.....	48	SXMP1B440.....	41	SXPP1A438.....	53
SXPP1B418.....	49	SXMP1B441.....	41	SXPP1A439.....	53
SXPP1B419.....	49	SXMP1B442.....	41	SXPP1A440.....	53
SXPP1B420.....	48	SXMP1B443.....	41	SXPP1A441.....	53
SXPP1B421.....	48	SXMP1B444.....	41	SXPP1A442.....	53
SXPP1B422.....	48	SXMP1B445.....	41	SXPP1A443.....	53
SXPP1B424.....	49	SXMP1B447.....	41	SXPP1A444.....	53
SXPP1B425.....	48	SXMP1B448.....	42	SXPP1A445.....	53
SXPP1B426.....	48	SXMP1B449.....	41	SXPP1A446.....	52
SXPP1B427.....	48	SXMP1B450.....	41	SXPP1A447.....	53
SXPP1B428.....	48	SXMP1B451.....	41		
SXPP1B429.....	48	SXMP1B452.....	41	Sutures: Absorbable	56-133
SXPP1B430.....	48	SXMP1B453.....	42	Coated VICRYL™ Plus Antibacterial	
SXPP1B431.....	48	SXMP1B454.....	42	(polyglactin 910) Suture	61-92
SXPP1B432.....	48	SXMP1B455.....	42	J104G.....	90
SXPP1B433.....	48	Unidirectional STRATAFIX™ Symmetric PDS™		J106T.....	91
SXPP1B450.....	48	Plus Knotless Tissue Control Device	52-55	J110T.....	90
SXPP1B452.....	49	SXPP1A100.....	55	J113T.....	91
SXPP1B453.....	48	SXPP1A101.....	55	J195H.....	80
SXPP1B454.....	48	SXPP1A200.....	54	J196H.....	80
SXPP1B455.....	48	SXPP1A201.....	54	J204G.....	92
SXPP1B456.....	48	SXPP1A202.....	54	J206G.....	92
SXPP1B457.....	49	SXPP1A203.....	54	J212H.....	64
SXPP1B458.....	49	SXPP1A204.....	54	J214H.....	65
		SXPP1A205.....	54	J219H.....	66
		SXPP1A205.....	54	J232H.....	66
		SXPP1A300.....	52	J258H.....	69
		SXPP1A301.....	52	J259H.....	69
		SXPP1A302.....	52	J260H.....	70
		SXPP1A303.....	52	J261H.....	71
		SXPP1A304.....	52	J267H.....	79
		SXPP1A305.....	52	J268H.....	79
		SXPP1A306.....	52	J269H.....	66
		SXPP1A400.....	53	J275H.....	72
		SXPP1A401.....	53	J281H.....	73
		SXPP1A402.....	53	J283G.....	92
		SXPP1A403.....	53	J284G.....	92
		SXPP1A404.....	53	J285G.....	92
		SXPP1A405.....	53	J304H.....	65
		SXPP1A406.....	53	J306H.....	65
		SXPP1A407.....	52	J310H.....	66
		SXPP1A408.....	52	J315H.....	67
		SXPP1A409.....	52	J316H.....	67
		SXPP1A410.....	52	J317H.....	68
		SXPP1A411.....	52	J322H.....	71
		SXPP1A412.....	53	J323H.....	71
		SXPP1A413.....	52	J327H.....	65
		SXPP1A414.....	52	J328H.....	65
		SXPP1A415.....	52	J329H.....	65
		SXPP1A416.....	53	J333H.....	66
		SXPP1A417.....	53	J335H.....	67
		SXPP1A418.....	53	J337H.....	69
		SXPP1A419.....	52	J338H.....	69
		SXPP1A420.....	53	J339H.....	70
		SXPP1A421.....	52	J340H.....	70
		SXPP1A422.....	52	J344H.....	69
		SXPP1A423.....	53	J345H.....	70
		SXPP1A424.....	54	J346H.....	71
		SXPP1A425.....	53	J347H.....	71
		SXPP1A426.....	54	J352H.....	73
		SXPP1A427.....	52	J353H.....	73
		SXPP1A428.....	52	J359H.....	73
		SXPP1A429.....	53	J363H.....	74
		SXPP1A430.....	52	J364H.....	74
		SXPP1A431.....	52	J365H.....	74
		SXPP1A433.....	53	J369H.....	74
		SXPP1A434.....	53	J370H.....	74
		SXPP1A435.....	53	J371H.....	74
		SXPP1A436.....	53		

J381H.....	75	J566G.....	89	J912G.....	91
J385H.....	76	J568H.....	81	J913G.....	91
J386H.....	76	J569H.....	81	J916H.....	76
J392H.....	76	J570G.....	88	J917H.....	76
J401G.....	64	J571G.....	88	J935H.....	85
J405G.....	64	J574G.....	89	J944H.....	69
J408G.....	76	J575G.....	89	J945H.....	70
J409G.....	76	J576G.....	89	J946H.....	71
J415H.....	67	J578G.....	88	J947H.....	71
J417H.....	68	J579G.....	88	J958H.....	73
J418H.....	68	J581G.....	75	J959H.....	73
J421H.....	76	J582G.....	75	J974G.....	89
J428H.....	84	J583G.....	75	J975G.....	89
J433H.....	64	J590G.....	88	J977H.....	74
J436H.....	72	J591G.....	88	J978H.....	74
J437H.....	72	J599G.....	89	J979H.....	74
J446G.....	88	J606H.....	90	J980H.....	74
J451H.....	77	J610H.....	83	JB725.....	63
J455H.....	77	J614H.....	90	JB841H.....	63
J458H.....	78	J633H.....	90	JB945H.....	62
J459H.....	78	J634H.....	90	JB947H.....	63
J460H.....	78	J635H.....	90	JB977H.....	63
J461H.....	78	J636H.....	91	JB978H.....	63
J463G.....	83	J643H.....	90	VCP103G.....	90
J464G.....	83	J644H.....	90	VCP104G.....	90
J467H.....	79	J645H.....	90	VCP105G.....	90
J468H.....	79	J646H.....	91	VCP106G.....	91
J471H.....	79	J649G.....	75	VCP107G.....	91
J473H.....	79	J650G.....	75	VCP109G.....	90
J474H.....	79	J656G.....	86	VCP110G.....	90
J479H.....	80	J657G.....	86	VCP111G.....	90
J480H.....	80	J662H.....	82	VCP112G.....	91
J485H.....	80	J670G.....	88	VCP195H.....	80
J486H.....	80	J671G.....	88	VCP196H.....	80
J488G.....	83	J682H.....	85	VCP197H.....	80
J489G.....	83	J683H.....	85	VCP205G.....	92
J490G.....	83	J693H.....	77	VCP207G.....	92
J492G.....	83	J699H.....	81	VCP208G.....	92
J492H.....	83	J707T.....	78	VCP213H.....	64
J493G.....	83	J726D.....	66	VCP214H.....	65
J493H.....	83	J727D.....	67	VCP215H.....	65
J494H.....	83	J739D.....	69	VCP218H.....	66
J495G.....	84	J740D.....	70	VCP219H.....	66
J495H.....	84	J741D.....	71	VCP232H.....	66
J496G.....	84	J756T.....	81	VCP259H.....	70
J496H.....	84	J757T.....	81	VCP260H.....	70
J497G.....	84	J774D.....	67	VCP261H.....	71
J497H.....	84	J775D.....	68	VCP265H.....	79
J499G.....	84	J783G.....	67	VCP266H.....	79
J500G.....	84	J801D.....	69	VCP267H.....	79
J506G.....	86	J808T.....	86	VCP268H.....	79
J507G.....	86	J813G.....	87	VCP269H.....	66
J510G.....	85	J814G.....	87	VCP270H.....	67
J517H.....	76	J823G.....	87	VCP275H.....	72
J518H.....	76	J823H.....	87	VCP280H.....	73
J519H.....	76	J824G.....	87	VCP281H.....	73
J523H.....	82	J833G.....	87	VCP284G.....	92
J531H.....	82	J839D.....	69	VCP285G.....	92
J532H.....	82	J843G.....	87	VCP286G.....	92
J544G.....	88	J844G.....	87	VCP287G.....	92
J546G.....	89	J845G.....	87	VCP303H.....	64
J547G.....	89	J864D.....	67	VCP304H.....	65
J548G.....	89	J868H.....	78	VCP305H.....	65
J55.....	88	J870H.....	79	VCP306H.....	65
J551G.....	83	J871H.....	79	VCP310H.....	66
J553G.....	88	J885G.....	92	VCP311H.....	66
J555G.....	88	J906G.....	91	VCP315H.....	67
J556G.....	88	J907G.....	91	VCP316H.....	68
J562G.....	88	J910T.....	90	VCP317H.....	68

VCP322H.....	71	VCP494H.....	83	VCP726D.....	66
VCP323H.....	71	VCP495G.....	84	VCP727D.....	67
VCP327H.....	65	VCP495H.....	84	VCP738D.....	69
VCP328H.....	65	VCP496ZG.....	84	VCP739D.....	69
VCP329H.....	65	VCP496ZH.....	84	VCP740D.....	70
VCP332H.....	66	VCP497H.....	84	VCP741D.....	71
VCP333H.....	66	VCP500G.....	84	VCP743D.....	64
VCP334H.....	67	VCP500H.....	84	VCP748T.....	77
VCP335H.....	67	VCP502H.....	77	VCP5749T.....	78
VCP338H.....	69	VCP503G.....	85	VCP750T.....	80
VCP339H.....	70	VCP504G.....	85	VCP751D.....	72
VCP340H.....	70	VCP506G.....	86	VCP752D.....	73
VCP341H.....	71	VCP507G.....	86	VCP753D.....	73
VCP344H.....	69	VCP511G.....	85	VCP754T.....	80
VCP345H.....	70	VCP516H.....	76	VCP756T.....	81
VCP346H.....	71	VCP517H.....	76	VCP757T.....	81
VCP347H.....	71	VCP518H.....	76	VCP761D.....	78
VCP351H.....	72	VCP519H.....	76	VCP762D.....	78
VCP352H.....	73	VCP527H.....	68	VCP764D.....	74
VCP353H.....	73	VCP533H.....	80	VCP765D.....	74
VCP356H.....	72	VCP534H.....	80	VCP771D.....	66
VCP357H.....	73	VCP535H.....	80	VCP772D.....	66
VCP358H.....	73	VCP568H.....	81	VCP773D.....	67
VCP359H.....	73	VCP569H.....	81	VCP774D.....	67
VCP363H.....	74	VCP581G.....	75	VCP775D.....	68
VCP364H.....	74	VCP583G.....	75	VCP784D.....	68
VCP365H.....	74	VCP588H.....	77	VCP785D.....	68
VCP369H.....	74	VCP593G.....	86	VCP787D.....	77
VCP370H.....	74	VCP594G.....	86	VCP789D.....	68
VCP371H.....	74	VCP596H.....	77	VCP790D.....	78
VCP375H.....	75	VCP602H.....	75	VCP801D.....	69
VCP376H.....	75	VCP603H.....	75	VCP822G.....	87
VCP391H.....	76	VCP607H.....	90	VCP823G.....	87
VCP392ZH.....	76	VCP608H.....	91	VCP823H.....	87
VCP393H.....	77	VCP615H.....	90	VCP824G.....	87
VCP415H.....	67	VCP616H.....	91	VCP834G.....	87
VCP416H.....	67	VCP617H.....	91	VCP835G.....	87
VCP417H.....	68	VCP618H.....	91	VCP838D.....	69
VCP418H.....	68	VCP634T.....	90	VCP839D.....	69
VCP421H.....	76	VCP635T.....	90	VCP840D.....	70
VCP422H.....	76	VCP636T.....	91	VCP841D.....	71
VCP423H.....	77	VCP645T.....	90	VCP844G.....	87
VCP426H.....	84	VCP646T.....	91	VCP845G.....	87
VCP427H.....	84	VCP649G.....	75	VCP849G.....	75
VCP428H.....	84	VCP656G.....	86	VCP864D.....	67
VCP433H.....	64	VCP662H.....	82	VCP868H.....	78
VCP434H.....	64	VCP663H.....	82	VCP869H.....	79
VCP436H.....	72	VCP682G.....	85	VCP870H.....	79
VCP437H.....	72	VCP683G.....	85	VCP871H.....	79
VCP441H.....	77	VCP683H.....	85	VCP880T.....	75
VCP442H.....	77	VCP694H.....	78	VCP285G.....	92
VCP443H.....	77	VCP695H.....	78	VCP916H.....	76
VCP451H.....	77	VCP698H.....	81	VCP917H.....	76
VCP452H.....	77	VCP699H.....	81	VCP931H.....	82
VCP458H.....	78	VCP700D.....	72	VCP935H.....	85
VCP459H.....	78	VCP701D.....	72	VCP936H.....	85
VCP460H.....	78	VCP702D.....	72	VCP940H.....	82
VCP463H.....	83	VCP706T.....	77	VCP944H.....	69
VCP464G.....	83	VCP708T.....	78	VCP945H.....	70
VCP466H.....	79	VCP710T.....	80	VCP946H.....	71
VCP467H.....	79	VCP711T.....	80	VCP947H.....	71
VCP468H.....	79	VCP713D.....	65	VCP956H.....	72
VCP472H.....	79	VCP714D.....	65	VCP957H.....	72
VCP474H.....	79	VCP716T.....	81	VCP958H.....	73
VCP478H.....	80	VCP717T.....	81	VCP959H.....	73
VCP480H.....	80	VCP718T.....	81	VCP977H.....	74
VCP486H.....	80	VCP719T.....	81	VCP978H.....	74
VCP493G.....	83	VCP723D.....	74	VCP979H.....	74
VCP493H.....	83	VCP724D.....	74	VCP980H.....	74

VCPB259H.....	62	MCP497G.....	102	Y489G.....	101
VCPB260H.....	63	MCP513G.....	102	Y490G.....	101
VCPB269H.....	62	MCP604H.....	100	Y492G.....	101
VCPB340H.....	63	MCP605H.....	100	Y493G.....	101
VCPB341H.....	63	MCP682H.....	102	Y494G.....	102
VCPB416H.....	62	MCP822G.....	103	Y495G.....	102
VCPB417H.....	62	MCP823G.....	103	Y496G.....	102
VCPB41D.....	63	MCP834G.....	103	Y497G.....	102
VCPB650G.....	64	MCP835G.....	103	Y510G.....	102
VCPB724D.....	63	MCP844G.....	103	Y513G.....	102
VCPB725D.....	63	MCP845G.....	103	Y523H.....	101
VCPB726D.....	62	MCP935H.....	102	Y527H.....	97
VCPB727D.....	62	MCP936H.....	102	Y604H.....	100
VCPB739D.....	62	MCP944H.....	98	Y605H.....	100
VCPB740D.....	63	MCP945H.....	98	Y606H.....	100
VCPB741D.....	63	MCP946H.....	99	Y682H.....	102
VCPB765D.....	63	MCPB316H.....	96	Y701D.....	99
VCPB774D.....	62	MCPB346H.....	96	Y732H.....	97
VCPB839D.....	62	MCPB416H.....	96	Y738D.....	98
VCPB840D.....	62	MCPB944H.....	96	Y739D.....	98
VCPB841D.....	63	MCPB945H.....	96	Y740D.....	98
VCPB864D.....	62	MCPB946H.....	96	Y752D.....	99
VCPB944H.....	62	Y213H.....	97	Y762H.....	98
VCPB945H.....	62	Y214H.....	97	Y801D.....	98
VCPB946H.....	62	Y215H.....	97	Y814G.....	103
VCPB947H.....	63	Y219H.....	97	Y817H.....	100
VCPB958H.....	63	Y227H.....	98	Y822G.....	103
VCPB977H.....	63	Y266H.....	100	Y823G.....	103
VCPB978H.....	63	Y267H.....	100	Y833G.....	103
VCPJ490G.....	83	Y303H.....	97	Y834G.....	103
VCPP31D.....	71	Y304H.....	97	Y835G.....	103
VCPP40D.....	71	Y305H.....	97	Y844G.....	103
VCPP41D.....	70	Y315H.....	97	Y845G.....	103
VCPP42D.....	70	Y316H.....	97	Y872D.....	97
VCPP71D.....	72	Y318H.....	97	Y908G.....	103
VCPP80D.....	69	Y333H.....	98	Y909G.....	103
VCPP81D.....	69	Y334H.....	98	Y910G.....	103
VCPPB31D.....	63	Y335H.....	98	Y911G.....	103
		Y338H.....	98	Y917H.....	100
		Y340H.....	98	Y931H.....	101
MONOCRYL™ Plus Antibacterial		Y341H.....	99	Y932H.....	101
(poliglecaprone 25) Suture.....	95-103	Y344H.....	98	Y935H.....	102
MCP213H.....	97	Y345H.....	98	Y936H.....	102
MCP214H.....	97	Y346H.....	98	Y944H.....	98
MCP266H.....	100	Y347H.....	99	Y945H.....	98
MCP304H.....	97	Y350H.....	99	Y946H.....	98
MCP315H.....	97	Y351H.....	99	Y947H.....	99
MCP316H.....	97	Y352H.....	99	YB257.....	96
MCP317H.....	97	Y353H.....	99	YB258.....	96
MCP333H.....	98	Y357H.....	99	YB416.....	96
MCP340H.....	99	Y358H.....	99	YB944.....	96
MCP344H.....	98	Y359H.....	99	YB946.....	96
MCP345H.....	98	Y377H.....	100	YB947.....	96
MCP346H.....	99	Y378H.....	100	YB978H.....	96
MCP347H.....	99	Y397H.....	99	YY30G.....	99
MCP358H.....	99	Y398H.....	99	YY31G.....	99
MCP377H.....	100	Y399H.....	99	YY70G.....	99
MCP415H.....	97	Y415H.....	97	YY71G.....	99
MCP416H.....	97	Y416H.....	97		
MCP417H.....	97	Y417H.....	97		
MCP426H.....	102	Y426H.....	102	PDS™ Plus Antibacterial (polydioxanone)	
MCP427H.....	102	Y427H.....	102	Suture.....	104-114
MCP463G.....	102	Y432H.....	96	PDP114H.....	107
MCP489G.....	101	Y433H.....	96	PDP117H.....	107
MCP490G.....	101	Y434H.....	96	PDP123H.....	107
MCP492G.....	101	Y436H.....	99	PDP126H.....	107
MCP493G.....	101	Y437H.....	99	PDP127H.....	107
MCP494G.....	102	Y463G.....	101	PDP137H.....	109
MCP495G.....	102	Y464G.....	102	PDP148H.....	107
MCP496G.....	102			PDP149H.....	107

PDP303H.....	107	PDPB370T.....	106	Z833G.....	114
PDP304H.....	107	PDPB371T.....	106	Z844G.....	114
PDP305H.....	107	PDPB879G.....	106	Z845G.....	114
PDP310H.....	108	PDPB880G.....	105	Z877G.....	109
PDP311H.....	108	PDPB990G.....	105	Z879G.....	108
PDP320H.....	107	PDPB991G.....	105	Z880G.....	108
PDP332H.....	108	Z00JG.....	109	Z881G.....	109
PDP333H.....	108	Z114H.....	107	Z968H.....	111
PDP334H.....	108	Z117H.....	107	Z969H.....	111
PDP370T.....	108	Z123H.....	107	Z970H.....	111
PDP371T.....	108	Z126H.....	107	Z990G.....	108
PDP420G.....	109	Z127H.....	107	Z991G.....	108
PDP421H.....	110	Z135H.....	107	Z997G.....	110
PDP422H.....	110	Z137H.....	109	ZB340.....	106
PDP423H.....	110	Z148H.....	107	ZB345.....	106
PDP432H.....	107	Z149H.....	107	ZB346.....	106
PDP441H.....	110	Z155H.....	107	ZB347.....	106
PDP442H.....	110	Z195T.....	111	ZB370.....	106
PDP443H.....	110	Z303H.....	107	ZB371.....	106
PDP463G.....	112	Z304H.....	107	ZB740.....	106
PDP466H.....	111	Z305H.....	107		
PDP467H.....	111	Z310H.....	108	Surgical Gut Suture – Chromic	115-127
PDP468H.....	111	Z311H.....	108	1637G.....	124
PDP487G.....	112	Z320H.....	107	1638G.....	124
PDP489G.....	112	Z332H.....	108	1731G.....	126
PDP493G.....	112	Z333H.....	108	1744G.....	126
PDP494G.....	112	Z334H.....	108	1745G.....	126
PDP495G.....	112	Z370T.....	108	1754G.....	126
PDP496ZG.....	112	Z371T.....	108	1757G.....	126
PDP497G.....	112	Z420G.....	109	1766G.....	126
PDP503G.....	113	Z421H.....	110	1771G.....	126
PDP504G.....	113	Z422H.....	110	1791G.....	126
PDP507G.....	113	Z423H.....	110	1792G.....	126
PDP513G.....	112	Z432H.....	107	1793G.....	126
PDP534T.....	111	Z441H.....	110	1797G.....	126
PDP535H.....	111	Z442H.....	110	1798G.....	126
PDP569T.....	111	Z443H.....	110	3802T.....	119
PDP582G.....	109	Z463G.....	112	3803T.....	119
PDP662H.....	112	Z466H.....	111	3811T.....	118
PDP663H.....	112	Z467H.....	111	3812T.....	118
PDP682G.....	113	Z468H.....	111	3813G.....	118
PDP683G.....	113	Z487G.....	112	44T.....	120
PDP800G.....	107	Z489G.....	112	47T.....	116
PDP802D.....	107	Z493G.....	112	4887G.....	116
PDP823G.....	114	Z494G.....	112	48G.....	116
PDP824G.....	114	Z495G.....	112	4915H.....	119
PDP833G.....	114	Z496G.....	112	4923H.....	118
PDP844G.....	114	Z497G.....	112	4924H.....	118
PDP845G.....	114	Z503G.....	113	4925H.....	118
PDP877G.....	109	Z504G.....	113	4995H.....	117
PDP880G.....	108	Z507G.....	113	5589T.....	118
PDP881G.....	109	Z513G.....	112	636H.....	122
PDP968H.....	111	Z534T.....	111	654H.....	123
PDP969H.....	111	Z535H.....	111	723G.....	122
PDP970H.....	111	Z569T.....	111	724G.....	122
PDP990G.....	108	Z583G.....	109	744G.....	122
PDP9910G.....	108	Z662H.....	112	751G.....	122
PDPB316H.....	105	Z663H.....	112	752G.....	122
PDPB317H.....	105	Z682G.....	113	792G.....	124
PDPB332H.....	105	Z683G.....	113	793G.....	124
PDPB333H.....	105	Z690G.....	108	794G.....	124
PDPB334H.....	105	Z695H.....	110	795G.....	124
PDPB339H.....	106	Z761D.....	111	796G.....	124
PDPB340.....	106	Z762D.....	111	798G.....	125
PDPB345.....	106	Z765D.....	108	800H.....	119
PDPB346.....	106	Z776D.....	108	801H.....	119
PDPB347H.....	106	Z807T.....	113	802H.....	119
PDPB358T.....	106	Z823G.....	114	803H.....	119
PDPB359T.....	106	Z824G.....	114	805H.....	123

810H.....	118	L11G.....	127	S102H.....	133
811H.....	118	L112G.....	127	S103H.....	133
812H.....	118	L113G.....	127	S104H.....	133
813H.....	118	L114G.....	127	S105H.....	133
814H.....	118	L115G.....	127	SG2T.....	133
816H.....	123	M760G.....	120	U207H.....	129
817H.....	123	N877H.....	120		
832H.....	123	N878H.....	120	VICRYL™ (polyglactin 910) Suture.....	93-94
833H.....	123	N879H.....	120	V402G.....	94
837H.....	121	S111H.....	127	V448G.....	94
838H.....	121	S112H.....	127	V449G.....	94
839H.....	121	S113H.....	127	V450G.....	94
862H.....	120	S114H.....	127	V5495G.....	94
863H.....	120	S115H.....	127	V950G.....	94
864H.....	120	S116H.....	127	V960G.....	94
865H.....	120	S117H.....	127	V966G.....	94
882H.....	117	SG13T.....	127	V980G.....	94
883H.....	117	SG14T.....	127		
884H.....	117	SG15T.....	127	VICRYL RAPIDE™ (polyglactin 910) Suture ..	57-60
885H.....	117	U202H.....	116	VR1695.....	58
886H.....	122	U203H.....	116	VR214.....	58
887H.....	122	U204H.....	116	VR2218.....	58
888H.....	122	U213H.....	122	VR416.....	58
892H.....	117	U233H.....	121	VR417.....	58
893H.....	117	U245H.....	120	VR426.....	59
903H.....	120	U246H.....	120	VR490.....	59
904H.....	120	U247H.....	120	VR493.....	59
905H.....	120	U255H.....	121	VR494.....	59
912H.....	119	U256H.....	121	VR496.....	59
913H.....	119			VR497.....	59
914H.....	119	Surgical Gut Suture – Plain	128-133	VR834.....	60
915H.....	119	1627G.....	132	VR835.....	60
922H.....	118	1735G.....	133	VR844.....	60
923H.....	118	1772G.....	133	VR845.....	60
924H.....	118	1775G.....	133	VR916.....	58
925H.....	118	1824H.....	131	VR917.....	58
942H.....	121	1828H.....	131	VR931.....	59
943H.....	121	1915G.....	132	VR932.....	59
944H.....	121	52T.....	129	VR935.....	59
945H.....	121	53T.....	129	VR944.....	58
975H.....	123	5503T.....	129	VR945.....	58
977H.....	123	612G.....	130		
982H.....	123	665H.....	131		
983H.....	123	742G.....	130	Sutures: Nonabsorbable	134-207
990H.....	119	770G.....	131	ETHIBOND EXCEL™ Polyester Suture.....	140-151
991H.....	119	772G.....	131	10X42N.....	142
CC01G.....	118	773G.....	131	10X52N.....	147
CC02G.....	118	774G.....	131	10X53N.....	147
CC03G.....	117	775G.....	131	10X66N.....	148
CC04G.....	117	842H.....	129	10X72N.....	148
CC38G.....	118	843H.....	129	10X82N.....	144
CC40G.....	118	844H.....	129	10X86N.....	144
CC41G.....	118	846H.....	130	B409T.....	149
G121H.....	117	849G.....	130	B499T.....	149
G122H.....	117	852H.....	129	CX12D.....	143
G123H.....	117	853H.....	129	CX13D.....	143
G124H.....	117	872H.....	129	CX14D.....	143
G126H.....	118	873H.....	129	CX21D.....	144
G127H.....	118	898H.....	130	CX22D.....	144
G128H.....	118	CP53G.....	129	CX26D.....	143
G133H.....	119	G321H.....	129	CX27D.....	143
G152H.....	121	G322H.....	129	CX30D.....	145
G172H.....	117	L101G.....	133	CX31D.....	145
G173H.....	117	L102G.....	133	CX41D.....	142
G181H.....	117	L103G.....	133	CX42D.....	142
G182H.....	117	L104G.....	133	CX45D.....	143
GL30G.....	120	N862H.....	129	CX46D.....	143
GL31G.....	120	N863H.....	129	CX74T.....	149
K895H.....	116	N864H.....	129	CX75T.....	149

MB46G.....	146	X306H.....	151	X934H.....	147
MB47G.....	151	X406T.....	149	X935H.....	147
MB66G.....	148	X411H.....	143	X936H.....	147
MPX17.....	147	X412H.....	143	X937H.....	147
MPX37.....	147	X423H.....	144	X942H.....	145
MX523.....	142	X424H.....	144	X951H.....	146
MX551.....	142	X425H.....	144	X963H.....	146
MX552.....	142	X444H.....	146	X964H.....	146
MX553.....	142	X496T.....	149	X966H.....	148
MX563.....	144	X513H.....	144	X975H.....	147
MX69G.....	148	X517H.....	149	X976H.....	147
MX763.....	143	X518H.....	149	X977H.....	147
MX823.....	143	X519H.....	149	X987H.....	147
MX833.....	143	X521H.....	144	X997G.....	148
MX927.....	146	X522H.....	144		
MX944.....	146	X523H.....	144	ETHILON™ Nylon Suture.....	152-161
MX977.....	147	X524H.....	144	1603G.....	158
PX17H.....	147	X538H.....	150	1611G.....	157
PX514G.....	144	X548H.....	150	1629H.....	155
PX54H.....	147	X550H.....	142	1647G.....	157
PX62H.....	148	X551H.....	142	1660G.....	158
PX64H.....	148	X556H.....	142	1661G.....	158
PX65H.....	148	X557H.....	142	1662G.....	158
PX82H.....	142	X558H.....	142	1663G.....	158
PX83H.....	142	X559H.....	142	1663H.....	158
PX937H.....	147	X562H.....	144	1665G.....	157
PX977.....	147	X563H.....	144	1666G.....	157
PXX22N.....	146	X581H.....	141	1666H.....	157
PXX41N.....	142	X582H.....	142	1667G.....	157
PXX43N.....	142	X588H.....	141	1667ZH.....	157
PXX50N.....	147	X622H.....	150	1668G.....	157
PXX52N.....	147	X663H.....	149	1668H.....	157
PXX76N.....	144	X664H.....	149	1669H.....	157
PXX77N.....	148	X681G.....	150	1670BH.....	155
PXX79N.....	148	X683G.....	148	1674BH.....	155
PXX80N.....	142	X692G.....	150	1683H.....	156
PXX82N.....	142	X695G.....	150	1691H.....	155
PXX85N.....	144	X698G.....	150	1696G.....	157
PXX86.....	144	X706H.....	141	1697H.....	156
PXX87N.....	144	X717H.....	141	1698G.....	157
PXX92SCN.....	143	X726H.....	141	1714G.....	159
SX17H.....	147	X727H.....	141	1716G.....	160
SX54H.....	147	X728H.....	145	1740G.....	160
SXP02.....	144	X761H.....	143	1845G.....	159
SXP07.....	147	X762H.....	143	1854G.....	158
SXP22.....	146	X763H.....	143	1855G.....	158
SXP52.....	147	X806H.....	147	1856G.....	158
SXP54.....	147	X807H.....	147	1864G.....	158
SXP76.....	144	X831H.....	143	1865G.....	158
SXP77N.....	148	X832H.....	143	1866G.....	158
SXP82.....	142	X833H.....	143	1893G.....	159
SXP86.....	144	X834H.....	143	1894G.....	159
SXX43N.....	142	X843H.....	144	1895G.....	159
SXX50N.....	147	X844H.....	144	1954G.....	158
SXX53.....	147	X864H.....	145	1956G.....	158
SXX54N.....	147	X865H.....	145	1966G.....	158
SXX86N.....	144	X870H.....	142	2800G.....	153
X114H.....	141	X871H.....	142	2808G.....	153
X182H.....	151	X872H.....	142	2809G.....	153
X183H.....	151	X873H.....	142	2810G.....	153
X184H.....	151	X881H.....	141	2813G.....	153
X185H.....	151	X889H.....	141	2815G.....	153
X186H.....	151	X890H.....	141	2819G.....	153
X190G.....	145	X905H.....	148	2820G.....	153
X204H.....	141	X916H.....	147	2822G.....	153
X205H.....	141	X917H.....	147	2829G.....	153
X303H.....	151	X925H.....	146	2830G.....	153
X304H.....	151	X926H.....	146	2849G.....	154
X305H.....	151	X927H.....	146	2850G.....	153

2860G.....	153	9065G.....	160	N124T.....	168
2870G.....	153	G666G.....	157	PERMA-HAND™ Silk Suture	171-181
2881G.....	153	G667G.....	157	1588H.....	179
2890G.....	159	G695G.....	157	1639G.....	178
2892G.....	154	G696G.....	157	1675H.....	178
2894G.....	154	G697G.....	157	1676G.....	178
2897G.....	154	G698G.....	157	1677G.....	178
2898G.....	154	L880G.....	154	1678H.....	178
460T.....	156	L886T.....	153	1679H.....	178
489T.....	170	MERSILENE™ Polyester Fiber Suture.....	162-166	1684G.....	178
490T.....	170	1722G.....	165	1687G.....	179
562H.....	155	1760G.....	165	1732G.....	179
585H.....	158	1761G.....	165	1780G.....	179
593H.....	157	1762G.....	165	1925G.....	179
626H.....	156	1764G.....	165	1984G.....	179
627H.....	156	1778G.....	165	413H.....	173
628H.....	156	1779G.....	165	423H.....	174
642G.....	156	740G.....	164	424H.....	174
661G.....	155	746G.....	164	443H.....	177
661H.....	155	MR21T.....	163	482T.....	177
662G.....	155	R424H.....	163	485T.....	177
662H.....	155	R495T.....	164	486H.....	179
663G.....	155	R633H.....	163	580H.....	179
664G.....	155	R647H.....	163	583H.....	178
664H.....	155	R665H.....	163	590H.....	179
667G.....	154	R690G.....	164	591H.....	177
668G.....	155	R691G.....	164	621H.....	178
669H.....	155	R749G.....	165	622H.....	178
689G.....	157	R750G.....	165	623H.....	178
690G.....	157	R756G.....	165	624H.....	178
691G.....	157	R770G.....	165	629G.....	176
697G.....	157	R832H.....	163	632G.....	177
697H.....	157	R833H.....	163	639G.....	178
698G.....	157	RS22.....	166	640G.....	178
698H.....	157	RS23.....	166	641G.....	178
699G.....	157	NUROLON™ Nylon Suture.....	167-170	673H.....	179
699H.....	157	5424H.....	169	675H.....	176
749G.....	159	5425H.....	169	677G.....	176
7707G.....	160	5665G.....	170	678G.....	176
7711G.....	160	5666G.....	170	679H.....	176
7717G.....	160	C502D.....	168	680H.....	176
7718G.....	160	C503D.....	168	682G.....	176
7731G.....	159	C512D.....	168	683G.....	176
7756G.....	160	C513D.....	168	684G.....	176
7757G.....	160	C514D.....	168	685G.....	176
7760G.....	160	C520D.....	169	706G.....	172
7770G.....	160	C521D.....	169	707G.....	172
824G.....	154	C522D.....	169	709H.....	172
825G.....	154	C526D.....	168	714G.....	176
830H.....	154	C527D.....	168	734G.....	177
9000G.....	161	C541D.....	168	735G.....	176
9001G.....	161	C542D.....	168	736G.....	176
9002G.....	161	C545D.....	168	737G.....	177
9003G.....	161	C546D.....	168	765G.....	177
9005G.....	161	C550D.....	169	766G.....	177
9006G.....	161	C551D.....	169	768G.....	177
9007G.....	161	C553D.....	168	769G.....	177
9007G.....	161	C554D.....	168	7729G.....	175
9010G.....	161	C563D.....	169	7733G.....	179
9011G.....	161	C566D.....	169	7734G.....	175
9011G.....	161	C570D.....	169	780G.....	178
9013G.....	161	C572D.....	170	783G.....	178
9013G.....	161	C575D.....	170	786G.....	178
9016G.....	161	C576D.....	170	789G.....	178
9030G.....	161	C584D.....	168	A182H.....	181
9031G.....	161	C595T.....	170	A183H.....	181
9032G.....	161	N104T.....	168	A184H.....	181
9033G.....	161			A185H.....	181
9036G.....	161				
9038G.....	161				
9060G.....	160				
9061G.....	161				

A186H.....	181	N272H.....	175	8455H.....	194
A187H.....	181	SA10G.....	180	8519H.....	191
A302H.....	181	SA11G.....	180	8521H.....	193
A303H.....	181	SA5H.....	180	8522H.....	193
A304H.....	181	SA63H.....	180	8523H.....	193
A305H.....	181	SA64T.....	180	8526H.....	193
A306H.....	181	SA65T.....	180	8533H.....	193
A307H.....	181	SA66G.....	180	8534H.....	193
C003D.....	173	SA6H.....	180	8554H.....	191
C004D.....	173	SA73H.....	180	85550H.....	192
C012D.....	174	SA74H.....	180	8556H.....	192
C013D.....	174	SA75H.....	180	8557H.....	192
C014D.....	174	SA76G.....	180	8558H.....	192
C016D.....	174	SA77G.....	180	8559H.....	192
C017D.....	174	SA7H.....	180	8571G.....	198
C018D.....	174	SA83H.....	180	8580H.....	191
C041D.....	173	SA84H.....	180	8581H.....	191
C045D.....	174	SA85H.....	180	8592G.....	197
C053D.....	173	SA86G.....	180	8603G.....	199
C054D.....	173	SA87G.....	180	8604G.....	199
C084D.....	173	SA8H.....	180	8605G.....	199
CO21D.....	174	SA9G.....	180	8606G.....	199
CO22D.....	174			8610H.....	189
CO26D.....	173	PROLENE™ Polypropylene Suture.....	182-201	8617G.....	200
CO27D.....	173	1710G.....	201	8618G.....	200
K533H.....	174	1713G.....	198	8619G.....	200
K534H.....	174	1754G.....	200	8621H.....	198
K570H.....	173	1755G.....	201	8622H.....	198
K571H.....	173	1757G.....	200	8623H.....	198
K572H.....	173	1771G.....	200	8629H.....	198
K580H.....	172	2775G.....	183	8630G.....	200
K581H.....	172	2790G.....	183	8631G.....	200
K801H.....	172	2793G.....	183	8632G.....	200
K802H.....	172	2794G.....	183	8634G.....	200
K809H.....	172	2795G.....	183	8635G.....	200
K831H.....	174	3844G.....	195	8636G.....	200
K832H.....	174	3846T.....	198	8648G.....	199
K833H.....	174	7740G.....	200	8649H.....	198
K834H.....	174	7761G.....	201	8655G.....	199
K835H.....	174	7796G.....	201	8656G.....	199
K843H.....	175	788G.....	195	8661G.....	197
K844H.....	175	8190G.....	195	8663G.....	199
K845H.....	175	8204H.....	192	8665G.....	197
K846H.....	175	8205ZH.....	192	8675H.....	198
K852H.....	175	8206H.....	191	8680G.....	199
K870H.....	173	8301H.....	188	8681G.....	199
K871H.....	173	8304H.....	188	8682G.....	199
K872H.....	173	8305H.....	188	8682H.....	199
K873H.....	173	8307H.....	190	8683G.....	197
K880H.....	172	8318H.....	190	8684G.....	198
K881H.....	172	8321H.....	190	8685H.....	198
K882H.....	172	8325H.....	190	8686G.....	199
K883H.....	172	8356H.....	192	8687H.....	199
K889H.....	172	8357H.....	192	8689H.....	198
K890H.....	172	8381H.....	191	8690H.....	198
K891H.....	172	8410H.....	193	8695G.....	199
K945G.....	175	8411H.....	193	8696G.....	199
K952H.....	175	8412H.....	193	8697G.....	199
K964H.....	175	8417H.....	193	8698G.....	199
LA53G.....	180	8418H.....	193	8699G.....	199
LA54G.....	180	8422H.....	194	8701H.....	188
LA55G.....	180	8423H.....	194	8702H.....	188
M104T.....	173	8424H.....	194	8703H.....	188
M13T.....	174	8425H.....	194	8704H.....	196
M3T.....	173	8433H.....	194	8705H.....	195
M72T.....	175	8434H.....	194	8706H.....	190
N266H.....	173	8435H.....	194	8707H.....	196
N267H.....	172	8444H.....	197	8708H.....	189
N271H.....	175	8454H.....	194	8709H.....	188

8710H.....	191	C826G.....	193	EPM8734.....	184
8711H.....	191	C827G.....	193	EPM8735.....	185
8713H.....	191	C841G.....	192	EPM8736.....	184
8714H.....	191	C845G.....	193	EPM8737.....	185
8717H.....	190	EP8610H.....	189	EPM8738.....	183
8718H.....	190	EP8701H.....	188	EPM8739.....	184
8719H.....	196	EP8702H.....	188	EPM8741.....	184
8720ZH.....	190	EP8703H.....	188	EPM8742.....	184
8721H.....	196	EP8704H.....	196	EPM8743.....	184
8725H.....	190	EP8705H.....	196	EPM8745.....	187
8726H.....	190	EP8706H.....	190	EPM8746.....	186
8730H.....	183	EP8707H.....	196	EPM8747.....	187
8732H.....	184	EP8708H.....	189	EPM8751.....	187
8734H.....	184	EP8709H.....	188	EPM8752.....	186
8735H.....	185	EP8711H.....	191	EPM8753.....	186
8740H.....	184	EP8714H.....	191	EPM8754.....	186
8741H.....	184	EP8718H.....	190	EPM8755.....	187
8746H.....	186	EP8719H.....	196	EPM8765.....	185
8747H.....	186	EP8726H.....	190	EPM8766.....	185
8752H.....	186	EP8730H.....	183	EPM8767.....	185
8753H.....	186	EP8732H.....	184	EPM8768.....	185
8754H.....	186	EP8734H.....	184	EPM8775.....	185
8755H.....	187	EP8735H.....	185	EPM8776.....	189
8756H.....	192	EP8740H.....	184	EPM8785.....	185
8757H.....	192	EP8741H.....	184	EPM8800.....	189
8762H.....	193	EP8746H.....	186	EPM8805.....	188
8766H.....	185	EP8747H.....	186	EPM8806.....	188
8776H.....	189	EP8752H.....	186	EPM8807.....	196
8800H.....	189	EP8753H.....	186	EPM8812SL.....	197
8805H.....	188	EP8754H.....	186	HS6822H.....	193
8806H.....	188	EP8755H.....	187	HS6855H.....	193
8807H.....	196	EP8766H.....	185	HS6856H.....	192
8824G.....	194	EP8767LG.....	185	HS6857H.....	192
8825G.....	194	EP8776H.....	189	HS6861H.....	192
8831H.....	193	EP8800H.....	189	HS8556H.....	192
8832H.....	193	EP8805H.....	188	M8203.....	188
8833H.....	193	EP8806H.....	188	M8205.....	188
8834H.....	193	EP8807H.....	196	M8206.....	190
8842H.....	194	EP8889H.....	189	M8220.....	190
8843H.....	194	EPH7477H.....	190	M8301.....	188
8845G.....	195	EPH8710H.....	191	M8303.....	188
8853H.....	194	EPH8716H.....	191	M8304.....	188
8860H.....	191	EPH8720H.....	190	M8305.....	188
8861H.....	191	EPH8725H.....	190	M8307.....	190
8862H.....	191	EPH8890H.....	189	M8318.....	190
8870H.....	192	EPH9702H.....	188	M8321.....	190
8871H.....	192	EPM8203.....	188	M8325.....	190
8872H.....	192	EPM8205.....	188	M8556.....	192
8881H.....	191	EPM8206.....	190	M8557.....	192
8889H.....	189	EPM8610.....	189	M8600.....	189
8890H.....	189	EPM8701.....	188	M8610.....	189
8925H.....	196	EPM8702.....	188	M8701.....	188
8926H.....	196	EPM8703.....	188	M8702.....	188
8934H.....	197	EPM8704.....	196	M8703.....	188
8935H.....	197	EPM8704SL.....	197	M8704.....	196
8936H.....	197	EPM8705.....	195	M8705.....	195
8937H.....	197	EPM8706.....	190	M8706.....	190
8975H.....	197	EPM8707.....	196	M8707.....	196
8976H.....	197	EPM8709.....	188	M8709.....	188
8977H.....	197	EPM8711.....	191	M8710.....	191
8978H.....	197	EPM8718.....	190	M8711.....	191
9090G.....	201	EPM8719.....	196	M8717.....	190
9091G.....	200	EPM8723.....	196	M8718.....	190
9556H.....	192	EPM8726.....	190	M8719.....	196
9702H.....	188	EPM8727.....	189	M8720.....	190
9706H.....	191	EPM8730.....	183	M8721.....	196
C811G.....	194	EPM8731.....	185	M8725.....	190
C821G.....	194	EPM8732.....	184	M8726.....	190
C822G.....	194	EPM8733.....	185	M8727.....	189

M8730.....	183	M3768.....	136	TPW92.....	206
M8731.....	185	M3770.....	136		
M8732.....	184	PN3305H.....	137		
M8733.....	185	PN3522H.....	138	Suture Cart.....	214
M8734.....	184	PN3523H.....	138		
M8735.....	185	PN3557H.....	138		
M8736.....	184	PN3558H.....	138		
M8737.....	185	PN3581H.....	138		
M8738.....	183	PN3701H.....	137	Suture Charts.....	214-216
M8739.....	184	PN3702H.....	137		
M8741.....	184	PN3703H.....	137		
M8742.....	184	PN3709H.....	137	Absorbable Sutures Chart.....	215
M8743.....	184	PN3717H.....	137	Nonabsorbable Sutures Chart.....	216
M8745.....	187	PN3720H.....	137		
M8746.....	186	PN3725H.....	137		
M8747.....	187	PN3805H.....	137		
M8752.....	186	PN3806H.....	137		
M8753.....	186	PN3976.....	139	True-to-Scale Needle Diagram.....	29-33
M8754.....	186	PNM3703.....	137		
M8755.....	187	PNM3706.....	137	Blunt Point.....	33
M8764.....	185	PNM3725.....	137	Cutting Edge.....	30
M8765.....	185	PNM3805.....	137	Endoscopy.....	33
M8766.....	185			TAPERCUT™.....	33
M8767.....	185	Surgical Stainless Steel Suture.....	202-207	Taperpoint.....	32
M8768.....	185	467G.....	204	Trocar.....	33
M8775.....	185	497G.....	204		
M8776.....	189	603G.....	204		
M8785.....	185	607G.....	204		
M8800.....	189	618G.....	204		
M8805.....	188	631G.....	204		
M8806.....	188	DS18.....	207		
M8807.....	196	DS20.....	207		
P8522H.....	193	DS22.....	207		
P8935H.....	197	DS23.....	207		
PB20H.....	183	DS24.....	207		
PB21H.....	183	DS25.....	207		
PB30H.....	183	DS26.....	207		
		DS28.....	207		
PRONOVA™ Poly (Hexafluoropropylene – VDF)		DS30.....	207		
Suture.....	135-139	DS32.....	207		
3304H.....	137	DS35.....	207		
3307H.....	137	M400G.....	204		
3321H.....	137	M414G.....	203		
3521H.....	138	M435G.....	203		
3526H.....	138	M436G.....	205		
PN3555H.....	138	M645G.....	206		
PN3556H.....	138	M649G.....	203		
3580H.....	138	M650G.....	203		
PN3706H.....	137	M651G.....	203		
3707H.....	139	M652G.....	205		
PN3710H.....	138	M653G.....	205		
PN3711H.....	138	M654G.....	205		
PN3726H.....	137	M655G.....	205		
PN3735H.....	136	M656G.....	203		
3762H.....	138	M657G.....	205		
3776H.....	137	M660G.....	203		
3842H.....	139	M682G.....	203		
3975H.....	139	SS26G.....	206		
3976H.....	139	SS28G.....	206		
PNM3702.....	137	SS30G.....	206		
M3709.....	137	TPW10.....	205		
M3720.....	137	TPW20.....	205		
M3735.....	136	TPW30.....	206		
M3737.....	136	TPW32.....	206		
M3739.....	136	TPW40.....	205		
M3743.....	136	TPW42.....	205		
M3747.....	136	TPW72.....	206		
M3766.....	136	TPW80.....	206		
M3767.....	136	TPW90.....	206		

Place an Order

Contact your Ethicon Representative
or call 1-877-ETHICON (384-4266) to order.

More Information

To find the right product for your wound closure needs, visit [Ethicon.com](https://www.ethicon.com).

Ethnicity • Culture • Language • Religion • Values • Attitudes • Beliefs • Behaviors