

# MICRO-TOUCH® NITRA-TEX™ E.P.

Performance Series

## Nitrile, powder-free, extended-cuff medical exam glove with strong barrier protection

- **Comfortable fit:** This non-sterile glove is made from a softer synthetic formulation for a more comfortable fit and feel, and offers high tear resistance for optimal protection
- **Easy donning:** The inner chlorinated surface guarantees easy donning
- **Extended protection:** With its extended cuff, this glove is suitable for scenarios where additional wrist protection is required
- **Chemical resistance:** These gloves have excellent resistance to a broad range of commonly used chemicals
- **Reduced allergy risk:** These nitrile gloves also protect from latex Type I allergies
- **Specialized protection:** Tested and approved for use with chemotherapy drugs\*

\*see product packaging for details

### Recommended For

- Extra protection over the wrist and arm
- Sustained, high-risk examination procedures
- Environmentally friendly facilities
- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients

### FEATURED TECHNOLOGIES



### MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH® NITRA-TEX™ E.P.

Performance Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Extended
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.10 / 3.9
Finger Thickness (mm/mil)	0.12 / 4.7	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	295 / 11.6
Product Segmentation	Performance Series		

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	700121	700122	700123	700124	700125

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10993-10:2002, ISO 9001, ISO 13485, Food Contact, UKCA

## Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)	Zinc Diethyldithiocarbamate (ZDEC)
------------------------------------	------------------------------------

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com