

MICRO-TOUCH® Nitrile (White)

Standard Thin

White nitrile medical exam glove with textured fingertips

- **Heightened tactility:** Due to its lightweight design, this glove provides enhanced tactility and dexterity
- **Reduced allergy risk:** Nitrile material prevents Type I (latex) allergies
- **Smart packaging:** 150-count packaging reduces waste
- **Secure grip:** Textured fingertips ensures secure handling of medical instruments
- **Easy donning and doffing:** The polymer inner coating reduces friction for easy donning and doffing
- **Specialized protection:** Tested and approved for use with chemotherapy drugs*




*see product packaging for details

Recommended For

- Tested for use with chemotherapy drugs
- Standard, moderate risk examination procedures
- Environmentally friendly facilities
- Protection from Type I latex allergy in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® Nitrile (White)

Standard Thin

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.07 / 2.8
Finger Thickness (mm/mil)	0.10 / 3.9	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Standard Thin

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	700151	700152	700153	700154	700155

PACKAGING AND STORAGE

Packaging	150 gloves per dispenser 10 dispensers per case 1500 gloves per case *120 gloves per dispenser (Size XL) * 1200 gloves per case (Size XL)
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 455 1-4, EN 421:2010, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600-US Food Contact Approved, ISO 9001, ISO 10993-10:2002, ISO 13485, Food Type: Rubber, UKCA

Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)	
Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.	

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com